

PrecisionWeb

Point of Care Data Management System

User's Manual

QC Manager 3.0



PrecisionWeb, Precision PCx and i-STAT are trademarks of the Abbott Group of Companies.
Microsoft, Windows, Microsoft Excel ® are either registered trademarks or trademarks of
Microsoft Corporation.

Abbott Diabetes Care Inc., 1360 South Loop Road, Alameda, CA 94502, USA
© 2009 Abbott
ART12275 Rev. B 04/09

 **Abbott**

| | |
|---|------------|
| 1. ABOUT THE PRECISIONWEB POINT OF CARE DATA MANAGEMENT SYSTEM | 1-1 |
| 2. SCREEN OVERVIEW | 2-1 |
| USER LOGON SCREEN | 2-1 |
| WARNING MESSAGES | 2-2 |
| MAIN SCREEN DETAILS | 2-3 |
| THE MENU BAR | 2-4 |
| THE LOCATION TREE | 2-5 |
| THE BREADCRUMB TRAIL | 2-8 |
| THE TOTAL ALARM COUNT | 2-9 |
| THE DETAILS WORKING AREA | 2-10 |
| USER MANUAL | 2-19 |
| 3. DASHBOARD | 3-1 |
| DASHBOARD OVERVIEW | 3-1 |
| DASHBOARD DISPLAY CATEGORIES | 3-1 |
| ALARM TYPES | 3-3 |
| LIS/TEST ALARMS | 3-3 |
| NOTIFICATION ALARMS | 3-5 |
| OPERATOR | 3-7 |
| OTHER DASHBOARD SCREEN FEATURES | 3-8 |
| MY VIEWS | 3-8 |
| INSTRUMENT SUMMARY | 3-8 |
| LOCATION SCOPE | 3-9 |
| SCREEN LINKS | 3-9 |
| LIS/TEST ALARMS | 3-10 |
| NOTIFICATION ALARMS | 3-10 |
| OPERATOR ALARMS | 3-11 |
| INSTRUMENT SUMMARY | 3-14 |
| MY VIEWS | 3-14 |
| 4. REVIEW | 4-1 |
| REVIEW OVERVIEW | 4-1 |
| SUMMARY SCREEN | 4-1 |
| SEARCH CRITERIA | 4-2 |
| TABLE | 4-3 |
| E-MAIL A SUMMARY RECORD | 4-3 |
| ALARMS SCREEN | 4-5 |
| SEARCH CRITERIA | 4-6 |
| EDIT ALARM RECORDS | 4-6 |
| ACKNOWLEDGE ALARM RECORDS | 4-9 |
| E-MAIL ALARM RECORDS | 4-10 |
| ALL RESULTS SCREEN | 4-11 |
| SEARCH CRITERIA | 4-12 |
| EXPORT CURRENT VIEW | 4-14 |
| EXCLUDE RESULTS | 4-19 |
| REVIEW RESULTS | 4-21 |
| EDIT RESULTS | 4-22 |
| VIEW RESULTS | 4-29 |
| RE-ROUTE RESULT DATA | 4-30 |
| ACKNOWLEDGE RESULTS | 4-32 |

| | |
|--|------------|
| E-MAIL RESULTS | 4-33 |
| PATIENT RESULTS SCREEN | 4-35 |
| SEARCH CRITERIA | 4-36 |
| E-MAIL PATIENT RESULT RECORDS | 4-37 |
| PATIENT LOOKUP | 4-38 |
| 5. REPORTS | 5-1 |
| REPORTS OVERVIEW | 5-1 |
| VIEW REPORTS SCREEN | 5-3 |
| EDIT A REPORT | 5-3 |
| VIEW A REPORT | 5-4 |
| PRINT A REPORT | 5-5 |
| E-MAIL A REPORT | 5-5 |
| SAVE A REPORT | 5-7 |
| REPORT IMPORT SCREEN | 5-8 |
| REPORTS | 5-9 |
| AUDIT TRAIL LOG | 5-12 |
| COMMENT CODE LOG | 5-13 |
| CONTROL DATA STATISTICS BY INSTRUMENT | 5-15 |
| CONTROL DATA STATISTICS BY REAGENT LOT | 5-17 |
| DATA LISTING | 5-19 |
| GLUCOSE METER NURSE MANAGER REPORT | 5-21 |
| GLUCOSE METER QI REPRESENTATIVE REPORT | 5-23 |
| INSTRUMENT CONFIGURATION REPORT | 5-24 |
| INSTRUMENT ERROR LOG | 5-25 |
| INSTRUMENT WORKLOAD STATISTICS | 5-26 |
| LEVEY-JENNINGS GRAPH | 5-27 |
| LEVEY-JENNINGS GRAPH BY REAGENT LOT | 5-29 |
| LINEARITY GRAPH | 5-30 |
| OPERATOR CERTIFICATION EXPIRATION REPORT | 5-31 |
| OPERATOR COMPETENCY REPORT | 5-33 |
| PATIENT OUTLIER LOG | 5-35 |
| PROFICIENCY DATA LISTING REPORT | 5-37 |
| QC OUTLIER LOG | 5-39 |
| I-STAT DATA LISTING | 5-41 |
| I-STAT INSTRUMENT EVENT | 5-43 |
| QUALITY CHECK CODES | 5-44 |
| QUALITY CHECK CODES BY OPERATOR | 5-45 |
| REAGENT USAGE | 5-46 |
| 6. OPERATORS | 6-1 |
| OPERATORS OVERVIEW | 6-1 |
| NEW OPERATOR SCREEN | 6-1 |
| ADD AND CERTIFY NEW OPERATOR | 6-1 |
| ADD A NEW OPERATOR (NO CERTIFICATION PERMISSION) | 6-5 |
| OPERATOR CERTIFICATION SCREEN | 6-8 |
| EDIT OPERATOR RECORDS | 6-10 |
| DELETE AN OPERATOR | 6-11 |
| RE(CERTIFY) AN OPERATOR | 6-13 |
| E-MAIL OPERATOR CERTIFICATION RECORDS | 6-14 |
| OPERATOR SEARCH SCREEN | 6-15 |

| | |
|---|------------|
| 7. INSTRUMENTS | 7-1 |
| INSTRUMENTS OVERVIEW | 7-1 |
| MOVING, RENAMING AND DELETING INSTRUMENTS | 7-1 |
| MOVING INSTRUMENTS | 7-2 |
| RENAMING INSTRUMENTS | 7-3 |
| DELETING INSTRUMENTS | 7-4 |
| WARNING MESSAGES | 7-5 |
| INSTRUMENT CONFIGURATION SCREENS | 7-6 |
| ALL TESTS OPTIONS | 7-10 |
| PATIENT TEST OPTIONS | 7-11 |
| CONTROL TEST OPTIONS | 7-12 |
| SECURITY OPTIONS | 7-14 |
| PROFILE OPTIONS | 7-15 |
| TEST TYPE OPTIONS | 7-16 |
| SYSTEM OPTIONS | 7-18 |
| INSTRUMENTS COMMENT CODE SCREENS | 7-21 |
| ADD A COMMENT CODE | 7-22 |
| EDIT A COMMENT CODE | 7-22 |
| DELETE A COMMENT CODE | 7-23 |
| WARNING MESSAGES | 7-23 |
| INSTRUMENT PORT ASSIGNMENT SCREENS | 7-26 |
| ADD A PORT ASSIGNMENT | 7-26 |
| EDIT A PORT ASSIGNMENT | 7-27 |
| DELETE A PORT ASSIGNMENT | 7-28 |
| WARNING MESSAGES | 7-28 |
| TRACKING LOG | 7-30 |
| ADD A USER-DEFINED EVENT: | 7-31 |
| EDIT A USER-DEFINED EVENT | 7-31 |
| DELETE A USER-DEFINED EVENT | 7-32 |
| UPLOAD LOG | 7-33 |
| 8. LOTS | 8-1 |
| LOTS OVERVIEW | 8-1 |
| CONTROL LOTS SCREEN | 8-1 |
| ADD A CONTROL LOT | 8-2 |
| EDIT A CONTROL LOT | 8-3 |
| DELETE A CONTROL LOT | 8-4 |
| USING THE FILTER OPTIONS | 8-5 |
| LINEARITY LOTS SCREEN | 8-8 |
| ADD A LINEARITY LOT | 8-9 |
| EDIT A LINEARITY LOT | 8-11 |
| DELETE A LINEARITY LOT | 8-12 |
| USING THE FILTER OPTIONS | 8-13 |
| PROFICIENCY LOTS SCREEN | 8-16 |
| ADD A PROFICIENCY LOT | 8-17 |
| EDIT A PROFICIENCY LOT | 8-20 |
| DELETE A PROFICIENCY LOT | 8-22 |
| USING THE FILTER OPTIONS | 8-22 |
| WARNING MESSAGES | 8-23 |
| REAGENT LOTS SCREEN | 8-25 |
| ADD A SINGLE ANALYTE REAGENT LOT | 8-26 |

| | |
|--|------------|
| EDIT A REAGENT LOT | 8-31 |
| DELETE A REAGENT LOT | 8-32 |
| USING THE FILTER OPTIONS | 8-33 |
| WARNING MESSAGES | 8-34 |
| 9. ADMINISTRATIVE | 9-1 |
| ADMINISTRATIVE OVERVIEW | 9-1 |
| ORGANIZATION SETUP SCREEN | 9-1 |
| NAME AN ORGANIZATION | 9-1 |
| EDIT ORGANIZATION NAME AND DETAILS | 9-2 |
| WARNING MESSAGES | 9-3 |
| FACILITY SETUP SCREEN | 9-4 |
| ADD A FACILITY NAME AND DETAILS | 9-4 |
| EDIT FACILITY NAME AND DETAILS | 9-5 |
| DELETE A FACILITY | 9-5 |
| WARNING MESSAGES | 9-6 |
| DEPARTMENT SETUP SCREEN | 9-8 |
| ADD A DEPARTMENT NAME AND DETAILS | 9-8 |
| EDIT DEPARTMENT NAME AND DETAILS | 9-9 |
| DELETE A DEPARTMENT NAME | 9-10 |
| WARNING MESSAGES | 9-10 |
| LOCATION SETUP SCREEN | 9-12 |
| ADD A LOCATION NAME AND DETAILS | 9-12 |
| EDIT LOCATION NAME AND DETAILS | 9-13 |
| DELETE A LOCATION NAME AND INFORMATION | 9-13 |
| ADD A TERMINAL SERVER TO A LOCATION | 9-15 |
| EDIT TERMINAL SERVER DETAILS | 9-16 |
| DELETE A TERMINAL SERVER | 9-17 |
| WARNING MESSAGES | 9-18 |
| FACILITY SETTINGS | 9-21 |
| NEW USER SETUP SCREEN | 9-23 |
| ADD A NEW USER | 9-23 |
| IDENTIFICATION | 9-24 |
| PERMIT ACCESS TO LOCATIONS | 9-24 |
| QCM3 MENU PERMISSIONS | 9-25 |
| QCM3 FUNCTION PERMISSIONS | 9-26 |
| WARNING MESSAGES | 9-27 |
| USER SEARCH SCREEN | 9-28 |
| EDIT USERS | 9-29 |
| MY PROFILE SCREEN | 9-30 |
| SELECT A STARTUP SCREEN | 9-31 |
| DEFAULT INSTRUMENT | 9-31 |
| ALARM PREFERENCES | 9-31 |
| SHOW ALARMS OF TYPE | 9-32 |
| PASSWORD | 9-33 |
| WARNING MESSAGES | 9-33 |
| DATABASE SCREENS | 9-35 |
| DATABASE IMPORT | 9-35 |
| IMPORTING DATA FROM AN EXCEL SPREADSHEET | 9-38 |
| PANELS SCREEN | 9-40 |
| ADD A TEST PANEL | 9-41 |

| | |
|--|-------------|
| EDIT A TEST PANEL | 9-42 |
| DELETE A TEST PANEL | 9-43 |
| WARNING MESSAGES | 9-43 |
| 10. TECHNICAL SUPPORT INFORMATION | 10-1 |

1. About the PrecisionWeb Point of Care Data Management System

Thank you for selecting the PrecisionWeb Point of Care Data Management System.

Intended Use

The PrecisionWeb Data Management System is a browser-based software application that allows you to manage your Point of Care (POC) program from any computer throughout your healthcare facility. Note that the PrecisionWeb Data Management system was previously called QC Manager 3.0 and is referred to as QCM3 throughout this manual.

The PrecisionWeb data management system includes configurable features to automatically or manually forward patient and Quality Control (QC) results to the Laboratory or Hospital Information System (LIS or HIS).

In addition, a bi-directional interface allows data to be uploaded from and downloaded to any supported instrument. The bi-directional interface helps you to manage data and comply with laboratory standards and policies.

This product is not for diagnostic use; all patient diagnostics should be based on results reported by the point of care instrument.

User Interface Screens

- The User Interface (UI) screens are designed to support the management of test results, operators, instruments and lots (reagent lots, control lots, linearity lots and proficiency lots).
- The UI screens allow you to set up the facilities within your organization and name the associated departments and locations.

Multiple Users

- Multiple users can access the system at one time.
- The system allows for customized levels of access. Each level of access will provide the user with visibility to the UI screens and functions specific to the user's needs. This is referred to as custom scope and permissions.
- Authorized users can access the system from any networked computer using Windows Internet Explorer (version 5.5 or later).

Reports

- The system allows you to select 23 different default report types that can be used to create individual customized reports for your facility. Report settings can also be saved for easy retrieval at a later time.
- Each default report has specific filters available that allow you to include certain types of data in your report, such as comment codes, instrument errors, and test result outlier data.
- When a customized report is viewed later, the current data from the database will appear in the report. This eliminates the need to re-create reports that are required on a weekly or monthly basis.
- The system also gives you the ability to e-mail links to reports saved as PDF files.

We are confident that the PrecisionWeb (QCM3) Point of Care Data Management System will help you manage your POC program more effectively and efficiently.

Please read this user manual for additional information and step-by-step instructions and illustrations on how to use the system.

2. Screen Overview


User Logon Screen

The Administrator or Point of Care Coordinator (POCC) is responsible for providing each user with a User ID and Password, and assign scope and permissions to access the system.

User scope and permissions affect the menu options available, the viewable Tree locations and the actions one can perform when accessing the system.

The **User Logon** screen is the initial screen viewed by all users.

To log on to the system:



1. Enter the **User ID** and **Password** you have been assigned. (See Chapter 9, Administrative, Section [My Profile Screen](#) for information on changing your password.)
2. Click the **Logon** button, or press **Enter** on your keyboard.
3. If a valid User ID and Password are entered, the system displays the QCM3 main screen.

Warning Messages

If the User ID or Password is not entered before clicking the Logon button, one of the following warning messages is displayed:




- Click the **OK** button to clear the message and enter the missing field.

If the User ID or Password entered is not valid the following warning message is displayed:




- Click the **OK** button to clear the message and enter a valid User ID or Password

If a QCM3 "database restore" is in progress, the following warning message is displayed:



- Click the **OK** button to clear the message and log on again later.


The QCM3 session will expire if there is no activity within the application for 20 minutes. (Except for the Summary screen; see Chapter 4, Section [Summary Screen](#)). After 20 minutes, the following message appears:



- Click the **OK** button and log on again.

Main Screen Details

The main screen has eight main features:



1. The **Menu Bar** — Displays a list of the available menu options across the top of the screen. When you move your cursor over the menu items, dropdown submenus appear.
2. The **Location Tree** — Displays an expandable and collapsible list of the available locations in the organization. The Tree can remain visible or be hidden to allow more information to be displayed on the Details Working area. The Tree is located on the left side of the screen.
3. The **Breadcrumb Trail** — Displays the current menu selection, followed by the current Tree selection (organization, facility, department or location).
4. The **Total Alarm Count** — Displays the total count of unacknowledged alarms for the Location Tree scope assigned to the user. The totals include all locations accessible to the user regardless of the current Tree selection.
5. The **Details Working Area** — Displays screens that will vary based on the selected menu option. The menu and location selected in the Tree will determine the information displayed. Search Criteria may be present to help locate or filter information.
6. **Action Buttons** — These buttons vary according to the selected screen. They take the labeled action (Exclude, Review, Edit, Delete, E-Mail, etc...) to items selected in the table.
7. **Current Logged In User** — Displays the name of the current user.
8. **Log Off** — Logs off the system

The Menu Bar

The horizontal **Menu Bar** displays the options available to you depending on the scope and permission level assigned to you by the Administrator or POCC.

The following main menu items are displayed:

| | |
|-----------------------|---|
| Dashboard | Clicking this menu item will display the Dashboard screen. |
| Review | Moving the cursor over this menu item will display the Review menu. |
| Reports | Moving the cursor over this menu item will display the Reports menu. |
| Operators | Moving the cursor over this menu item will display the Operators menu. |
| Instruments | Moving the cursor over this menu item will display the Instruments menu. |
| Lots | Moving the cursor over this menu item will display the Lots menu. |
| Administrative | Moving the cursor over this menu item will display the Administrative menu. |
| Help | Moving the cursor over this menu item will display a list of available User Manuals. The table of contents will provide links to chapters and sections of the selected User Manual. If the PDF viewer is not installed, a message will ask you to install the PDF viewer. |

Selecting a Menu Option

To **view** menu options:

1. Place the mouse cursor over a menu item until the cursor changes from an arrow to a hand and a list of associated options appears.
2. If an option has additional choices, an arrow will appear to the right of the initial option.
3. Move the mouse cursor over the menu option to view the list of associated options in the submenu displayed to the right of the initial list of options.

To **select** a menu option:

1. Place the mouse cursor over a menu item to highlight the associated menu options.
2. Move the cursor (now a hand) over the list of options that appear below the menu item until the option you would like to select is highlighted.

- Click the left mouse button to select that option.

To **select** a submenu option from the **Administrative**>>Organization Setup, Administrative>>User Setup, Administrative>>Database or **Lots**>>**Reagent** options (these options have an additional submenu):

- Move the mouse cursor over **Administrative** or **Lots** to display the menu options.
- Move the mouse cursor over the **Reagent** option under **Lots** or the **Organization Setup**, **User Setup** or **Database** option under **Administrative** until the submenu appears.
- Click on the desired menu option.

The Location Tree

The Location Tree is made up of one or more Facilities and associated Departments and Locations. Each level of the Tree is linked to the level above. This is known as the Parent-Child relationship. The Organization level is a Parent Level. The Location level is a child of the Department level. Facilities and Departments can be either a Parent or Child.

| | |
|---------------|---|
| Organization: | Parent to all Facilities within the Organization |
| Facility: | Parent to all Departments within the Facility, Child to the Organization |
| Department: | Parent to all Locations within the Department, Child to the associated Facility |
| Location: | Child to the associated Department |






The **Location Tree** displays the locations available to you depending on the scope and permission level assigned to you by the Administrator or POCC.

- If you have not been assigned access to a Facility, Department or Location the corresponding location is not displayed in the Tree area.
- The Tree displays an expandable and collapsible diagram of the location structure of the organization.
- Locations are created by the Administrator or POCC and appear in the Tree in the following order:



Organization Name

The named Organization

| | | |
|---|------------------------|---|
|  | Facility Name | Associated hospitals and satellite clinics (Facilities) in the Organization |
|  | Department Name | Associated Departments in a Facility |
|  | Location Name | Associated Locations in a Department |
|  | Instrument Type | Instruments that have been added to a Location |
|  | Instruments | Associated Instruments under an Instrument Icon in a Location |

Note: *There are no limits on the number of Facilities, Departments, Locations or Instruments that can be listed in the Tree.*





Selecting from the Tree Area

To **select** a location in the Tree:

1. Click the left mouse button on Organization, Facility, Department or Location.
2. Your selection is highlighted in blue in the Tree and displayed in the Breadcrumb Trail.

For example, the screens below show:

- The highlighted Organization with the associated Facilities
- The highlighted Facility with the associated Departments
- The highlighted Department with the associated Locations
- The highlighted Location with the associated Instrument Icons and Instruments



To **expand** and **collapse** the locations in the Tree:

- Click the left mouse button on the "+" symbol to display the contents of that location.
- Locations with a "+" symbol are collapsed.
- Click the left mouse button on the "-" symbol to collapse the contents of that location.
- Locations with a "-" symbol are expanded.


To **refresh** the Tree:

1. Click the left mouse button on the **Refresh** button to refresh the Tree with the current information from the database.
2. When should you refresh the Tree? When you want to see new information, such as newly added instruments or new locations added by others.

To **hide** the Tree:

1. Click on the **Hide Tree** button in the top right corner of the Tree.
2. Clicking on the **Hide Tree** button will hide the Tree and replace it with a thin, blue vertical bar.
 - A **Show Tree** button will be displayed at the top of the blue bar.

- When the Tree is hidden, the display area to the right (Details Working Area) is enlarged and encompasses the area where the Tree had been displayed.
- Press the **Show Tree** button to display the Tree again.



To **pin** the Tree:

- When the Tree is hidden, move the cursor into the small, minimized sidebar on the left of the screen (do not click the **Show Tree** button).
 - The Tree will re-emerge with the **Pin Tree** option showing rather than the Hide Tree option.




- Pinning the Tree restores a stable view of the Tree.
- Clicking on the **Pin Tree** button reverts the button back to **Hide Tree**.

The Breadcrumb Trail


A locator line is displayed below the horizontal menu and above the Details Working Area that shows the selected menu option in teal-colored capital letters, followed by the Facility, Department and/or Location you have chosen from the Location Tree.

- The current Tree selection (Facility, Department or Location) is shown at the end of the **Breadcrumb Trail** in bold.



- If the current Tree selection is a Facility, the Breadcrumb Trail will indicate that you are looking at data from "All Departments>>All Locations."

- Similarly, if the current Tree Selection is a Facility and a Department, but not a Location, the Breadcrumb Trail will note that you are looking at data from "All Locations."




The Total Alarm Count

An **Alarm Count** status bar is displayed in the upper right-hand corner to the right of the Breadcrumb Trail and below the User ID identification.

- The **Alarm Count** is displayed on all screens except for the Logon screen.

The **Alarm Count** displays the total count of unacknowledged alarms for each of the following categories:

- Notifications (No)
- LIS/Test Alarms (LIS)
- Operator Alarms (Op)



Each **Alarm Count** will display a maximum of five digits (max: 99,999) and is calculated as the sum of all unacknowledged alarms, including all instruments for each alarm category.

If there are any unacknowledged alarms, the number is displayed in white text on red background. If there are no unacknowledged alarms, the number "0" is displayed with black text on white background.

If the number of alarms exceeds 99,999, this field will still display "99,999."

- The **Alarm Counts** that are displayed include all alarms for all the locations to which the logged-on user has access.

Notes:

- The counts are independent of the location selected in the **Location Tree**.
- The **Total Alarm Count** bar will be displayed only if you checked **Show Total Alarm Counts** under **Alarm Preferences** in the **My Profile** screen under the **Administrative** option on the horizontal menu bar.
- The alarm counts shown reflect your selection of alarm types under **Show Alarms of Type** in the **My Profile** screen.
- An explanation of alarm types can be found in Chapter 4, Review, Section [Alarms Screen](#).

The Details Working Area

The Details Working Area displays screens that will vary based the selected menu option. The menu and location selected in the Tree will determine the information displayed. Search Criteria may be present to help locate or filter information.

The screen you see immediately after you log on is determined by the option chosen from the dropdown list in the **Startup Screen** field on the **My Profile** screen, under the **Administrative** horizontal menu option.

To choose a **Startup** screen:

1. Move the mouse cursor over **Administrative** in the Menu Bar to display the menu options.
2. Move the mouse cursor down the options to **My Profile**.
3. Click on the **My Profile** option.
4. In the **My Profile** screen, under **Profile Options**, click on the arrow to the right of the **Startup Screen** field to choose from the dropdown list of screen choices.

- Click the **Save** button at the bottom of the screen.

MY PROFILE:

Profile Options

Startup Screen: (No Selection)

Default Instrument: (No Selection)

Alarm Preference

Show Alarms of T: Review > Summary, Review > Alarms, Review > All Results, Review > Patient Results, Review > Patient Lookup, Reports > View, Reports > Suites, Reports > Import, Operators > New, ☒ LIS Transfer/Generated Alarms, ☒ Operator Alarms, ☒ Test Alarms, ☒ Lot Alarms

Password

☐ Change Password

Old:

New:

Confirm:

Save

Note: The remaining choices on the My Profile screen will be addressed in Chapter 9, Administrative, Section [My Profile Screen](#).

If a user has not specified a startup screen, the **Details Working Area** displays the startup graphic.



The screen displayed in the **Details Working Area** changes depending on which menu option you have chosen from the dropdown lists on the horizontal **Menu Bar**.

When any menu option is selected from the horizontal menu, the **Details Working Area** displays:

- **Search Criteria** specific to the selected menu option (on the upper half of the screen in blue-gray background), and/or
- A system default table with fields and action buttons, which vary depending on the selected menu option and the user's assigned Function Permissions (on the bottom half of the screen).


Notes

- The **Auto Refresh** option on the **Summary** screen allows you to select how often the information in the Table is automatically updated.
- When **Auto Refresh** is selected, the Auto Logout option is overridden.
- See Chapter 4, Review: Section [Summary Screen](#) for more information on Auto Refresh.

Filtering Results

To filter the results displayed on a screen (use the **Review>>Summary** screen as an example):

1. Move your cursor to the **Search Criteria** area,
2. Select **Instrument** and **Criteria** from the dropdown lists,
3. Select a **Location** (Facility, Department or Location) from the Tree,
4. Click the **Display Records** button.



PrecisionWeb Dashboard | Review | Reports | Operators | Instruments | Lots | Administrative | Help | admin | Log Off

SUMMARY: Eastern Medical Center » All Departments » All Locations No: **134** LIS: **58** Op: **29**

Search Criteria

Instrument: **All** Criteria: **Full listing** Auto Refresh: **5 Minutes**

Display Records

2 Records Found

| Location | Instrument Type | Instrument Serial No. | Instrument Name | Connect Status | LIS Transfer | Last Good QC Date/Time | Last Upload Location | Last Upload Date/Time |
|---------------------------------|-----------------|-----------------------|-----------------|----------------|--------------|------------------------|----------------------|-----------------------|
| <input type="checkbox"/> All | | | | | | | | |
| <input type="checkbox"/> Ward A | PCx | M08990068 | M08990068 | N/A | N/A | -- | Ward A | 07/07/2004 10:01AM |
| <input type="checkbox"/> Ward B | PCx | M06210023 | M06210023 | N/A | N/A | 01/30/2003 9:15AM | Ward B | 07/07/2004 9:20AM |

2 Records Found

The information in the Table is updated based on the search criteria and location selected.

Table Features

The **Action Buttons** appear on a dark gray background along both the top and bottom of the main table in the **Details Working Area**.

Action Buttons listed on the Action Bar allow you to act on the information displayed in the table depending on the permissions assigned to you by the Administrator or POCC.

- For example, the Action Bar on the All Results table displays the **Exclude, Review, Edit, View, Re-Route, Acknowledge** and **E-Mail** action buttons.

If your User ID does not have permission to perform a certain function (i.e. edit a result record), the corresponding action button is not displayed in the action bar.

The screenshot displays the PrecisionWeb application interface. At the top, there is a navigation bar with tabs: Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. The user is logged in as 'admin' with a 'Log Off' link. Below the navigation bar, there is a search criteria section with various filters including Instrument, Patient ID, Patient Name, Patient Date Of Birth, Patient Data Source, Test, Reagent Type, Reagent Lot, Lot, Level, From, To, and Period. The main table displays 3292 records, showing columns for Location, Test Date/Time, Instrument Name, Instrument Serial No., Result Type, Reviewed, Review Comment, Operator Name, Operator ID, and Instrument Comment. The table is filtered by 'All' locations. The action bar at the top of the table includes buttons for Exclude, Review, Edit, View, Re-Route, Acknowledge, and E-Mail.


| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Result Type | Reviewed | Review Comment | Operator Name | Operator ID | Instrument Comment |
|------------------------------|--------------------|-----------------|-----------------------|---------------|----------|----------------|---------------|-------------|--------------------|
| <input type="checkbox"/> All | | | | | | | | | |
| <input type="checkbox"/> L1 | 10/16/2000 9:01AM | -- | 299999 | Electronic QC | No | -- | -- | -- | -- |
| <input type="checkbox"/> L1 | 05/07/2004 8:28AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 10:38AM | -- | 299999 | Patient | No | -- | -- | 123 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 11:21AM | -- | 299999 | Patient | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 12:00PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 12:13PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 12:36PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/07/2004 1:37PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/11/2004 8:07AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- |
| <input type="checkbox"/> L1 | 05/12/2004 9:04AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- |

Columns

Each screen table has a different set of columns that may be displayed. You can select which columns to display on most screens by using the **Customize** button.

Selecting Columns to Display

1. Click on the **Customize** button (to the right of each **Action Bar** either above or below the table).
 - This will bring up a **Select Columns to Display** pop-up screen with the list of columns available for the table associated with the menu item selected.
 - All columns currently displayed in the table will be checked in the pop-up screen.



2. Choose which columns to display in the table. You can select one or more columns. Columns not checked will not be displayed in the table.

Note: If all the columns have been selected, the **Customize** button may be to the far right of the table and not initially visible on the screen. You may need to scroll to the right using the scroll bar on the bottom of your screen to see the **Customize** button.

3. Click **OK** to close the pop-up and save the selection. Selected columns are saved as defined for each Menu option per user.
 - If you check **All**, all columns will be displayed on the screen when the pop-up closes. If you check **All** again, it will uncheck all selections and give you the error message: "Please select at least one column to display."
 - Clicking the **Cancel** button will close the pop-up without saving any selections.

Notes:

- *Selections will be saved only for the current user.*
- *This customization is available for each menu item's associated table.*

Sorting Data Columns

- Column headers with underlined hypertext can be sorted.
- Click on the underlined hypertext of a column header to sort the data in the Table in ascending order by the selected column.
- You can sort in descending order by clicking on the column header a second time.
- The sort order is indicated with an up or down arrow to the right of the column name.
- The column header field that displays an arrow is the column the table is currently sorted by.
- Selected column sort is saved as defined for each Menu option per user.

ALL RESULTS: F3 » All Departments » All Locations No: **175** LIS: **123** Op: **0**

Search Criteria

| | | |
|---|--|--|
| Instrument: <input type="text" value="PCx"/> | Patient ID: <input type="text"/> | Reagent Type: <input type="text" value="All"/> |
| Instrument Serial No.: <input type="text"/> | Patient Name: <input type="text"/> | Reagent Lot: <input type="text" value="All"/> |
| Result Type: <input type="text" value="All"/> | Patient Date Of Birth: <input type="text" value="31"/> | Lot: <input type="text" value="All"/> |
| Result Based Alarm: <input type="text" value="(No Selection)"/> | Patient Data Source: <input type="text" value="(No Selection)"/> | Level: <input type="text" value="All"/> |
| Alarm Status: <input type="text" value="All"/> | Test: <input type="text" value="All"/> | From: <input type="text" value="02/14/2006"/> |
| Status: <input type="text" value="(No Selection)"/> | Comment Status: <input type="text" value="(No Selection)"/> | To: <input type="text" value="02/14/2007"/> |
| Operator ID/Name: <input type="text"/> | Comment Code: <input type="text"/> | Period: <input type="text" value="Today"/> |

6 Records Found

| <input type="button" value="Exclude"/> | <input type="button" value="Review"/> | <input type="button" value="Edit"/> | <input type="button" value="View"/> | <input type="button" value="Re-Route"/> | <input type="button" value="Acknowledge"/> | <input type="button" value="E-Mail"/> | | | | | | |
|--|---------------------------------------|-------------------------------------|-------------------------------------|---|--|---------------------------------------|---------------|-------------|--------------------|--|--|--|
| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Result Type | Reviewed | Review Comment | Operator Name | Operator ID | Instrument Comment | | | |
| <input type="checkbox"/> All | | | | | | | | | | | | |
| <input type="checkbox"/> L3 | 01/09/2007 9:56AM | -- | M06210023 | Patient | No | -- | -- | 223 | -- | | | |
| <input type="checkbox"/> L3 | 01/10/2007 11:19AM | -- | M06210023 | Control | Yes | edit | -- | 787 | -- | | | |
| <input type="checkbox"/> L3 | 01/10/2007 11:21AM | -- | M06210023 | Control | Yes | edit | -- | 787 | -- | | | |
| <input type="checkbox"/> L3 | 01/10/2007 11:22AM | -- | M06210023 | Patient | No | -- | -- | 787 | -- | | | |
| <input type="checkbox"/> L3 | 01/10/2007 11:23AM | -- | M06210023 | Patient | No | -- | -- | 787 | -- | | | |
| <input type="checkbox"/> L3 | 01/10/2007 4:00PM | -- | M06210023 | Patient | No | -- | -- | 787 | -- | | | |

6 Records Found

Notes:

- The system stores the column sort for each table you have selected during a logged-on session.
- The next time you log onto the system the tables are displayed with your last preferred setup.

Field Highlight Display

Highlighted text will be displayed with white letters on red background in result data columns when the data is related to an alarm. Highlighted text will be displayed with black letters on yellow background in result data columns when the data is related to a previously acknowledged alarm. Highlighted text may appear on the following screens: **Summary, All Results, Edit Result, View Result and Operator Certification.**

Record Count

- A **record count** appears above the Table that displays how many records meet the chosen search criteria.
- A **page counter** appears above the Table if the search criteria selected returns more records than 100 records (up to 100 records will be displayed per page).
- For example, the table below displays the text "2178 Records Found. Displaying Page 1 of 22 << Previous 1 2 3 4 5 Next >>"

ALL RESULTS: BUMC » All Departments » All Locations No: **561** LIS: **752** Op: **5764**

Search Criteria

| | | |
|---|--|--|
| Instrument: <input type="text" value="PCx"/> | Patient ID: <input type="text"/> | Reagent Type: <input type="text" value="All"/> |
| Instrument Serial No.: <input type="text"/> | Patient Name: <input type="text"/> | Reagent Lot: <input type="text" value="All"/> |
| Result Type: <input type="text" value="All"/> | Patient Date Of Birth: <input type="text" value="31"/> | Lot: <input type="text" value="All"/> |
| Result Based Alarm: <input type="text" value="(No Selection)"/> | Patient Data Source: <input type="text" value="(No Selection)"/> | Level: <input type="text" value="All"/> |
| Alarm Status: <input type="text" value="All"/> | Test: <input type="text" value="All"/> | From: <input type="text" value="03/01/2006"/> |
| Status: <input type="text" value="(No Selection)"/> | Comment Status: <input type="text" value="(No Selection)"/> | To: <input type="text" value="03/01/2006"/> |
| Operator ID/Name: <input type="text"/> | Comment Code: <input type="text"/> | Period: <input type="text" value="Today"/> |

2178 Records Found Displaying Page 1 of 22 << Previous 1 2 3 4 5 Next >>

| Exclude | Review | Edit | View | Re-Route | Acknowledge | E-Mail | Customize | | |
|--------------------------------------|--------------------|-----------------------|-------------|-------------|--------------|---------------------|-----------|-----|--------------|
| Location | Test Date/Time | Instrument Serial No. | Patient ID | Reagent Lot | LIS Transfer | Patient Data Source | GLU | EQC | Simulator ID |
| <input type="checkbox"/> All | | | | | | | | | |
| <input type="checkbox"/> 4 Truett | 03/01/2006 12:00AM | M61490187 | 01055312003 | 67zk2g | N/A | -- | 112 mg/dL | -- | -- |
| <input type="checkbox"/> 4 Truett | 03/01/2006 12:00AM | M61630068 | 01059209001 | 67zk2g | N/A | -- | 119 mg/dL | -- | -- |
| <input type="checkbox"/> 13 Roberts | 03/01/2006 12:00AM | M61110108 | 00571220011 | 67zk2g | N/A | -- | 280 mg/dL | -- | -- |
| <input type="checkbox"/> CCU | 03/01/2006 12:01AM | M61660107 | -- | 67zk2g | N/A | -- | 48 mg/dL | -- | -- |
| <input type="checkbox"/> CCU | 03/01/2006 12:01AM | M40900045 | 01060527001 | 67zk2g | N/A | -- | 103 mg/dL | -- | -- |
| <input type="checkbox"/> CCU | 03/01/2006 12:01AM | M61660107 | -- | 67zk2g | N/A | -- | 261 mg/dL | -- | -- |
| <input type="checkbox"/> 4 North ICU | 03/01/2006 12:02AM | M61660011 | 01059155001 | 67zk2g | N/A | -- | 76 mg/dL | -- | -- |
| <input type="checkbox"/> 2 ICU | 03/01/2006 12:02AM | M61890144 | 01060061001 | 67zk2g | N/A | -- | 138 mg/dL | -- | -- |
| <input type="checkbox"/> CCU | 03/01/2006 12:03AM | M40900045 | 01060905001 | 67zk2g | N/A | -- | 201 mg/dL | -- | -- |
| <input type="checkbox"/> 11 Roberts | 03/01/2006 12:04AM | M61890017 | 00945612005 | 67zk2g | N/A | -- | 125 mg/dL | -- | -- |

To change pages:

- Click on a page number to display the group of data associated with the selected page.
- Click on **Previous** or **Next** to display the previous or next group of pages, respectively.

- Results are displayed in groups of five pages.
- If there is no next page (after the last page), the Next field is disabled.
- If there is no previous page (before the first page), the Previous field is disabled.
- The group of pages displayed will only contain numbers up to the highest page. For example, if there are three pages, the Page Counter line will read: "Previous 1 2 3 Next" with both the Previous and Next fields disabled.

User Manual

This user manual is accessible as a PDF file from the **Help** menu option on the right side of the horizontal menu. The PDF document will allow you to access specific topics from the **Table of Contents** via a hyperlink.

3. Dashboard

Dashboard Overview

The **Dashboard** screen displays an overview of summary and alarm information. The Dashboard screen provides a quick summary of test result exceptions, operator certification and instrument and connectivity status. The Dashboard screen is available only if your account was granted this permission by the Administrator or POCC. The information that appears is based on a User's assigned scope and permissions.

- The alarm counts displayed in the **Dashboard** screen reflect the selection in the **Location Tree**.
- Changing your selection in the Location Tree will change the display of alarms to those for the newly selected location.
- The Tree will only display locations for which you have been assigned permission to view.

Dashboard Display Categories

The information in the **Dashboard** screen is grouped under the following separate category headers.

LIS/TEST ALARMS

NOTIFICATIONS

OPERATOR

INSTRUMENT SUMMARY

MY VIEWS

- The **Alarm Type** and **Alarm Count** (total unacknowledged alarms up to 99,999) for each alarm type will be displayed in the rows under each category.
- Only the rows for alarm types that you have selected to view in **My Profile** under the **Administrative** menu option will be displayed.
- Even if an **Alarm Type** is turned off from the My Profile screen alarms will continue to be created as applicable, but the alarms will not be displayed under the respective category on the Dashboard screen.
- For example, if you did not select the **Test Alarms** in the My Profile screen, the Invalid QC – Inst, Out-of-Range QC, and Out-of-Range Result, Test Error and Patient ID Not Confirmed - Inst alarm rows will not be displayed under the **LIS/Test Alarms** category header, however these alarms will still be generated (for other users), just not displayed for you.

- If you do not have Review Menu Permissions, the **My Views** list is not displayed. The Review menu permission is required to create and save custom Views.

The screenshot displays the PrecisionWeb Dashboard interface. The top navigation bar includes links for Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. The user is logged in as 'admin' with a 'Log Off' option. The dashboard title is 'DASHBOARD: F3 » All Departments » All Locations'. On the left, there is an 'Organization' tree with nodes for F1, F2, D3, L3, I-STAT, and PCx. The main content area is divided into several sections:

- LIS/TEST ALARMS:** A table listing various alarm types and their counts.

| | |
|-------------------------------|----|
| Out-of-Range Result-LIS | 1 |
| Invalid Operator-LIS | 1 |
| Invalid Patient-LIS | 46 |
| Repeat Test | -- |
| Scripting Error | 1 |
| Out-of-Range QC | -- |
| Out-of-Range Result | -- |
| Test Rejected-LIS | 1 |
| Invalid QC-Inst | 1 |
| Result Edit Not Forwarded | 2 |
| Patient ID Not Confirmed-Inst | -- |
- NOTIFICATIONS:** A table listing various notification types and their counts.

| | |
|------------------------|----|
| New Reagent Lot | 29 |
| Expired Reagent Lot | -- |
| New Control Lot | 3 |
| Expired Control Lot | -- |
| New Linearity Lot | -- |
| Expired Linearity Lot | -- |
| New Instrument | 4 |
| Test Error | -- |
| Instrument Moved | -- |
| New Proficiency Lot | 39 |
| Reagent Exp. Warning | -- |
| Control Exp. Warning | -- |
| Linearity Exp. Warning | -- |
| Invalid Operator-Inst | 3 |
| Invalid Reagent-Inst | 3 |
| Rejected Result Alert | No |
- MY VIEWS:** A section with a 'Delete View' button and a large empty box for displaying custom views.
- OPERATOR:** A table listing operator-related notifications and counts.

| | |
|---------------------------|----|
| New Operator | 5 |
| Expired Op. Certification | -- |
| Operator Exp. Warning | -- |
- INSTRUMENT SUMMARY:** A table listing instrument-related notifications and counts.

| | |
|-------------------------|----|
| Not Connected | 4 |
| Upload Required | 4 |
| Upload Location Warning | -- |

- If there are no unacknowledged alarms for an alarm type under the respective alarm category the count column to the right of each alarm type is empty and "--" is displayed.
- If there are unacknowledged alarms, the alarms are displayed in bold white text on red background with the corresponding Alarm Type displayed as a hyperlink (blue type, underlined).

Alarm Types

In the lists below, the alarm types chosen in **My Profile** are indicated next to the alarm.

LIS/Test Alarms

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|----------------------------------|-------------|---|
| Out-of-Range Result - LIS | LIS | Alarm is generated based on information returned from the LIS/HIS. Test result submitted is out of range (patient Action/Critical Range or QC Range). |
| Invalid Operator - LIS | LIS | Alarm is generated based on information returned from the LIS/HIS. The Operator ID associated to test result was not accepted by LIS/HIS. The Operator ID should be corrected and result resubmitted. |
| Invalid Patient - LIS | LIS | Alarm is generated based on information returned from the LIS/HIS. Patient ID associated to test result was not accepted by LIS/HIS. Patient ID should be corrected and result resubmitted. |
| Repeat Test | LIS | Alarm is generated based on information returned from the LIS/HIS. The test result submitted was a repeat test. |
| Scripting Error | LIS | Alarm is generated when the interface script times out before the test result is sent to LIS/HIS. Result should be resubmitted via QCM. |
| Out-of-Range QC | Test | Alarm is generated based on the QC Range set at time of test. Alarm is created when the test result is uploaded. |
| Out-of-Range Result | Test | Alarm is generated based on Action/Critical Range set at time of test. Alarm is created when the test result is uploaded. |
| Test Rejected - LIS | LIS | Alarm is generated based on information returned from the LIS/HIS. Test result submitted was not accepted by the LIS/HIS. |
| Invalid QC - Inst | Test | Alarm is generated when the QC requirement is set to Allow Test or Warn. By default Allow test is set to an interval of 24 hours, therefore, when QC is not completed every 24 hours, the system perceives that |

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|--------------------------------------|-------------|--|
| | | the QC requirement is invalid or not satisfied even though the meter does not prompt to complete QC. A new alarm is created each time the instrument is uploaded until the QC requirement is satisfied. This alarm is applicable only to Abbott glucose monitors. |
| Result Edit Not Forwarded | LIS | Alarm is generated when a result has been modified and saved, but not resubmitted to LIS/HIS (Forward to LIS box not checked on Edit Result Record screen). The original alarm is acknowledged and the 'Result Edit Not Forwarded' alarm is created to alert that result record was changed but not sent to the LIS/HIS. If multiple changes are made this alarm is generated for each change, however, the associated alarms can be acknowledge in one step from the Edit Result Record screen. |
| Custom Alarms | LIS | Up to 10 custom alarms can be created. Alarms are created by Abbott personnel so that appropriate setup can take place to generate the respective alarm. Custom alarm(s) appear just below the Result Edit Not Forwarded alarm. |
| Patient ID Not Confirmed-Inst | Test | The instrument was unable to confirm the patient information entered, but based on configuration settings the test was allowed. This alarm is applicable only to instruments that have the Patient ID Confirmation capability. |

Notification Alarms

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|----------------------------|-------------|---|
| New Reagent Lot | New | Alarm is generated when a new reagent lot is used. Alarm is created upon initial upload if the reagent lot is not currently listed in QCM. If the lot is in QCM no alarm will be generated, as this alarm is generated only once per lot. |
| Expired Reagent Lot | Lot | Alarm is generated when a reagent lot used to run a test is past its expiration date. Alarm is created based on the expiration date entered in QCM. |

Note: Abbott glucose test strip packaging is bar-coded. If the strip lot scanned is past its expiration date, the glucose monitor will prompt "Invalid Lot, Date Expired" and thus not allow use of the expired strip lot.

| | | |
|----------------------------|-----|---|
| New Control Lot | New | Alarm is generated when a new control lot is used. Alarm is created upon initial upload if the control lot is not currently listed in QCM. If the lot is in QCM no alarm will be generated, as this alarm is generated only once per lot. |
| Expired Control Lot | Lot | Alarm is generated when a QC lot used to run a test is past its expiration date. Alarm is created based on the expiration date entered in QCM. Remember: QC bottle for glucose monitor is barcoded. If the scanned lot is used past its expiration date, glucose monitor will not allow use of the expired QC lot. |

Note: Abbott control lot packaging is bar-coded. If the control lot scanned is past its expiration date, the glucose monitor will prompt "Invalid Lot, Date Expired" and thus not allow use of the expired control lot.

| | | |
|------------------------------|------|---|
| New Linearity Lot | New | Alarm is generated when a new linearity lot is used. Alarm is created upon initial upload if the linearity lot is not currently listed in QCM. If the lot is in QCM no alarm will be generated, as this alarm is generated only once per lot. |
| Expired Linearity Lot | Lot | Alarm is generated when a linearity lot used to run a test is past its expiration date. Alarm is created based on the expiration date entered in QCM. |
| New Instrument | New | Alarm is generated when a new instrument is uploaded to QCM (instrument serial number is not currently in the active database). |
| Test Error | Test | Alarm is generated when an instrument test error is uploaded to QCM. One alarm is created for every test |

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|-------------------------------|-------------|---|
| | | error uploaded. This alarm is applicable only to Abbott glucose monitors. |
| Instrument Moved | New | Alarm is generated when an instrument is uploaded with a new test result from a different location than it was previously assigned. For this alarm to be generated the Configuration Download assignment under Instruments > Configuration must be set to Download Location. If configuration Download is set to Assigned Location or if no new test results exist when the instrument is uploaded an alarm will not be created. |
| New Proficiency Lot | New | Alarm is generated when a new proficiency lot is used. Alarm is created upon initial upload if the proficiency lot is not currently listed in QCM. If the lot is in QCM no alarm will be generated, as this alarm is generated only once per lot. |
| Reagent Exp. Warning | Warning | Alarm is generated when a reagent lot is due to expire. Alarm is created based on the expiration date entered in QCM under Lots > Reagent > Single Analyte or Multi Analyte and the Expiration Date Notification (days in advance) setting under Administrative > Facility Settings. The system generates only one alarm based on this criteria; additional alarms will not be created. This alarm is facility specific. Alarm will not be created unless expiration date is defined in QCM. |
| Control Exp. Warning | Warning | Alarm is generated when a control lot is due to expire. Alarm is created based on the expiration date entered in QCM under Lots > Control and the Expiration Date Notification (days in advance) setting under Administrative > Facility Settings. The system generates only one alarm based on this criteria; additional alarms will not be created. This alarm is facility specific. Alarm will not be created unless expiration date is defined in QCM. |
| Linearity Exp. Warning | Warning | Alarm is generated when a linearity lot is due to expire. Alarm is created based on the expiration date entered in QCM under Lots > Linearity and the Expiration Date Notification (days in advance) setting under Administrative > Facility Settings. The system generates only one alarm based on these criteria; additional alarms will not be created. This alarm is facility specific. |

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|--------------------------------|-------------|--|
| | | Alarm will not be created unless expiration date is defined in QCM. |
| Invalid Operator - INST | Operator | Alarm is generated when the Operator Not Certified requirement is set to Allow Test or Warn. When an Operator completes a test because Lockout is not enabled and the Operator has either an expired instrument certification or is not on the current instrument Operator list an alarm is generated. When the meter is uploaded an alarm is created for each test completed by the non-certified/expired operator. This alarm is applicable only to Abbott glucose monitors. |
| Invalid Reagent - INST | Lot | Alarm is generated when the Strip Lot Not on List requirement is set to Allow Test or Warn. When a strip lot is used because Lockout is not enabled and the strip lot is not on the approved strip lot download list (Download to Instruments option under Lots, Reagent menu) an alarm is generated. Upon upload an alarm is created for every instance this strip lot was used to test. This alarm is applicable only to Abbott glucose monitors. |

Note: If download list exceeds the maximum number of strip lots accepted by the glucose monitor this alarm will be created for all lots added to the download list over and above the maximum number. Expired strip lots or lots no longer in use should be removed from the download list to free up a space for new strip lots and avoid unnecessary alarms.

| | | |
|------------------------------|-----------|--|
| Rejected Result Alert | Yes or No | Alert is set to 'Yes' when inbound rejected results (uploaded test results) are present. If test result upload is successful alert will remain as 'No'. |
|------------------------------|-----------|--|

Operator

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|----------------------------------|-------------|--|
| Expired Op. Certification | Operator | Alarm is generated when individual operator certification has expired for a specific instrument type. Alarm is generated on the day the certification status changes from 'due to expire' (Op Exp Warning) to 'expired' because the certification was not updated. The 'Op Exp Warning' alarm is canceled for this Operator and the certification status changed to expire. One alarm is created per Operator for the expired instrument type certification. |

| <u>Alarm</u> | <u>Type</u> | <u>Description</u> |
|------------------------------|-------------|--|
| New Operator | New | Alarm is generated when a test is completed using an Operator ID not listed in QCM (i.e. this can occur if glucose monitor is set to Allow Test or Warn for the Operator Certification requirement). Alarm is created upon upload if the Operator ID is not currently on the Operator list. If the operator ID exists no alarm will be generated. |
| Operator Exp. Warning | Warning | Alarm is generated when individual Operator certification is due to expire for a specific instrument type. Alarm is generated based on Expiration Date Notification setting for Operator Certification under Administration > Facility Settings. Once the Operator is expired, the status changes from due to expire to expire, and an 'Expired Op Certification' alarm is created for this Operator. This alarm cancels out the 'Expiration Warning' alarm, as the Operator certification is no longer due to expire, the Operator certification has expired, thus a duplicate alarm will not exist regarding the status of this Operator's certification. One alarm is created per Operator for the instrument type certification due to expire. |

Other Dashboard Screen Features

My Views

A custom View is specific search criteria saved as a link for quick viewing from the Dashboard. Views are user specific and created from the All Results screen. Saved views will be displayed in a scrollable list box in the **My Views** section of the Dashboard screen.

Instrument Summary

A total **Alarm Count** for each **Alarm Type** will be displayed in the rows under the Instrument Summary category. Alarm Types are as follows:

| <u>Alarm</u> | <u>Description</u> |
|------------------------|--|
| Not Connected | Total number of failed terminal server or Data Repeater connections (i.e. alarms generated when there is no response from a terminal server). |
| Upload Required | Total numbers of instruments that have exceeded the Meter Download Warning Interval defined under Facility Settings and as such require upload (i.e. last instrument upload is older than the number of hours specified for the Meter Download Warning). |

| <u>Alarm</u> | <u>Description</u> |
|--------------|--------------------|
|--------------|--------------------|

Note: Meter Download Warning Interval is independent of instrument Upload Interval setting.

Upload Location Warning Total number of instruments uploaded from a location other than their assigned home location.

Location Scope

The location selected in the Location Tree will affect the alarm counts displayed for the following types of alarms:

Test Location

- Any LIS/Test Alarm plus:

| | |
|-------------------------|------------------------|
| Invalid Operator – Inst | Invalid Reagent – Inst |
| New Instrument | Instrument Moved |
| Test Error | New Reagent Lot |
| Expired Reagent Lot | New Control Lot |
| Expired Control Lot | New Linearity Lot |
| Expired Linearity Lot | New Proficiency Lot |
| New Operator Alarm | |

Operator Home Location

Operator Expiration Warning
Expired Operator Warning

Operator Upload Facility

New Operator

Reagent Location, Department or Facility

Reagent Expiration Warning

Lot Facility

Control Expiration Warning
Linearity Expiration Warning

Instrument Location

Instrument Summary Alarms

Screen Links

If there are unacknowledged alarms for an Alarm Type, and you have permission to access the linked page (All Results, Alarms, Summary and/or Operators Certification screens), the **Alarm Type** will provide a hyperlink (blue type, underlined) that will display information for the selected alarm(s) by linking to certain screens depending on the alarm category.

If there are no unacknowledged alarms or you do not have permission to view the associated screens, the **Alarm Type** name will not provide a link.

Clicking on an **Alarm Type** hyperlink will display information for the selected Alarm with the following links according to the category:

LIS/Test Alarms

Clicking on any of the hyperlinked LIS/Test Alarm rows will display the **All Results** screen with the following search criteria applied automatically:

| | |
|---------------------------|-------------------------------|
| Result-Based Alarm | Will display alarm type |
| Alarm Status | Will display "Unacknowledged" |
| All Dates | No date filter entry |

Note: The program assumes you do not want to set a date range filter, since you want to see which alarm caused the link to appear on the dashboard. If you were to filter by date range, you may not see all of the alarms. The date range should remain blank unless you want to see a specific date-range subset of the alarms.

- When clicking the LIS/Test Alarm hyperlink, the All Results screen is displayed. The filtered data is listed on the bottom half of the screen just below the Search Criteria. The selected LIS/Test Alarm for each record in the table is unacknowledged and highlighted in red. The data displayed is filtered for the facility, department or location selected in the Tree. If multiple alarms are present for a specific record, more than one item may be highlighted in red.

Notification Alarms

- Clicking on any of the linked **Notifications** rows will display the **Alarms** screen with the **Alarm** field automatically filtered for the alarm type you have selected from the **Dashboard** screen. The filtered data is listed on the bottom half of the screen just below the Search Criteria.
- The selected Alarm for each record in the table is unacknowledged. The data displayed is filtered for the facility, department or location selected in the Tree.
- There will be no link from the **Rejected Result Alert** row. Rejected Result Alert will display 'No' or 'Yes'. If the Rejected Result Alert displays 'Yes' contact Technical Support for assistance.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

ALARMS: Eastern Medical Center » All Departments » All Locations No:134 LIS:58 Op:29

Search Criteria

Instrument: All From: 31 Source: All

Alarm: New Reagent Lot To: 31

[Display Records](#)

6 Records Found

| Date/Time | Alarm | Source | Description |
|---|-----------------|--------|--|
| <input type="checkbox"/> 07/01/2004 11:54AM | New Reagent Lot | QCM | Reagent Lot:155255 InstType:PCx ReagentType:Test Strip Operator:U28/56 Location:Pediatrics,Ward A Test:U/01/2004 11:48AM |
| <input type="checkbox"/> 07/01/2004 11:54AM | New Reagent Lot | QCM | Reagent Lot:149735 InstType:PCx ReagentType:Test Strip Operator:1 Location:Pediatrics,Ward A Test:04/01/2004 10:59AM |
| <input type="checkbox"/> 07/01/2004 12:41PM | New Reagent Lot | QCM | Reagent Lot:625123 InstType:PCx ReagentType:Test Strip Operator:2032 Location:Pediatrics,Ward A Test:07/01/2004 3:42PM |
| <input type="checkbox"/> 07/01/2004 12:41PM | New Reagent Lot | QCM | Reagent Lot:625124 InstType:PCx ReagentType:Test Strip Operator:2034 Location:Pediatrics,Ward A Test:07/01/2004 3:47PM |
| <input type="checkbox"/> 07/01/2004 12:41PM | New Reagent Lot | QCM | Reagent Lot:625126 InstType:PCx ReagentType:Test Strip Operator:2038 Location:Pediatrics,Ward A Test:07/01/2004 3:54PM |
| <input type="checkbox"/> 07/01/2004 12:41PM | New Reagent Lot | QCM | Reagent Lot:625127 InstType:PCx ReagentType:Test Strip Operator:2040 Location:Pediatrics,Ward A Test:07/01/2004 3:57PM |

6 Records Found

Operator Alarms

- Clicking on any of the linked **Operator** rows (New Operator, Expired Op. Certification and Operator Exp. Warning) will display the **Operator Certification** screen.
- The **Status** field is automatically filtered for the alarm type you have selected from the **Dashboard** screen.
- The **Instrument** field will be pre-selected to "All" in order to include all instrument types.
- The filtered data is listed on the bottom half of the screen just below the Search Criteria.
- The selected Alarm for each record in the table is unacknowledged. The data displayed is filtered for the facility, department or location selected in the Tree.

New Operator Alarm

- Clicking on this link will display the **Operator Certification** screen showing new operators who have never been certified for any instrument type.
- Duplicate rows will be present for each operator; one row for each instrument type certification. For an operator to be certified for a specific instrument type the operator certification information must be entered via Edit Certification screen by selecting the (Re)Certify button. If the operator has not been certified for a specific instrument type, certification details will remain blank and the operator is not considered certified for the specified instrument type. An operator may be certified for one or more instrument types.
- The Search Criteria in the **Operator Certification** screen are automatically filled in as follows:

| | |
|--------------------|-----------------|
| Instrument | All |
| Status | Never Certified |
| From (date) | Blank |
| To (date) | Blank |

- Select the Instrument from the Search Criteria (the screen will automatically refresh, this may take a few seconds) then check the box next to the Operator(s) that require certification. Click the **Edit** button, and enter demographic information for each Operator, including a home location by selecting the respective location in the Tree.
- The demographic information will only need to be entered once per operator (even if there are multiple instrument type rows for an operator).

OPERATOR CERTIFICATION: F2 » All Departments » All Locations No: 173 LIS: 123 Op: 0

Search Criteria

Operator ID: Instrument: PCx From:

Last Name: Status: Never Certified To:

Display Records

16 Records Found

| <input type="checkbox"/> All | <input type="checkbox"/> Edit | <input type="checkbox"/> Delete | <input type="checkbox"/> (Re)Certify | <input type="checkbox"/> E-Mail | <input type="button" value="Customize"/> | | | | | |
|--|-------------------------------|---------------------------------|--------------------------------------|---------------------------------|--|-----------------|-----------|-------------|--------------|---------------|
| Operator Name | Operator ID | Email | Instrument Type | Original Cert. Date | Cert. Date | Recert Interval | Exp. Date | Auto Recert | Last Good QC | Last Pat Test |
| <input type="checkbox"/> -- | 3 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | OPERATOR1 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 2 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 3622 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 999 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 223 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 123 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | TESTOPERATOR | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 456 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 777 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 45454 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> -- | 56824 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> blondin, Dwight D | 888 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> broad, Dustin A | 666 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> ButterCup, Jeff B | 234 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |
| <input type="checkbox"/> Morris, Joe B | 111 | -- | PCx | -- | -- | -- | -- | -- | -- | -- |

Customize

16 Records Found

Expired Op. Certification Alarm

- Clicking on this link will display the **Operator Certification** screen showing operators whose certification has expired for a specific instrument type.
- There may be duplicate rows per operator if the operator has not been certified for a specific instrument type.
- The Search Criteria in the **Operator Certification** screen are automatically filled in as follows:

| | |
|--------------------|--------------|
| Instrument | All |
| Status | Expiration |
| From (date) | Blank |
| To (date) | Current date |

- Select the Instrument from the Search Criteria (the screen will automatically refresh, this may take a few seconds) then check the box next to the Operator(s) that require Re-certification, Click the **(Re)Certify** button, and enter the updated certification information.
- Operator demographics cannot be updated via the (Re)Certify button. Select the Edit button to update Operator demographics. To change the Operator's home location select a new location in the Tree. The Operator ID field is the only field that cannot be modified.

OPERATOR CERTIFICATION: Eastern Medical Center » All Departments » All Locations No: 134 LIS: 38 Op: 29

Search Criteria

Operator ID: Instrument: All From: 31
 Last Name: * Status: Expiration To: 10/04/2004 31
 Display Records

5 Records Found

| Operator Name | Operator ID | Department | Location | Email | Instrument Type | Original Cert. Date | Cert. Date | Recert Interval | Exp. Date | Auto Recert | Last Good QC | Last Pat Test |
|---|-------------|------------|----------|------------------|-----------------|---------------------|------------|-----------------|------------|-------------|--------------|---------------|
| <input type="checkbox"/> Lang, David | 0287561 | Pediatrics | Ward A | dlang@emc.org | PCx | 06/01/2004 | 06/01/2004 | 3 Mths | 09/01/2004 | No | -- | -- |
| <input type="checkbox"/> Smith, Larry | 5718262 | Pediatrics | Ward A | lsmith@emc.org | PCx | 06/01/2004 | 07/01/2004 | Specified | 10/01/2004 | No | -- | 07/06/29:37AM |
| <input type="checkbox"/> Manning, Renee | 329657 | Pediatrics | Ward B | rmanning@emc.org | PCx | 07/01/2004 | 07/01/2004 | -- | Expired | No | -- | 07/07/28:33AM |
| <input type="checkbox"/> Hamilton, Tara | 123456 | Pediatrics | Ward A | -- | PCx | 09/13/2003 | 09/13/2004 | -- | Expired | No | -- | -- |
| <input type="checkbox"/> Doe, Jane | 1234 | Pediatrics | Ward A | -- | PCx | 09/13/2003 | 09/13/2004 | -- | Expired | No | -- | -- |

5 Records Found

Operator Exp. Warning Alarm

- Clicking on this link will display the **Operator Certification** screen showing operators whose certification will expire for a specific instrument type within the amount of time specified for the facility (Operator Expiration Date Notification defined under Facility Settings).
- There may be duplicate rows per operator if the operator certification is expired for more than one instrument.
- The Search Criteria in the **Operator Exp. Warning** screen are automatically filled in as follows:

| | |
|--------------------|---|
| Instrument | All |
| Status | Expiration |
| From (date) | Next day's date |
| To (date) | Future date corresponding to Facility setting for Operator Expiration Date Notification (Days in Advance) |

- Select the Instrument from the Search Criteria (the screen will automatically refresh, this may take a few seconds) then check the box next to the Operator(s) that require Re-certification, Click the **(Re)Certify** button, and enter the updated certification information.
- Operator demographics cannot be updated via the (Re)Certify button. Select the Edit button to update Operator demographics. To change the Operator's home location select a new location in the Tree. The Operator ID field is the only field that cannot be modified.

Instrument Summary

Clicking on the **Not Connected**, **Upload Required**, or the **Upload Location Warning** rows under the **Instrument Summary** category header will display the current Instrument Summary screen with default search criteria.

The Search Criteria in the **Instrument Summary** screen are automatically filled in as follows:

Instrument All
Criteria Errors Only
Auto Refresh 5 Minutes

SUMMARY: Test Facility » All Departments » All Locations

Search Criteria

Instrument: Criteria: Auto Refresh:

[Display Records](#)

5 Records Found

| E-Mail | | Customize | | | | | | | |
|--------------------------------------|-----------------|-----------------------|-----------------|----------------|--------------|------------------------|----------------------|-----------------------|--|
| Location | Instrument Type | Instrument Serial No. | Instrument Name | Connect Status | LIS Transfer | Last Good QC Date/Time | Last Upload Location | Last Upload Date/Time | |
| <input type="checkbox"/> All | | | | | | | | | |
| <input type="checkbox"/> LOC01-ISTAT | I-STAT | 299999 | -- | 10.10.90.21 | OK | -- | LOC01-ISTAT | 03/19/2007 4:24PM | |
| <input type="checkbox"/> LOC01-ISTAT | I-STAT Strip | 299999 | -- | 10.10.90.21 | OK | -- | LOC01-ISTAT | 03/19/2007 4:24PM | |
| <input type="checkbox"/> LOC02-PCX2 | PCx | M10000004 | -- | 169.254.240.27 | OK | -- | LOC02-PCX2 | 03/19/2007 4:11PM | |

E-Mail [Customize](#)

My Views

Double-clicking a **View** will display the **All Results** screen with the name of the selected view at the beginning of the **Breadcrumb Trail** in place of the words "All Results." Creating a View will save the current table column configuration, and search criteria. However, the current Tree selection will not be saved as part of a View. (See Chapter 4, Review for more information on [My Views](#).)

Note: Date range search criteria will be saved when 'Period' is selected and a range chosen from the drop down box. Date range search criteria will not be saved if a range is entered into the 'From' and 'To' fields.

To select a View:

1. Double-click on the underlined text (hyperlink).
2. The name of the selected View from the **My Views** category will be displayed at the beginning of the **Breadcrumb Trail** at the top of the new screen (In the screen shot below it is **PCx – Patient Results**).

PrecisionWeb

Dashboard | Review | Reports | Operators | Instruments | Lots | Administrative | Help

admin | Log Off

Refresh | Hide Tree

ALARMS: Eastern Medical Center » All Departments » All Locations

No: 134 LIS: 58 Op: 80

Search Criteria

Instrument: PCx From: 05/04/2004 To: 10/04/2004 Source: All

Alarm: All

Display Records

71 Records Found

| Date/Time | Alarm | Source | Description |
|---|-----------------------|------------|--|
| <input type="checkbox"/> All | | | |
| <input type="checkbox"/> 07/01/2004 11:54AM | Invalid Operator-Inst | Instrument | Lin Lot: 1 Level: 1 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 3 Test: 04/01/2004 11:08AM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3555 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2032 Test: 07/01/2004 3:42PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3556 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2032 Test: 07/01/2004 3:44PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3557 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2032 Test: 07/01/2004 3:46PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3558 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2034 Test: 07/01/2004 3:47PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3559 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2034 Test: 07/01/2004 3:49PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3561 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2036 Test: 07/01/2004 3:52PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3562 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2038 Test: 07/01/2004 3:54PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3563 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2038 Test: 07/01/2004 3:56PM Comment: 16 |
| <input type="checkbox"/> 07/01/2004 12:41PM | Invalid Operator-Inst | Instrument | PatientId: 3564 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2040 Test: 07/01/2004 3:57PM Comment: 16 |
| <input type="checkbox"/> 07/06/2004 9:31AM | Invalid Operator-Inst | Instrument | PatientId: 3556 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2035 Test: 07/01/2004 8:28PM Comment: 16 |
| <input type="checkbox"/> 07/06/2004 9:31AM | Invalid Operator-Inst | Instrument | PatientId: 3557 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2035 Test: 07/01/2004 8:29PM Comment: 16 |
| <input type="checkbox"/> 07/06/2004 9:31AM | Invalid Operator-Inst | Instrument | PatientId: 3558 Location: Pediatrics, Ward A Instr: PCx M08990068 Operator: 2035 Test: 07/01/2004 8:31PM Comment: 16 |

To edit a View:


1. Modify the **Search Criteria** in the **All Results** screen.
2. Click the **Save To My Views** button.
3. A pop-up will appear that will display the name of the previously saved view.
4. Click the **Save** button in the pop-up to save the modified criteria under the previously saved View name.

Notes:

- If you modify the **View** name in the pop-up before saving, a new View will be saved and added into the Dashboard's list of Views.
- If you modify the table column selections using the **Customize** button at the far right end of the Action Bar prior to clicking on **Save To My Views**, the changes to the column selections will not only be displayed in the newly saved view, but whenever you return to the **All Results** screen.
- If you would like to modify the column selections for a particular View, but not for the All Results screen, after saving the View, you must click on the **Customize** button once again (after having saved the view) and change the column selections back to their original settings.
- If you do not have the **Review All Results** Menu Permissions, the **My Views** category is not displayed.

Delete a View:

The **Delete View** button under the **My Views** category header displays a **pop-up** that allows you to delete selected views.



To **delete** a View:

1. Click on the **Delete View** button.
2. Click on **All** to select all views for deletion, or
3. Select one or more options one at a time for deletion. Deselecting **All** clears all selection options.
4. Click the **OK** button to delete the selected view(s) and close the pop-up or Click the **Cancel** button to close the pop-up without deleting the selected views.

4. Review

Review Overview

The **Review** menu option provides access to five screens for use in managing instruments, test results and alarms. The **Summary** screen provides information on instrument status. The **Alarms** screen allows you to manage alarms and notifications. The **All Results** screen allows you to review and manage both patient and QC test results. The **Patient Results** screen provides a simple view of patient test results only. The **Patient Lookup** screen allows you to retrieve current patient data for any Facility in your organization. Each screen is available only if your account was granted this permission by the Administrator or POCC.

Summary Screen

Review>>Summary

The **Summary** screen allows you to view instrument summary information on the connection status, the last good "QC" performed and the last upload time and location.

Summary information can be viewed by Facility, Department or Location.

To **view** summary information:

1. Move the mouse cursor over **Review** in the Menu Bar to display the menu options.
2. Click on the Summary option.

The **Summary** screen displays the Summary table based on the location highlighted in the Tree and the system default Search Criteria.

To **view** summary information from a different location click on a **Facility, Department** or **Location** in the Tree.

SUMMARY: Test Facility » All Departments » All Locations

Search Criteria

Instrument

All

Criteria

Full listing

Auto Refresh

5 Minutes

Display Records

5 Records Found

E-Mail

Customize

| Location | Instrument Type | Instrument Serial No. | Instrument Name | Connect Status | LIS Transfer | Last Good QC Date/Time | Last Upload Location | Last Upload Date/Time |
|--------------------------------------|-----------------|-----------------------|-----------------|----------------|--------------|------------------------|----------------------|-----------------------|
| <input type="checkbox"/> All | | | | | | | | |
| <input type="checkbox"/> LOC01-ISTAT | i-STAT | 299999 | -- | 10.10.90.21 | OK | -- | LOC01-ISTAT | 03/19/2007 4:24PM |
| <input type="checkbox"/> LOC01-ISTAT | i-STAT Strip | 299999 | -- | 10.10.90.21 | OK | -- | LOC01-ISTAT | 03/19/2007 4:24PM |
| <input type="checkbox"/> LOC02-PCX2 | PCx | M10000004 | -- | 169.254.240.27 | OK | -- | LOC02-PCX2 | 03/19/2007 4:11PM |

E-Mail

Customize

Search Criteria

The available search criteria options on the **Summary** screen are:

| | | |
|-------------------|--------------|--|
| Instrument | All | System displays summary information for all instrument types. |
| | Instrument X | System only displays summary information for the selected Instrument type. |
| Criteria | Full listing | System displays all instrument records. |
| | Errors only | System only displays instruments with error information. |

Click the **Display Records** button to display the records based on the selected search criteria options.

The **Auto Refresh** option on the **Summary** screen allows you to select how often the information in the Summary table is automatically updated.

The screenshot shows the 'SUMMARY: F2' screen with search criteria set to 'Instrument: PCx' and 'Criteria: Full listing'. The 'Auto Refresh' dropdown menu is open, showing options: 'Off', '5 Minutes', '10 Minutes', and '15 Minutes'. The '5 Minutes' option is selected. Below the search criteria, a table displays 1 record found. The table has columns: Location, Instrument Type, Instrument Serial No., Instrument Name, Connect Status, LIS Transfer, Last Good QC Date/Time, Last Upload Location, and Last Upload Date/Time. The record shows 'L2' location, 'PCx' instrument type, 'M08990020' serial number, and a last upload date of '01/10/2007 3:55PM'.

The system default setting for Auto Refresh is set to update information every **5** minutes.

To select a different time interval or to turn the auto refresh option off, select from the dropdown list (Off, 5, 10 or 15 minutes).

If **Off** is selected, information is only refreshed if the **Display Records** button is clicked or when a different location in the Tree is selected.

If **Auto Refresh** is enabled, the information on the Summary screen is automatically refreshed at the specified interval in the **Auto Refresh** field. The auto refresh only occurs when the Summary screen is displayed. If the Auto Refresh option is enabled, then the application will not automatically log the user off due to inactivity.

- Selecting the Display Records button will update the information displayed on the Summary screen.
- The information displayed on the Summary screen will update in real-time on the Dashboard screen.

Table

Column Header Notes

Connect Status column displays:

- **OK** if the terminal server or data repeater terminal is responding to network requests.
- **IP Address** with red background if the terminal server or data repeater terminal failed to respond to network requests.
- **N/A** if the instrument communicates via a local serial COM port connection.

LIS Transfer column displays:

- **OK** if the last data transfer to the LIS was successful.
- **FAIL** with red background if the last data transfer to the LIS failed.
- **N/A** if there is no interface to the LIS.

Last Good QC Date/Time column displays:

- The date and time of the last successful (passing) QC done on the instrument.
- A red highlighted date of the last known successful QC if the specified QC interval requirement is exceeded.

Last Upload Location column displays:

- The last upload location column displays the name of the last upload location for that instrument.
- The location name will appear in red if the instrument's last upload location is different than its home location in the Tree.

Last Upload Date /Time column displays:

- The date and time of the last upload.
- A red highlighted date and time if the specified time interval for the Meter Download Warning in Facility Settings screen is exceeded.

E-Mail a Summary Record

To **e-mail** a single summary record or multiple summary records:


1. Select one or more records in the table by click on the box(es) in the Location column.
2. Click the **E-Mail** action button.

To **e-mail** all summary records:

1. Click on the box marked **All** in the Location column header.
2. Click the **E-Mail** action button.

The default e-mail client opens a new mail message and includes the menu option title in the subject field.

The body of the message includes a hyperlink to an .htm file that contains the selected records. Linked files will expire after 30 days.



Alarms Screen

Review>>Alarms

The **Alarms** screen allows you to view, edit, acknowledge and e-mail alarm records.

If a test result has more than one alarm, each alarm is displayed in a separate row.

To **view** alarm records:

1. Move the mouse cursor over **Review** in the Menu bar to display the menu options.
2. Click on the **Alarms** option.

The **Alarms** screen displays the Alarms table based on the location highlighted in the Tree and the system default Search Criteria.

To display alarm records from a different location using the same search criteria, click on a **Facility, Department or Location** in the Tree.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin: LogOff

ALARMS: F2 » All Departments » All Locations No: 175 LIS: 128 Op: 9

Search Criteria
 Instrument: PCx From: 02/15/2006 To: 02/15/2007 Source: All
 Alarm: All Display Records

7 Records Found

| Date/Time | Alarm | Source | Description |
|--------------------|-----------------------|------------|---|
| 01/10/2007 3:55PM | Invalid Operator-Inst | Instrument | PatientId:10000 Location:D2,L2 Instr:PCx M08990020 Operator:787 Hackett, Brian A Test:01/10/2007 3:54PM Comment: |
| 01/10/2007 10:55AM | Invalid Operator-Inst | Instrument | Control Lot:17747 Level:3 Location:D2,L2 Instr:PCx M08990020 Operator:787 Hackett, Brian A Test:01/10/2007 10:50AM Comment: |
| 01/10/2007 10:55AM | Invalid Operator-Inst | Instrument | Control Lot:17693 Level:1 Location:D2,L2 Instr:PCx M08990020 Operator:787 Hackett, Brian A Test:01/10/2007 10:47AM Comment: |
| 01/10/2007 10:55AM | New Control Lot | QCM | Control Lot:17747 InstType:PCx ReagentType:Test Strip Operator:787 Hackett, Brian A Location:D2,L2 Test:01/10/2007 10:50AM |
| 01/10/2007 10:55AM | New Control Lot | QCM | Control Lot:17693 InstType:PCx ReagentType:Test Strip Operator:787 Hackett, Brian A Location:D2,L2 Test:01/10/2007 10:47AM |
| 01/10/2007 10:55AM | New Instrument | QCM | New Instrument: PCx M08990020 Location: D2,L2 DateTime: 01/10/2007 10:55AM |
| 01/10/2007 10:55AM | New Reagent Lot | QCM | Reagent Lot:68ok3g InstType:PCx ReagentType:Test Strip Operator:787 Hackett, Brian A Location:D2,L2 Test:01/10/2007 10:53AM |

7 Records Found

Search Criteria

The available search criteria options on the **Alarms** screen are:

| Criteria | Options | Description |
|--------------------------|--------------|--|
| Instrument | All | Allow you to view records by instrument type. System displays alarm records from all instrument types. |
| | Instrument X | System displays alarm records only for the selected Instrument type. |
| Alarm | | Allow you to view records by a specific alarm type in the Alarm dropdown list. |
| | All | System displays all records of all alarm types from the Notifications list. System displays records of a single alarm type from the Notifications list if a single alarm is selected from the dropdown list. |
| Date Range Filter | | Pop-up calendar allows you to filter the alarms by date range in the From and To fields. The date currently associated with the alarm and displayed in the records table (upload date or test date) will be used for the filter. Alarms with dates on or between the entered dates will be displayed. You can also enter dates into the fields without using the calendar pop-ups. |
| Source | | Allow you to select alarms generated by all sources or by a single source. |
| | All | System displays alarms generated by all sources. |
| | Instrument | System only displays alarms generated by instruments. |
| | QCM | System only displays alarms generated by QCM3. |

Click the **Display Records** button to display the records based on the selected search criteria options.

Edit Alarm Records

To **edit** a single alarm record or multiple alarm records:

1. Select one or more records in the table by click on the box(es) in the Date/Time column.
2. Click the **Edit** action button.

Note: Only result-based alarms can be edited.

To **edit** all alarm records:

1. Click on the box marked **All** in the Date/Time column header to select all records.
 2. Click the **Edit** action button.
- The **Edit Result Record** screen appears with the data row of the first selected record to be edited. If multiple records are selected each record is displayed on a separate page.
 - This same screen is accessed via the **All Results** screen (see Chapter 4, Review, Section [Edit Results](#)).

EDIT RESULT RECORD: Eastern Medical Center No:134 LIS:58 Op:33

Revised Data

Operator ID: 073928 Patient ID: []

Department: *Pediatrics Location: Ward B

Result Type: Control Sample Source: []

* Reagent Lot: 655k1e * Lot: 13552

Instrument: (No Selection) * Level: Low

Accession: No

Change Comment: []

Review Information

☐ Mark As Reviewed ☐ Mark As Excluded

Review Comment: []

Predefined: []

Recent: ppp

Custom: []

User ID: admin Review Date/Time: 09/10/

Result Actions

☐ Forward to LIS ☐ Bypass Error Check

Alarms Generated By This Record

Invalid Operator-Inst ☐ Acknowledge Alarm

Invalid Reagent-Inst ☐ Acknowledge Alarm

Current Data

| | | | |
|----------------------------------|------------------------------|---------------------------------------|----------------------------|
| Operator ID: 073928 | Operator Name: Ayers, Steven | Instrument Serial No.: M06210023 | Instrument Name: M06210023 |
| Patient ID: -- | Patient Name: -- | Test Date/Time: 01/01/2003 9:01AM EST | Instrument Comment: -- |
| Department: Pediatrics | Location: Ward B | LIS Transfer Status: N/A | Instrument Text 1: -- |
| Result Type: Control | Sample Source: -- | LIS Receive Time: -- | Instrument Text 2: -- |
| Reagent Type: Glucose Test Strip | Reagent Lot: 655k1e | Accession Number: -- | Instrument Text 3: -- |
| Lot: 13552 | Level: Low | LIS Message: -- | Alarms: Yes |
| Reviewed: Yes | Excluded: Yes | Review Comment: ppp | |

Analyte Results

Glucose 51 mg/dl

To view the selected result records, click on the **>>** right arrow to view the **next** record or click on the **<<** left arrow to view the **previous** record.

The following entry fields can be edited:

| Field | Description |
|--------------------|--|
| Operator ID | Enter the correct Operator ID number. Entry field accepts up to a maximum of 30 characters. |
| Patient ID | Enter the correct Patient ID number. Entry field accepts up to a maximum of 30 characters. Field is disabled if the record is not associated with a Patient test. |
| Department | Click on a Department in the Tree. The read-only field for Department automatically displays the selected department. |

| Field | Description | | | | | | | | |
|---------------------------|---|----------------|---|----------------|---|------------------|---|--------------------|------------------------------------|
| Location | Click on a Location in the Tree. The read-only field for Location automatically displays the selected location. | | | | | | | | |
| Result type | <table> <tr> <td>Patient</td><td>Enables the Patient ID, Reagent Lot and Sample Source fields for editing.</td></tr> <tr> <td>Control</td><td>Enables the Lot and Level fields for editing.</td></tr> <tr> <td>Linearity</td><td>Enables the Lot and Level fields for editing.</td></tr> <tr> <td>Proficiency</td><td>Enables the Lot field for editing.</td></tr> </table> | Patient | Enables the Patient ID, Reagent Lot and Sample Source fields for editing. | Control | Enables the Lot and Level fields for editing. | Linearity | Enables the Lot and Level fields for editing. | Proficiency | Enables the Lot field for editing. |
| Patient | Enables the Patient ID, Reagent Lot and Sample Source fields for editing. | | | | | | | | |
| Control | Enables the Lot and Level fields for editing. | | | | | | | | |
| Linearity | Enables the Lot and Level fields for editing. | | | | | | | | |
| Proficiency | Enables the Lot field for editing. | | | | | | | | |
| Sample Source | Enabled only if Patient result record is selected. Select sample source from the dropdown list. | | | | | | | | |
| Reagent Lot | Enabled only if Patient result record is selected. Select reagent lot from the dropdown list. | | | | | | | | |
| Lot | Enabled only if QC, Proficiency or Linearity result record is selected. Select lot from the dropdown list. | | | | | | | | |
| Instrument Comment | Select an instrument comment code from the dropdown list. | | | | | | | | |
| Level | Enabled only if QC or Linearity result record is selected. Select Low, Normal or High for QC. Select level numbered L1 through L10 for Linearity. | | | | | | | | |
| Change Comment | Enter a comment related to the reason for editing the record. | | | | | | | | |
| Review Information | <input type="checkbox"/> Mark as Reviewed <input type="checkbox"/> Mark as Excluded | | | | | | | | |

Note: When a result record is marked as excluded, QCM3 automatically marks the record as reviewed.

- | | |
|-------------------|--|
| Predefined | Select a predefined comment from the dropdown list. This selection disables the custom comment field. |
| Custom | Enter a custom comment (up to 60 alphanumeric characters) in the entry field. This selection disables the predefined comment field. |

Note: If a custom comment is entered it will be added to the custom comment selection list. The list supports up to ten custom comments. If an eleventh comment is added, the first comment entered is removed from the selection list.

- | | | |
|-----------------------|---|---|
| Result Actions | <input type="checkbox"/> Forward to LIS | <input type="checkbox"/> Bypass Error Check |
|-----------------------|---|---|

- Click the **Save** button to save the edited record. The next selected record is displayed, or Click the **Cancel** button to cancel the edit action without saving.

Note: The original record is saved intact for audit purposes. The edited record is saved and stamped with the User ID and the date and time the record was edited.

Acknowledge Alarm Records

To **acknowledge** a single alarm record or multiple alarm records:


1. Select one or more records in the table by click on the box(es) in the Date/Time column.
2. Click the **Acknowledge** action button.

To **acknowledge** all alarm records:

1. Click on the box marked **All** in the Date/Time column header.
2. Click the **Acknowledge** action button.

Warning Message

When the **Acknowledge** action button is clicked the following warning message is displayed:



- Click the **OK** button to mark the selected record(s) as acknowledged in the database.
- Click the **Cancel** button to cancel the acknowledge action.

E-Mail Alarm Records

To **e-mail** a single alarm record or multiple alarm records:


1. Select one or more records in the table by clicking on the box(es) in the Date/Time column.
2. Click the **E-Mail** action button.

To **e-mail** all alarm records:

1. Click on the box marked **All** in the Date/Time column header.
2. Click the **E-Mail** action button.

The default e-mail client opens a new mail message and includes the menu option title in the subject field.

The body of the message includes a hyperlink to an .htm file that contains the selected records. Linked files will expire after 30 days.



All Results Screen

Review>>All Results

The **All Results** screen allows you to exclude, review, edit, view, re-route, acknowledge and e-mail patient, linearity, proficiency, control test results and electronic QC. Records that meet the specified search criteria may be viewed or exported. The search criteria may also be saved to the Dashboard My Views sections.

To **view** result records:

1. Move the mouse cursor over **Review** in the Menu Bar to display the menu options.
2. Click on the **All Results** option.

ALL RESULTS: F1 » D1 » L1 No: 460 LIS: 15 Op: 559

Search Criteria

| | | |
|------------------------------------|-------------------------------------|-------------------|
| Instrument: I-STAT | Patient ID: | Reagent Type: All |
| Instrument Serial No.: | Patient Name: | Reagent Lot: All |
| Result Type: All | Patient Date Of Birth: 31 | Lot: All |
| Result Based Alarm: (No Selection) | Patient Data Source: (No Selection) | Level: All |
| Alarm Status: All | Test: All | From: 31 |
| Status: (No Selection) | Comment Status: (No Selection) | To: 03/27/2007 31 |
| Operator ID/Name: | Comment Code: | Period: Today |

Export Save To My Views Display Records

3292 Records Found Displaying Page 1 of 33 «Previous 1 2 3 4 5 Next»

| Exclude | Review | Edit | View | Re-Route | Acknowledge | E-Mail | | | | |
|------------------------------|--------------------|-----------------|-----------------------|---------------|-------------|----------------|---------------|-------------|--------------------|--|
| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Result Type | Reviewed | Review Comment | Operator Name | Operator ID | Instrument Comment | |
| <input type="checkbox"/> All | | | | | | | | | | |
| <input type="checkbox"/> L1 | 10/16/2000 9:01AM | -- | 299999 | Electronic QC | No | -- | -- | -- | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 8:28AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 10:38AM | -- | 299999 | Patient | No | -- | -- | 123 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 11:21AM | -- | 299999 | Patient | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 12:00PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 12:13PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 12:36PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/07/2004 1:37PM | -- | 299999 | Linearity | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/11/2004 8:07AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- | |
| <input type="checkbox"/> L1 | 05/12/2004 9:04AM | -- | 299999 | Electronic QC | No | -- | -- | 1 | -- | |

- The **From** and **To** date fields are filled in by default with today's date. To display records you must first set the **From** Date.
- The **All Results** screen displays the Table based on the location highlighted in the Location Tree and the system default Search Criteria.

To display result records from a *different* location using the same search criteria, click on a **Facility**, **Department** or **Location** in the Location Tree.

- Use the Search Criteria options to display specific result records.

- The information in the Table in the bottom half of the screen is refreshed based on the Search Criteria selected and location highlighted in the Location Tree.
- Click the **Display Records** button to display the results based on the selected search criteria options. You can also display the results that meet the selected criteria by hitting the "Enter" key on your keyboard after changing the search criteria.

Search Criteria

The available **Search Criteria** options on the **All Results** screen are:

| Criteria | Option | Description |
|------------------------------|----------------------------------|---|
| Instrument | All | System displays all instrument results. |
| | Instrument X | System displays results only from the selected Instrument type. The screen refreshes automatically when you change instrument types, but does not display new records until the Display Records button is clicked. |
| Instrument Serial No. | | Entry field for instrument serial number. System only displays test results from the instrument that matches, or is similar to, the serial number entered. |
| Result type | All | System displays results from patient, control, Electronic QC, linearity and proficiency tests. |
| | Various Result Types | System only displays test results with the selected result type. |
| Result-Based Alarm | No Selection | Selection System displays all test results. |
| | All | System only displays results with associated "result-based" alarms. |
| | Various Alarm Types | System displays results based on the selected alarm type. |
| Alarm Status | All | Selection System displays status of both acknowledged and unacknowledged alarms. |
| | Acknowledged | System only displays acknowledged alarms. |
| | Unacknowledged | System only displays unacknowledged alarms. This control is disabled if "No Selection" is specified in the Result-based Alarm field. |
| Status | No Selection | System displays results regardless of their status: Pending Review, Reviewed, Excluded, Not Excluded, LIS Transferred, or LIS Not Transferred. |
| | Various Status Selections | System only displays results with the selected status type. |

| Criteria | Option | Description |
|------------------------------|---------------------------|---|
| Operator ID/Name | | Entry field for Operator ID or Name System only displays results that match or are similar to the ID or Name entered. |
| Patient ID | | Entry field for Patient ID System only displays results that match or are similar to the ID entered. |
| Patient Name | | Entry field for Patient Name System only displays results that match or are similar to the name entered. |
| Patient Date of Birth | | Entry field for Patient Date of Birth System only displays results that match or are similar to the date of birth entered. |
| Patient Data Source | No Selection | Default |
| | QCM | Indicates that the Patient ID was updated in QCM from the Edit Result screen. |
| | Instrument | Indicates that Admission, Discharge, and Transfer data for a test is set by the instrument (this Information is originally uploaded to the instrument by QCM). |
| | LIS | Indicates that Admission, Discharge, and Transfer data for a test was updated with information from the LIS. |
| Test | All | System displays all test results. |
| | Various Test Types | System only displays results from selected test type. |
| Comment Status | No Selection | Default |
| | No Comment | |
| | Any Commented | |
| | Specific | System allows you to enter a comment code number into entry field below. |
| Comment Code | | If "Specific" was chosen under Comment Status (above), you can enter a comment code number here. The other choices under Comment Status will disable the Comment Code field. |
| Reagent Type | All | System default, which displays results from any reagent product. |
| | Reagent X | System displays results only from the selected reagent. The choices available vary depending on Instrument type selected. |

| Criteria | Option | Description |
|--------------------|---------------------------------------|---|
| Reagent Lot | All | System default |
| | Various Reagent Lot Selections | Select from the dropdown list or Name entered. The choices available vary depending on Reagent Type selected. |
| Lot | | The Lot list is enabled if the Result Type search criterion selected is Control, Linearity or Proficiency. |
| Level | | The Level list is enabled if the Result Type search criterion selected is Control or Linearity. Control: Choices are All, High, Low, or Normal. Linearity: Choices are L1 to L10. |
| Date Range | | System default displays today's date. Enter or select a From and a To date to display test results that were performed between the selected dates. |
| | From | Click on the calendar icon to display the pop-up calendar and click on a date. The From date field automatically displays the selected date. |
| | To | Click on the calendar icon to display the pop-up calendar and click on a date. The To date field automatically displays the selected date – required . |

To Export Current View

- Click the **Export** button to export all records that match the current search criteria.
 - All data fields associated with each record will be exported into a comma-separated value (CSV) file. This CSV file may then be accessed with a spreadsheet program such as Microsoft Excel.
- An Internet Explorer dialogue box will appear. Click **Save** to save the data file or **Open** to open it directly in the default application for the CSV file type.
 - The default application for CSV files is set by the Windows Folder Options control panel. To set the default application, open 'Folder Options' in the Windows Control Panel.
 - Select Start > Settings > Control Panel.
 - Open the 'Folder Options' file folder.
 - Select the 'File Types' tab and scroll down to the CSV file type.
 - Verify that Microsoft Excel is selected as the default application for CSV files. If not, select the 'Change' button to set Microsoft Excel as the default application and click 'OK'.

5. Select the 'Advanced' button and verify that the "confirm open after download" option is not checked. If this box is checked, uncheck it and click 'OK'.
 6. Click 'Close'
- The default name selected is 'Export.csv'. If a file with the same name exists, then Internet Explorer functionality will append '[x]' to the file name (where 'x' is a number that makes the file name unique). These brackets in the file name may cause an error if the file is opened then saved using Microsoft Excel. A unique name should be provided to all export files.
 - The data displayed in the All Results grid and the data exported will match. Pressing the Export button will export data for the current search criteria selection and also refresh the data displayed in the All Results grid.
 - The first row of the export file will contain column headings. Test data will be exported starting on the second row.
 - Each test record in QCM will be exported in one row. For each unique result associated with a test in the dataset, additional columns will be appended. The first column heading will contain the analyte name. The result value will be exported in this column. The second column heading will contain the analyte name followed by the static text "Unit". The unit associated with the result will be exported in this column. If a test is not associated with a particular analyte, then the column values will be left empty for that analyte.
 - Each row of data contains the following thirty (30) record-specific columns:
 - *Operator ID*
 - *Operator Name*
 - *Instrument Serial No.*
 - *Instrument Name*
 - *Patient ID*
 - *Patient Name*
 - *Patient Date of Birth*
 - *Patient Gender*
 - *Patient Data Source*
 - *Test Date/Time*
 - *Accession Number*
 - *LIS Transfer Status*
 - *Instrument Comment*
 - *Department*
 - *Location*
 - *LIS Receive Time*
 - *Instrument Text 1*
 - *Instrument Text 2*
 - *Instrument Text 3*
 - *Result Type*
 - *Sample Source*
 - *Alarms*
 - *Reagent Type*
 - *Reagent Lot*
 - *Excluded*
 - *Lot*
 - *Level*
 - *Reviewed*
 - *Review Comment*
 - *LIS Message*

If the csv file is opened in Microsoft Excel the display will appear as follows:

| | A | B | C | D | E | F | G | H | I |
|---|-------------|---------------------|-----------------------|-----------------|------------|--------------|-----------------------|----------------|---------------------|
| 1 | Operator ID | Operator Name | Instrument Serial No. | Instrument Name | Patient ID | Patient Name | Patient Date of Birth | Patient Gender | Patient Data Source |
| 2 | 28365 | Unser, Nicole/NICUB | 300836 | | 0 | | | | |

| | J | K | L | M | N | O | P | Q |
|---|-------------------|------------------|---------------------|--------------------|----------------------|----------|------------------|-------------------|
| 1 | Test Date/Time | Accession Number | LIS Transfer Status | Instrument Comment | Department | Location | LIS Receive Time | Instrument Text 1 |
| 2 | 03/03/2006 1:16PM | | N/A | | Hoblitzelle Building | NICU | | |

| | R | S | T | U | V | W | X | Y | Z | AA | AB | AC |
|---|-------------------|-------------------|-------------|---------------|--------|--------------|-------------|----------|-----|-------|----------|----------------|
| 1 | Instrument Text 2 | Instrument Text 3 | Result Type | Sample Source | Alarms | Reagent Type | Reagent Lot | Excluded | Lot | Level | Reviewed | Review Comment |
| 2 | | | Patient | Unknown | No | G3+ | unspecified | No | | | No | |


| | AD | AE | AF | AG | AH | AI | AJ | AK | AL | AM |
|---|-------------|----------------|------|--------------|-----------------|---------------------|----------------|------------|-----------|-------------|
| 1 | LIS Message | BatteryVoltage | Uses | Software | PreferencesName | PreferencesRevision | ID | MethodName | PanelCode | StaroutCode |
| 2 | | 8.40V | 8968 | JAMS119J-A10 | 062099HB | 5 | 00 300836 8968 | i-STAT | 0C | 0 |

| | AN | AO | AP | AQ | AR | AS | AT | AU | AV | AW | AX | AY | AZ | BA |
|---|-----------|---------|-------|---------|--------------|----------|----------|----------|----------|------|-----------|----------------|------------|------------|
| 1 | FluidCode | ISTATID | pH | pH Unit | pH ValueType | pH RefLo | pH RefHi | pH ActLo | pH ActHi | TCO2 | TCO2 Unit | TCO2 ValueType | TCO2 RefLo | TCO2 RefHi |
| 2 | 0 | 13 | 7.293 | | Measured | 7.35 | 7.45 | 7.2 | 7.6 | 26 | mmol/L | Calculated | 24 | 29 |

| | BB | BC | BD | BE | BF | BG | BH | BI | BJ | BK | BL | BM | BN |
|---|-----|----------|---------------|-----------|-----------|-----------|------|-----------|----------------|------------|------------|------------|------------|
| 1 | PO2 | PO2 Unit | PO2 ValueType | PO2 RefLo | PO2 RefHi | PO2 ActLo | PCO2 | PCO2 Unit | PCO2 ValueType | PCO2 RefLo | PCO2 RefHi | PCO2 ActLo | PCO2 ActHi |
| 2 | 62 | mmHg | Measured | 75 | 100 | 50 | 51.2 | mmHg | Measured | 35 | 45 | 20 | 60 |

Warning Message

If you attempt to export a file that contains more than 256 columns or 65536 rows, then a dialog with the following message will appear:



- Click **OK** to continue the export.
- Click **Cancel** to close the dialog box and terminate the export.

To Save Current View

The screenshot shows the PrecisionWeb interface with the following search criteria:

- Instrument: PCx
- Instrument Serial No.:
- Result Type: All
- Result Based Alarm: (No Selection)
- Alarm Status: All
- Status: (No Selection)
- Operator ID/Name:
- Patient ID:
- Patient Name:
- Patient Date Of Birth: 31
- Patient Data Source: (No Selection)
- Test: All
- Comment Status: (No Selection)
- Comment Code:
- Reagent Type: All
- Reagent Lot: All
- Lot: All
- Level: All
- From: 02/15/2006
- To: 02/15/2007
- Period: Today

Buttons: Export, **Save To My Views** (circled in red), Display Records

6 Records Found

| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Result Type | Operator Name | Operator ID | Patient ID | Patient Name | LIS Transfer | Patient Data Source | GLU |
|------------------------------|--------------------|-----------------|-----------------------|-------------|---------------|-------------|------------|-------------------------|--------------|---------------------|-----------|
| <input type="checkbox"/> All | | | | | | | | | | | |
| <input type="checkbox"/> L3 | 01/09/2007 9:56AM | -- | M06210023 | Patient | -- | 223 | 10000 | -- | N/A | -- | 74 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:19AM | -- | M06210023 | Control | -- | 787 | -- | -- | Yes | -- | 49 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:21AM | -- | M06210023 | Control | -- | 787 | -- | -- | Yes | -- | 323 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:22AM | -- | M06210023 | Patient | -- | 787 | 10000 | Zappa10000, Frank1.0000 | Fail | LIS | 82 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:23AM | -- | M06210023 | Patient | -- | 787 | 10009 | WhoIs10009, Pat10009 | Fail | LIS | 85 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 4:00PM | -- | M06210023 | Patient | -- | 787 | 10235 | Smith10235, Betty10235 | Fail | LIS | 94 mg/dL |

6 Records Found

1. Click the **Save To My Views** button to save the table column configuration and search criteria. The date range filter and Location Tree selection will not be saved.

The dialog box titled "SAVE TO MY VIEWS" contains the following text:

Type a name for your view and click save to add it to your dashboard.

View Name:

Buttons: Save, Cancel

2. Enter a descriptive name in the View Name field of the pop-up.
3. Click **Save** to save as a view for display in the **My Views** section of the **Dashboard** screen, or Click the **Cancel** to return to the All Results screen without saving the view.

The Table displays the list of records that meet the Selection Criteria you have chosen.

Available Columns

The following columns will be available for display selection:

- Location
- Test Date/Time
- Instrument Name
- Instrument Serial Number
- Result Type

- Reviewed
- Review Comment
- Operator Name
- Operator ID
- Instrument Comment
- Patient ID
- Patient Name
- Patient Date of Birth
- Reagent Type
- Reagent Lot
- Lot
- Level
- Alarms
- Excluded
- LIS Transfer
- Patient Data Source
- Glucose
- EQC
- Simulator ID

All the columns are shown by default until you customize the table using the **Customize** button at the far right of the table's Action Bar.

Alarm Highlighting

Columns with an associated alarm will be highlighted according to alarm type.

- If a test has multiple alarms, multiple fields may be highlighted.
- A field will only be highlighted if you have selected the Alarm Type in the My Profile screen.
- The highlight will be white text on red background for unacknowledged alarms and black text on yellow background for acknowledged alarms.
- If a field caused both an acknowledged and an unacknowledged alarm, the unacknowledged alarm will take precedence and the field will be highlighted in red.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin: Log Off

ALL RESULTS: F3 » All Departments » All Locations No:175 LIS:123 Op:0

Show Tree

Search Criteria

Instrument: PCx Patient ID:
 Instrument Serial No.: Patient Name:
 Result Type: All Patient Date Of Birth: 31
 Result Based Alarm: (No Selection) Patient Data Source: (No Selection)
 Alarm Status: All Test: All
 Status: (No Selection) Comment Status: (No Selection)
 Operator ID/Name: Comment Code:

Reagent Type: All
 Reagent Lot: All
 Lot: All
 Level: All
 From: 02/15/2006 31
 To: 02/15/2007 31
 Period: Today

Export Save To My Views Display Records

6 Records Found

| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Result Type | Operator Name | Operator ID | Patient ID | Patient Name | LIS Transfer | Patient Data Source | GLU |
|------------------------------|--------------------|-----------------|-----------------------|-------------|---------------|-------------|------------|------------------------|--------------|---------------------|-----------|
| <input type="checkbox"/> All | | | | | | | | | | | |
| <input type="checkbox"/> L3 | 01/09/2007 9:56AM | -- | M06210023 | Patient | -- | 223 | 10000 | -- | N/A | -- | 74 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:19AM | -- | M06210023 | Control | -- | 787 | -- | -- | Yes | -- | 49 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:21AM | -- | M06210023 | Control | -- | 787 | -- | -- | Yes | -- | 323 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:22AM | -- | M06210023 | Patient | -- | 787 | 10000 | Zappa10000, Frank10000 | Fail | LIS | 82 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 11:23AM | -- | M06210023 | Patient | -- | 787 | 10009 | WhoIs10009, Pat10009 | Fail | LIS | 85 mg/dL |
| <input type="checkbox"/> L3 | 01/10/2007 4:00PM | -- | M06210023 | Patient | -- | 787 | 10235 | Smith10235, Betty10235 | Fail | LIS | 94 mg/dL |

6 Records Found

Action Bar

The Action Bar is displayed above and below the Table. The Action Bar includes the following Action buttons:

- **Exclude**
- **Review**
- **Edit**
- **View**
- **Re-Route**
- **Acknowledge**
- **E-Mail**

Exclude Results

Clicking this button allows you to exclude records from statistical reports.

To **exclude** a single result record or multiple result records from statistical reports:

1. Select one or more records in the table by click on the box(es) in the Location column.
2. Click the **Exclude** action button.

To **exclude** all result records from statistics:

1. Click on the box marked **All** in the Location column header.
2. Click the **Exclude** action button.

The **Exclude From Statistics** screen displays the User ID and the Review Date/Time in read-only fields and allows you to change the predefined comment or add a custom comment to the records. Result records that are excluded from statistics will not appear on reports containing statistics.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

EXCLUDE FROM STATISTICS: Test Facility1 » All Departments » All Locations No:12 LIS:4 Op:2

User ID: admin

Review Date/Time: 05/02/2007 04:03 PM

Review Comment:

Predefined ☐ 1

Most Recent ☐

Custom ☒

Save Cancel

| Location | Test Date/Time | Instrument Serial No. | Result Type | Reviewed | Operator ID | Patient ID | Reagent Type | Reagent Lot | Excluded | LIS Transfer | GLU |
|----------|--------------------|-----------------------|-------------|----------|-------------|------------|--------------------|-------------|----------|--------------|------------|
| LOC1 | 04/18/2007 11:58AM | XP0001-0005 | Control | No | 01 | -- | Glucose Test Strip | 44C15H | No | N/A | 7.7 mmol/L |
| LOC1 | 04/18/2007 11:59AM | XP0001-0005 | Control | No | 01 | -- | Glucose Test Strip | 44C15H | No | N/A | 7.7 mmol/L |
| LOC1 | 04/18/2007 12:00PM | XP0001-0005 | Patient | No | 01 | 10019 | Glucose Test Strip | 44C15H | No | Fail | 7.7 mmol/L |

Review Comment

Select from the following three choices:

| Field | Description |
|--------------------|--|
| Predefined | <p>Select from a list of 10 predefined comments from the dropdown list. Predefined comments are entered in the Facility Settings screen; see Chapter 9, Administrative, Section Facility Settings.</p> <p>This selection disables the Custom and Most Recent comment fields.</p> |
| Most Recent | <p>Select a comment from the dropdown list. Up to 10 comments will appear in this list.</p> <p>The comments listed are those most recently entered in the Custom comment field, below.</p> <p>The most recently entered comment in the Custom comment field will appear as the default entry in the Most Recent comment field until changed.</p> <p>This selection disables the Custom and Predefined comment fields.</p> |
| Custom | <p>Enter a custom comment (up to 60 alphanumeric characters) in the entry field.</p> <p>This selection disables the Predefined and Most Recent comment fields.</p> <p>When a comment is entered it will be added to the Most Recent comment selection list upon saving.</p> <p>The list supports up to ten custom comments. If an eleventh comment is added, the first comment entered is removed from the selection list.</p> |

- Click the **Save** button to mark the selected records as excluded from statistics or Click the **Cancel** button to cancel the exclude action without saving.

Review Results

Clicking this button allows you to add a comment to one or more result records. Selected results will be marked as reviewed and stamped with your PrecisionWeb User ID (electronic signature) and date and time of result review.

To **review** a single result record or multiple result records:

1. Select one or more records in the table by click on the box(es) in the Location column.
2. Click the **Review** action button.

To **review** all result records:

1. Click on the box marked **All** in the Location column header.
2. Click the **Review** action button.

The **Review Results** screen displays the User ID and Review Date/Time in read-only fields as well as the following fields that allow you to change the Review Comment or add a new comment to the records.

REVIEW RESULTS: Test Facility 1 > All Departments > All Locations

User ID: admin
 Review Date/Time: 05/02/2007 04:00 PM

Review Comment:

Predefined: ☒ 1
 Most Recent: ☐
 Custom: ☐

Save Cancel

| Location | Test Date/Time | Instrument Serial No. | Result Type | Reviewed | Operator ID | Patient ID | Reagent Type | Reagent Lot | Excluded | LIS Transfer | GLU |
|----------|--------------------|-----------------------|-------------|----------|-------------|------------|--------------------|-------------|----------|--------------|------------|
| LOC1 | 04/18/2007 11:58AM | XP0001-0005 | Control | No | 01 | -- | Glucose Test Strip | 44C15H | No | N/A | 7.7 mmol/L |
| LOC1 | 04/18/2007 11:59AM | XP0001-0005 | Control | No | 01 | -- | Glucose Test Strip | 44C15H | No | N/A | 7.7 mmol/L |
| LOC1 | 04/18/2007 12:00PM | XP0001-0005 | Patient | No | 01 | 10019 | Glucose Test Strip | 44C15H | No | Fail | 7.7 mmol/L |

Review Comment

Select from the following three choices:

| Field | Description |
|--------------------|--|
| Predefined | Select from a list of 10 predefined comments from the dropdown list. Predefined comments are entered in the Facility Settings screen; see Chapter 9, Administrative, Section Facility Settings . This selection disables the Custom and Most Recent comment fields. |
| Most Recent | Select a comment from the dropdown list. Up to 10 comments will appear in this list. |

| Field | Description |
|---------------|--|
| | <p>The comments listed are those most recently entered in the Custom comment field, below.</p> <p>The most recently entered comment in the Custom comment field will appear as the default entry in the Most Recent comment field until changed.</p> <p>This selection disables the Custom and Predefined comment fields.</p> |
| Custom | <p>Enter a comment (up to 60 alphanumeric characters) in the entry field.</p> <p>This selection disables the Predefined and Most Recent comment fields.</p> <p>When a comment is entered it will be added to the Most Recent comment selection list upon saving.</p> <p>The list supports up to ten custom comments. If an eleventh comment is added, the first comment entered is removed from the selection list.</p> <ul style="list-style-type: none"> Click the Save button to mark the selected records as reviewed or Click the Cancel button to cancel the review action without saving. |

Edit Results

Clicking this button allows you to edit data for one or more result records.

To **edit** a single result record or multiple result records:

1. Select one ore more records in the table by click on the box(es) in the Location column.
2. Click the **Edit** action button.

To **edit** all result records:

1. Click on the box marked **All** in the Location column header to select all records.
2. Click the **Edit** action button.

The **Edit Result Record** screen appears with the data of the first selected record to be edited.

If multiple records are selected each record is displayed on a separate page.

EDIT RESULT RECORD: Fletcher Allen No:1712 LIS:849 Op:52

| | | | |
|--|--|--|--|
| Revised Data Operator ID: 000 Patient ID: 99 * Department: d1 * Location: l2 Result Type: Patient Sample Source: arterial/capillary * Reagent Lot: 999 Lot: Instrument Comment: 0 Level: Instrument Text 1: CHAND Text 2: Accession No: Text 3: Change Comment: | | Review Information <input type="checkbox"/> Mark As Reviewed <input type="checkbox"/> Mark As Excluded Review Comment: Predefined: Recent: Custom: User ID: Review Date/Time: Result Actions <input type="checkbox"/> Forward to LIS <input type="checkbox"/> Bypass Error Check Alarms Generated By This Record Result Edit Not Forwarded <input type="checkbox"/> Acknowledge Alarm Result Edit Not Forwarded <input type="checkbox"/> Acknowledge Alarm | |
|--|--|--|--|

<< 1 of 1 >>
 Save Cancel

| | | | |
|--------------------------|-----------------------------------|-----------------------------------|--|
| Current Data | | | |
| Operator ID: 000 | Operator Name: -- | Instrument Serial No.: 304085 | Instrument Name: -- |
| Patient ID: 99 | Patient Name: -- | Patient Date Of Birth: -- | Patient Gender: -- |
| ID Confirmation: -- | Test Date/Time: 08/24/2020 5:09PM | Accession Number: -- | LIS Transfer Status: Pending |
| Instrument Comment: 0 | Department: d1 | Location: l2 | LIS Receive Time: -- |
| Instrument Text 1: CHAND | Result Type: Patient | Sample Source: arterial/capillary | Alarms: Yes |
| Instrument Text 2: -- | Reagent Type: EC8+ | Reagent Lot: 999 | Excluded: No |
| Instrument Text 3: -- | Lot: -- | Level: -- | Reviewed: No |
| Review Comment: -- | | | |

LIS Message: --

| | | | | | |
|------------------------|-----------------|----------------|--------------|----------------|--------------|
| Analyte Results | | | | | |
| Cl 114 mmol/L | K 3.6 mmol/L | Na 153 mmol/L | BUN 56 mg/dL | pH >8.200 | GLU 53 mg/dL |
| TCO2 <> mmol/L | AnGap <> mmol/L | PCO2 <5.0 mmHg | BE <> mmol/L | HCO3 <> mmol/L | Hb <> g/dL |
| Hct <10 %PCV | | | | | |

| | | | | | |
|---------------------------|---------------|-------------------|-------------------|--------------------|--|
| Test Attributes << | | | | | |
| BatteryVoltage: 8.19V | Uses: 228 | Software: 115-A95 | ID: 00 304085 228 | MethodName: i-STAT | |
| Field1: CHAND | PanelCode: 08 | StaroutCode: 00 | FluidCode: 80 | ISTATID: 10 | |
| >CommentText: 0 | | | | | |

To view the selected result records, click on the  right arrow button to view the **next** record or click on the  left arrow button to view the **previous** record.

The screen is divided into three sections:

- **Revised Data** – upper left corner of the screen
- **Review Information, Result Action and Alarms Generated By This Record** – upper right corner
- **Current Data, Analyte Results and Test Attributes** – bottom half of the screen

Revised Data

The following entry fields can be edited under the **Revised Data** section:

| Field | | Description |
|---|------------------------------|--|
| Operator ID | | Enter the correct Operator ID number Entry field accepts up to a maximum of 30 characters. |
| Note: The updated Operator ID must be in the system or the following error will occur: Error 1150: Assigned operator missing from Operators table. | | |
| Patient ID | | Entry field accepts up to a maximum of 30 characters. Field is disabled if the record is not associated with a Patient test. |
| Department | | Click on a Department in the Tree. The read-only field for Department automatically displays the selected Department. |
| Location | | Click on a Location in the Tree. The read-only field for Location automatically displays the selected location. |
| Result type: | Patient | Enables the Patient ID, Reagent Lot and Sample Source fields for editing. |
| | Control | Enables the Lot, Level and Reagent Lot fields for editing. |
| | Electronic QC | Enables Reagent Lot. |
| | Linearity Proficiency | Enables the Lot, Level and Reagent Lot fields for editing. Enables the Lot and Reagent Lot fields for editing. |
| Sample Source | | Enabled only if "Patient" result record is selected. Select sample source from dropdown list. |
| Reagent Lot | | Select reagent lot from the dropdown list. |
| Lot | | Enabled only if Control, Proficiency or Linearity result record is selected from the dropdown list for the Result Type field . Select lot from the dropdown list. These lot numbers must be predefined in the Lots screens (control, linearity, etc.) See Chapter 8, Lots . |
| Instrument Comment | | Select an instrument comment code from the dropdown list. 0 - 8 are predefined comments for PCx only; additional comments can be added in the Instrument Comment Codes screen. |

| Field | Description |
|--|---|
| Level | <p>Enabled only if Control or Linearity result record is selected from the dropdown list for the Result Type field.</p> <p>Select Low, Normal or High for QC. Select level numbered L1 through L10 for Linearity.</p> |
| Instrument Text 1, Text 2, Text 3 | <p>Displays the instrument freeform text field 1, 2 or 3 and allows entry.</p> <p>The value of any Instrument Test Attribute with the name "Field 1," "Field 2," "Field 3," respectively, will be displayed. If the instrument type does not support this attribute, the label and entry field will not be displayed.</p> |

Note: The **Edit Result** screen will not display the Instrument Text fields when **PCx** is selected from the Instrument field, as the PCx meter does not support free text entry.

| | |
|-----------------------|---|
| Accession No. | Displays the LIS/HIS Accession Number. This field can be edited. |
| Change Comment | Enter a comment related to the reason for editing the record. This comment will be forwarded to the LIS along with the entire record. |

Review Information

The following entry fields can be edited under the **Review Information** section:

☐ Mark as Reviewed

☐ Mark as Excluded

Review Comment entry fields

| Field | Description |
|-------------------------|--|
| Predefined | <p>Select from a list of 10 predefined comments from the dropdown list.</p> <p>Predefined comments are entered in the Facility Settings screen; see Chapter 9, Section Facility Settings.</p> <p>This selection disables the Custom and Recent comment fields.</p> |
| Recent | <p>Select a predefined comment from the dropdown list. Up to 10 comments will appear in this list.</p> <p>The comments listed are the most recently entered in the Custom comment field, below.</p> <p>The most recently entered comment in the Custom comment field or predefined comment will appear as the default entry in the Recent comment field until changed.</p> <p>This selection disables the Custom and Predefined comment fields.</p> |
| Custom | <p>Enter a comment (up to 60 alphanumeric characters) in the entry field.</p> <p>This selection disables the Predefined and Recent comment fields.</p> <p>When a comment is entered it will be added to the Recent comment selection list upon saving.</p> <p>The list supports up to ten custom comments. If an eleventh comment is added, the first comment entered is removed from the selection list.</p> |
| User ID | Your user ID is automatically inserted in this field. This field cannot be modified. |
| Review Date/Time | The date and time you are editing the data is automatically inserted in this field. This field cannot be modified. |

Result Actions

The following entry fields can be edited under the **Result Actions** section:

☐ Forward to LIS ☐ Bypass Error Check

- Bypass Error Check overwrites the interface's upfront error check (if enabled). Enabling the Bypass Error Check setting will forward result record to the LIS/HIS even if the result contains an error that would otherwise fail the defined interface error check.
- Select the Forward to LIS box to resend the result record to the LIS/HIS.

- If a result record has been change, but the **Forward to LIS** box remains unchecked upon saving the record, a "Result Edit Not Forwarded" alarm will be generated.
- The "Result Edit Not Forwarded" alarm will be displayed in the **Dashboard** under the LIS/Test category. Selecting the alarm will display the **All Results** screen showing the Result Edit Alarms.

Notes:

- *If you have selected the **Forward to LIS** box, the **Change Comment** field and **Bypass Error Check** flag you have specified will be forwarded to the LIS/HIS except in QCM2.x forwarding mode.*
- *The **Bypass Error Check** box is disabled if **Forward to LIS** is not checked. The box is also disabled if the LIS forwarding status is "Pending" or "N/A."*





Alarms Generated By This Record

This section displays all Alarm Types associated with the result you are currently editing.

If the **Edit Result** screen you are viewing is displayed coming from the **All Results** screen with a specific Alarm Type filter, the corresponding Alarm Type will be checked by default.

- Unacknowledged alarms display a selectable checkbox.
- Acknowledged alarms display a grayed and checked checkbox.
- Selected alarms will be acknowledged upon saving the record.


Page Counter

The  left arrow button allows you to move to previous records; the  right arrow button allows you to move to subsequent records. The position of the record you are viewing is displayed in relation to other records found that fit the Selection Criteria, for example "2 of 32," is displayed between the  left arrow and  right arrow symbols.

- Click the **Save** button to save the edited record. The next selected record is then displayed or Click the **Cancel** button to cancel the edit action without saving.

Warning Message

If you click the **Save** button without entering a comment in the Change Comment field, the following Warning Message is displayed:



- Click **OK** if you want to go on to the next record without adding a comment, or
- Click **Cancel** if you would like to add a comment to the Change Comment field.

Current Data Table

The **Current Data** table shows the current data you are editing in the Revised Data section of the Edit Result Record screen.

The current data is in read-only fields.

When an edit is made to the Patient ID field in the Revised Data section, the Patient Data Source field will be set to "QCM." The Patient ID associated with the test result will be updated.

The data displayed in the Patient ID, Patient Name, Patient Date of Birth, Patient Gender, and Patient Data Source fields will be read from Admission, Discharge, and Transfer data associated with the test result.

Field text related to alarms will be highlighted with white letters on red background for unacknowledged alarms and black text on yellow background for acknowledged alarms.

Analyte Results


Analyte result values are displayed below the **Current Data** table.

Only Analytes for the Reagent Type will be displayed.

Analyte result values that caused an **Out-of-Range** alarm will be displayed highlighted with white text on red background in the Analyte value field.

Note: All analyte values, except for Glucose (GLU), will be displayed only in the **Edit Results, View Result Details** and **Patient Results** screens.

Test Attributes

The read-only Instrument **Test Attributes** table is only displayed by clicking the  button. The Test Attributes table appears only for some instruments, including i-STAT and i-STAT Strip. Various lists of Test Attributes will be shown, depending on the instrument type and the particular test type.

When the attributes are displayed, the button reads: 

Clicking the button will hide the **Test Attributes** table and re-display the  button.

Attributes that may be displayed in this table include:

| | | |
|--------------------|----------------------------|-------|
| SimulatorID | Temperature | POT0P |
| Pressure | CommentText | POT0N |
| ProbeDelta | Run | POT0Z |
| CPBApplied | MethodName | POT1P |
| FIO2 | ID | POT1N |
| Low Range | ActLo | POT1Z |
| High Range | ActHi | POT2P |
| Field1 | RefLo | POT2N |
| Field2 | RefHi | POT2Z |
| Field3 | AMP0Z | POT3P |
| PanelCode | AMP1P | POT3N |
| StaroutCode | AMP1N | POT3Z |
| FluidCode | AMP1Z | POT4P |
| ISTATID | AMP2P | POT4N |
| PreferencesName | AMP2N | POT4Z |
| PreferencesVersion | AMP2Z | POT5P |
| Uses | OriginalTestType | POT5N |
| Software | RHI | POT5Z |
| PanelCode | RLO | POT6P |
| BatteryVoltage | Sample Source | POT6N |
| | CplItem2 | POT6Z |
| | CplItem3 | AMP0P |
| | CplItem4 | AMP0N |
| | CpSampleType | |
| | CpPtTemp | |
| | CplItem6 | |
| | CpCPB | |
| | CplItem9 through CplItem54 | |

View Results

Clicking this button displays the **View Results Record** screen, which shows detailed current data for the record(s) chosen. Result records cannot be modified from View Results screen.

The following aspects of the screen are the same as Current Data table in the **Edit Result** screen:

- **Forward** and **Backward** page navigation buttons
- **Cancel** button
- **Current Data** table
- **Analyte Results** table
- **Test Attributes** table

VIEW RESULT RECORD: F3 No: 175 LIS: 123 Op: 9

<< 1 of 1 >>
Cancel

Current Data

| | | | |
|-------------------------|-----------------------------------|-----------------------------------|--------------------------|
| Operator ID: 223 | Operator Name: -- | Instrument Serial No.: M06210023 | Instrument Name: -- |
| Patient ID: 10000 | Patient Name: -- | Patient Date Of Birth: -- | Patient Gender: -- |
| Patient Data Source: -- | Test Date/Time: 01/09/2007 9:56AM | Accession Number: -- | LIS Transfer Status: N/A |
| Instrument Comment: -- | Department: D3 | Location: L3 | LIS Receive Time: -- |
| Instrument Text 1: -- | Result Type: Patient | Sample Source: arterial/capillary | Alarms: Yes |
| Instrument Text 2: -- | Reagent Type: Glucose Test Strip | Reagent Lot: 178365 | Excluded: No |
| Instrument Text 3: -- | Lot: -- | Level: -- | Reviewed: No |
| Review Comment: -- | | | |

LIS Message: --

Analyte Results

GLU 74 mg/dL


Re-Route Result Data

Clicking this button allows you to acknowledge and forward selected results to the LIS/HIS.

In order to activate the Re-Route action you must choose an alarm type from the **Result-Based Alarm** dropdown list in the Selection Criteria.

A number of **dialogue boxes** are triggered when using the **Re-Route** action:

If you choose "All" from the Result-Based Alarm list, the following message will be displayed:



- Click **OK** to acknowledge all the alarms for the selected results. Only result-based alarms are acknowledged. Click **Cancel** to return to the All Results screen.

Once you have chosen an alarm type from the Result-based Alarm dropdown list, and clicked on the **Re-Route** button, the following message will be displayed:



- Click **OK** to continue or **Cancel** to return to the All Results screen.


Similar checks as in the **Edit Result** screen will be performed before forwarding a result to the LIS:

- The facility must be configured to forward results to the LIS.
- Tests with an LIS transfer status of "Pending" cannot be forwarded to the LIS.
- A result type must be configured for the facility to be forwarded to the LIS (Patient and/or QC tests only).
- The test date/time must be greater than the LIS Start Forwarding date for the facility.

If any of these (above) conditions fail, a dialogue box will inform you of the error and ask whether to continue:


- If you click OK and continue, only the results that pass the check will be forwarded, and all the alarms will be acknowledged.

If you forward one or more results to the LIS that have already been successfully forwarded, an additional dialogue box is displayed to warn the user.



- Click **OK** to continue or **Cancel** to return to the All Results screen.

If you have not selected an alarm type from the Result-Based Alarm field, the following warning message is displayed:



- Click **OK** and make a selection from the Result-Based Alarm dropdown list.

Note: The Re-Route button will be displayed only for users with Edit Results and Forward to LIS permissions.

Acknowledge Results

Clicking this button allows you to acknowledge alarms connected with the selected results of the displayed Alarms Selection Criteria.

In order to activate the **Acknowledge** action you must choose an alarm type from the **Result-Based Alarm** dropdown list in the Selection Criteria.


A number of dialogue boxes are triggered when using the **Acknowledge** action:

If you choose "All" from the Result-Based Alarm list, the following message will be displayed:




- Click **OK** to acknowledge all the alarms for the selected results (including result-based and non-result-based alarm categories).
- Click **Cancel** to return to the All Results screen.

Once you have chosen an alarm type from the Result-based Alarm dropdown list and clicked on the **Acknowledge** button, the following message will be displayed:



- Click **OK** to proceed.
- Click **Cancel** to return to the All Results screen.

If **(No Selection)** is selected from the Result-Based Alarm list, the following warning message is displayed:



- Click **OK** and select an alarm type.

When you click the **Acknowledge** button, any alarms of the selected results that are of the Alarm Type selected in the Alarms search criteria will be acknowledged.

- For example, if the Invalid Operator alarm is selected, three results are selected and you click on the **Acknowledge** button, only the three Invalid Operator Alarms for the three results will be acknowledged. Any additional alarms that any of the three results may have will not be acknowledged.

After acknowledging results, the screen will be refreshed to reflect the acknowledged alarms by changing the field text in affected columns from white text on red background to black text on yellow background.

Note:

- The **Acknowledge** button will be displayed only for users with Edit Results permissions.
- Only previously unacknowledged alarms and alarms within the permission scope of the user will be acknowledged.


E-Mail Results

To **e-mail** a single result record or multiple result records:

1. Select one or more records in the table by click on the box(es) in the Location column.
2. Click the **E-Mail** action button.

To **e-mail** all result records:

1. Click on the box marked **All** in the Location column header.
2. Click the **E-Mail** action button.
 - The default e-mail opens a new mail message and includes the menu option title in the subject field.
 - The body of the message includes a hyperlink to an .htm file that contains the selected records.
 - The e-mailed result table will include the result details for each result in a table format.
 - The e-mailed screen will be the **View Result Record** screen.



Patient Results Screen

Review>>Patient Results

The **Patient Results** screen allows you to view and e-mail patient test results and to view test panel information for multi analyte instruments.

To **view** patient results records:

1. Move the mouse cursor over **Review** to display the menu options.
2. Click on the **Patient Results** option.

| Location | Test Date/Time | Instrument Name | Instrument Serial No. | Reviewed | Review Comment | Operator Name | Operator ID | Instrument Comment | Patient |
|----------|-------------------|-----------------|-----------------------|----------|----------------|---------------|-------------|--------------------|---------|
| L1 | 01/03/2006 1:54PM | -- | 302026 | No | -- | -- | -- | -- | 1 |
| L1 | 01/03/2006 2:00PM | -- | 302026 | No | -- | -- | -- | -- | 2 |
| L1 | 01/03/2006 2:02PM | -- | 301701 | No | -- | -- | -- | -- | 3 |
| L1 | 01/03/2006 2:05PM | -- | 302026 | No | -- | -- | -- | -- | 4 |
| L1 | 01/03/2006 2:06PM | -- | 301701 | No | -- | -- | -- | -- | 5 |
| L1 | 01/03/2006 2:10PM | -- | 301701 | No | -- | -- | -- | -- | 6 |
| L1 | 01/03/2006 2:16PM | -- | 301701 | No | -- | -- | -- | -- | 7 |
| L1 | 01/03/2006 2:19PM | -- | 301701 | No | -- | -- | -- | -- | 9 |
| L1 | 01/03/2006 2:19PM | -- | 302026 | No | -- | -- | -- | -- | 8 |
| L1 | 01/03/2006 2:30PM | -- | 301701 | No | -- | -- | -- | -- | 11 |
| L1 | 01/03/2006 2:30PM | -- | 302026 | No | -- | -- | -- | -- | 10 |
| L1 | 01/03/2006 2:33PM | -- | 302026 | No | -- | -- | -- | -- | 12 |
| L1 | 01/03/2006 2:34PM | -- | 301701 | No | -- | -- | -- | -- | 13 |
| L1 | 01/03/2006 2:38PM | -- | 301701 | No | -- | -- | -- | -- | 15 |
| L1 | 01/03/2006 2:38PM | -- | 302026 | No | -- | -- | -- | -- | 14 |
| L1 | 01/03/2006 2:45PM | -- | 299999 | No | -- | -- | 147 | -- | 55 |
| L1 | 01/06/2005 6:52PM | -- | 299999 | No | -- | -- | -- | -- | unkno |
| L1 | 01/06/2005 7:00PM | -- | 299999 | No | -- | -- | -- | -- | unkno |

- The **From** and **To** date fields are filled in by default with today's date. To display records you must first set the **From** Date.
- The **Patient Results** screen displays the Table based on the location highlighted in the Location Tree and the system default Search Criteria.
- To display patient results records from a *different* location using the same search criteria, click on a **Facility**, **Department** or **Location** in the Location Tree.
- Use the Search Criteria options to display specific result records.

Search Criteria

To **find** a specific patient results record:

1. Enter a **Patient ID** in the entry field.
2. Click the **Display Records** button.
 - Results that match or are similar to the entered Patient ID will be displayed.

To **filter** the patient records displayed:

1. Select from the **Search Criteria** dropdown list available.
2. Click the **Display Records** button.

Available filter options:

| Filter | Options | Description |
|-------------------|---------------------|--|
| Panel | No Selection | |
| | Selected | System displays results based on the panel created from the Administrative >> Panels screen. |
| Test | All | |
| | GLU | System only displays glucose test results. |
| | Selected | System displays results with selected Test Type. |
| Alarm | No Selection | System displays results with and without associated alarms. |
| | All | System only displays results with associated alarms. |
| | Selected | System only displays results with the selected alarm type. |
| Date Range | | System default displays today's date. Enter or select a From and a To date to display test results records from tests that were performed between the selected dates. |
| | From | Click on the calendar icon to display the pop-up calendar and click on a date. The From date field automatically displays the selected date. |
| | To | Click on the calendar icon to display the pop-up calendar and click on a date – required . The To date field automatically displays the selected date. |

Click the **Display Records** button to display the results based on the selected Search Criteria options.

E-Mail Patient Result Records

To **e-mail** a single patient result record or multiple patient result records:


1. Select one or more records in the table by click on the box(es) in the Location column.
2. Click the **E-Mail** action button.

To **e-mail** all patient result records:

1. Click on the box marked **All** in the Location column header.
2. Click the **E-Mail** action button.

The default e-mail client opens a new mail message and includes the menu option title in the subject field.

The body of the message includes a hyperlink to an .htm file that contains the selected records. Linked files will expire after 30 days.



Patient Lookup

Review>> Patient Lookup

The **Patient Lookup** screen allows you to retrieve current patient data for any Facility in your organization. By default, existing users are not granted access to the Patient Lookup screen. To obtain access to this screen refer to [Menu Permissions](#) in Chapter 9.

| Current Data | | 1 of 7001 | |
|-----------------------|---------------------|------------------------|------------|
| Timestamp | 02/07/2007 09:41 AM | Patient First Name | Frank10000 |
| Patient ID | 10000 | Patient Last Name | Zappa10000 |
| Patient Date Of Birth | 06/17/1953 | Patient Middle Initial | A |
| Admission Date | 01/02/2006 07:01 AM | Patient Gender | M |

To **view** this screen:

1. Move the mouse cursor over **Review** in the Menu Bar to display the menu options.
2. Click on the **Patient Lookup** option.
3. Fill in the desired information. You may look up specific records by using any of the following search criteria:
 - Patient ID
 - Patient Date of Birth
 - Patient First Name
 - Patient Last Name
 - Admission Date
4. Click **Display Records**. The **Current Data** table displays the data for the patient(s) you have selected. Patient Records will be retrieved for the selected Facility in the Location Tree. The **Timestamp** field displays the time the data was last updated in the LIS/HIS.

Note: The Organization, Department, and Location nodes in the Tree are disabled in this screen.

The « left arrow button allows you to move to previous records; the » right arrow button allows you to move to subsequent records. The position of the record you are viewing is displayed in relation to other records found that fit the Selection Criteria, for example "1 of 9" is displayed between the « left arrow and » right arrow symbols. To limit the number of patient records recalled, use specific criteria to narrow the search.

5. Reports

Reports Overview

The **Reports** menu option provides access to formatted reports suitable for archiving on paper or online. There are three screens available. The **View** screen allows you to edit, view, save or e-mail a report or suite of reports. The **Suites** screen allows you to add, edit or delete suites of reports that can then be managed as a group. The **Import** screen allows you to add customized reports to the system. Each screen is available only if your account was granted this permission by the Administrator or POCC.

There are 23 standard report types, grouped into 2 suites:

QCM3 Default Reports

1. Audit Trail Log
2. Comment Code Log
3. Control Data Statistics by Instrument
4. Control Data Statistics by Reagent Lot
5. Data Listing
6. Glucose Meter Nurse Manager Report
7. Glucose Meter QI Representative Report

Note: *Glucose meter reports are specific to Abbott glucose monitors only and do not report data for other instrument types.*

8. Instrument Configuration
9. Instrument Error Log
10. Instrument Workload Statistics
11. Levey-Jennings Graph
12. Levey-Jennings Graph by Reagent Lot
13. Linearity Graph
14. Operator Certification Expiration
15. Operator Competency
16. Patient Outlier Log
17. Proficiency Data Listing
18. QC Outlier Log

i-STAT Default Reports (Reports are specific i-STAT instrument types only)

19. i-STAT Data Listing
20. i-STAT Instrument Event
21. Quality Check Codes
22. Quality Check Codes by Operator
23. Reagent Usage

Each report in the Default Suites can be edited and saved into customized Report Suites to create batched reports that meet the needs of your facility.

VIEW REPORTS:

Selection Criteria

Suite Name

18 Records Found

| Edit | View | Print | E-Mail |
|------------------------------|--|-----------------------|------------------------|
| Report Name | | | |
| <input type="checkbox"/> All | | | |
| <input type="checkbox"/> | Audit Trail Log | | |
| <input type="checkbox"/> | Comment Code Log | | |
| <input type="checkbox"/> | Control Data Statistics by Instrument | | |
| <input type="checkbox"/> | Control Data Statistics by Reagent Lot | | |
| <input type="checkbox"/> | Data Listing | | |
| <input type="checkbox"/> | Glucose Meter Nurse Manager Report | | |
| <input type="checkbox"/> | Glucose Meter QI Representative Report | | |
| <input type="checkbox"/> | Instrument Configuration | | |
| <input type="checkbox"/> | Instrument Error Log | | |
| <input type="checkbox"/> | Instrument Workload Statistics | | |
| <input type="checkbox"/> | Levey-Jennings Graph | | |
| <input type="checkbox"/> | Levey-Jennings Graph by Reagent Lot | | |
| <input type="checkbox"/> | Linearity Graph | | |
| <input type="checkbox"/> | Operator Certification Expiration | | |
| <input type="checkbox"/> | Operator Competency | | |
| <input type="checkbox"/> | Patient Outlier Log | | |
| <input type="checkbox"/> | Proficiency Data Listing | | |
| <input type="checkbox"/> | QC Outlier Log | | |
| Edit | View | Print | E-Mail |

18 Records Found

View Reports Screen

View>>Reports

The **View Reports** screen allows you to view available report suites and to select reports for editing, viewing, printing and e-mailing.

To **view** a report suite:

1. Move the mouse cursor over **Reports** in the Menu Bar to display the menu options.
2. Click on the **View** option.

Your location in the Tree is irrelevant at this time as you will select a Tree location once a report has been selected from a Suite.
3. Select a suite name from the dropdown list in the **Suite Name** field of the **Selection Criteria**.

The default suites are QCM3 default and i-STAT default.

The reports associated with the **Suite Name** you have chosen in the Selection Criteria are displayed in the Table on the bottom half of the screen.

There are five available action buttons in the **Action Bar** above the list of reports: **Edit**, **View**, **Delete** (not available for default reports), **Print** and **E-mail**.

Edit a Report

To **edit** a single report or multiple reports in a suite:

1. Click the box(es) to the left of the Report Name.
2. Click the **Edit** action button.

To **edit** all reports in a suite:

1. Click the box marked **All** in the Report Name column header.
2. Click the **Edit** action button.

EDIT REPORT:

Audit Trail Log

Location Choice All Facilities » All Departments » All Locations

Instrument PCx

Comment Code

Patient ID ☐ Suppress

Patient Last Name

User ID ...

Operator ID ☐ Suppress

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

View Print E-Mail Save Cancel

The **Edit Report** screen allows you to customize and save reports that include:

- The Organization, Facility, Department or Location specific data.
- The selected filter options, date range options and print options.
- The selected database (Live or Archive) as the source for the report.
- The sort option for the report.

The next time the saved report is viewed or printed, the application automatically updates the information with the current data from the database eliminating the need to re-create the report.

If a user-defined date range is selected for the saved report, you can edit the date range before viewing, printing or e-mailing the report.

On reports that include the Operator ID / Operator Name and/or the Patient ID / Patient Name, , the application provides the options to include or exclude the ID and Name from the report. This option is available for security reasons.

See Section [Default Reports](#) below for more information on selecting filter options in the **Edit Report** screen.

View a report

To **view** a single report or multiple reports in a suite:

1. Click the box(es) to the left of the Report Name in the **View Reports** screen and click the **View** action button, or
2. Click the **View** action button on the **Edit Report** screen to view that specific report.

Print a Report

To **print** a single report or multiple reports in a suite:

1. Click the box(es) to the left of the Report Name in the **View Reports** screen and click the **Print** action button, or
2. Click the **Print** action button on the **Edit Report** screen to print that specific report.

To **print** all reports in a suite:

1. Click the box marked **All** in the Report Name column header.
2. Click the **Print** action button.


E-Mail a Report

To **e-mail** a link to a single report or multiple reports:

1. Click the box(es) to the left of the Report Name in the **View Reports** screen and click the **E-Mail** action button, or
2. Click the **E-Mail** action button on the **Edit Report** screen to e-mail a link to the specific report.

To **e-mail** links to all reports in a suite:

1. Click the box marked **All** in the Report Name column header.
2. Click the **E-Mail** action button.
 - The default e-mail client opens a new mail message and includes the menu option title in the subject field. The body of the message includes a hyperlink to a PDF file that contains the selected report.
 - Linked report files will be deleted after 30 days or if a maximum directory size specified on the registry has been reached. The default directory size is 200MB.



Save a Report

Click the **Save action** button to save the selected filter choices and sort options for the report.

Note: A default report cannot be edited and saved to a default suite. An edited report with the same name as a default report can be saved under a new customized suite name.

The **Save Report** screen displays the current Report Name in the **New Name** entry field.

- To save the report with a different name, enter a Report Name in the **New Name** entry field.
- To save the report to an existing customized suite, click on the **Existing** button and select a Suite Name from the dropdown options.
- To save the report to a new customized suite, click the **New** button and enter a New Suite Name.

The screenshot shows the PrecisionWeb application interface. On the left is a navigation tree with 'Medical Center Inc' expanded, showing sub-items: Eastern Medical Center, Pediatrics, Southern Medical Center, and Western Medical Center. The main content area is titled 'SAVE REPORT: All Facilities » All Departments » All Locations'. It contains two main sections: 'Report Name' with a text input field containing '* New Name QC Outlier Log', and 'Suite Name' with two radio buttons: 'Existing' (which is selected) and 'New'. Below the 'Existing' radio button is a dropdown menu showing 'test'. At the bottom of the form are 'Save' and 'Cancel' buttons.

The next time the saved report is viewed or printed, the application automatically updates the information with the current data from the database eliminating the need to re-create the report.

* The New Name is a **required** entry field.

Note: If a user-defined date range is selected for the saved report, edit the report to update the date range before viewing or printing the report.

Report Import Screen

Reports>>Import

The Report Import Screen allows you to import reports from other files on your system.

To **import** a report:

1. In the Import Details box enter the **Report Name**.
2. Select a **Suite Name** from the dropdown list of available suites.

Note: Default suites are not available from the suite name dropdown list. To import a report, you must select from an existing suite or create a new suite prior to attempting the import.

3. Select the **File** location of the report file by clicking on the **Browse** button.
 - All three fields are **required** fields.
4. Click the **Import** button.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

REPORT IMPORT: No: 12 LIS: 1 Op: 2

Show Tree

Input Details

* Report Name

* Suite Name

* File Browse...

Import

The imported report will now be available in the **View Report** screen under the suite name you selected during the importing process.

Reports

The filter options for all 23 default reports are listed in the following pages.

See Section [Edit a Report](#) above for additional information on editing the reports.

To **edit** a single report or multiple reports in a suite:

1. Click the box(es) to the left of the Report Name.
2. Click the **Edit** action button.

To **edit** all reports in a suite:

1. Click the box marked **All** in the Report Name column header.
2. Click the **Edit** action button.

Once in the **Edit Report** screen:

1. Click on a **Facility, Department** or **Location** in the Tree to customize the report for specific location.
2. Select the filter options available on the **Edit Report** screen for the report you have chosen to edit.
3. After customizing the report, click on the **View, Print, E-Mail, Save** or **Cancel** action buttons. See sections [View a report](#), [Print a Report](#), [E-Mail a Report](#) and [Save a Report](#) above for more information on these actions.

Notes on common filter options

All reports except for **Instrument Configuration** and **Operator Certification Expiration** have the same **Test Date Range, Print Options** and **Additional Criteria** filter options as listed below:

Note: *Instrument Configuration has **no** Test Date Range and Operator Certification Expiration has **Expiration Date Range** instead of Test Date Range, which will be explained below in Section [Operator Certification Expiration Report](#).*

Test Date Range

Range button

Click on the **Range** button to display the calendar icon in the **From** and **To** entry fields (this option disables the **Period** button).

- The **From** and **To** date fields automatically display today's date unless changed.
- To change the date range:

- Click on the **31** calendar icon next to the **From** field to display the pop-up calendar and select an appropriate date.
- Click on the **31** calendar icon next to the **To** field to display the pop-up calendar and select an appropriate date.



Period button

- Click on the **Period** button and select a period from the dropdown list.

This option disables the **Range** button

EDIT REPORT:

Audit Trail Log

Location Choice All Facilities » All Departments » All Locations

Instrument

Comment Code

User ID ...

Operator ID ☐ Suppress

Patient ID ☐ Suppress

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☒ Suppress Cover Page

☒ Suppress Page Breaks

Additional Criteria

Data Source

Sort By

« 1 of 1 »

Print Options

Click on the checkboxes to exclude the cover page and/or the page breaks in the printed report.

☐ **Suppress Cover Page**


☐ **Suppress Page Breaks**

Additional Criteria

Click the **More Options >>** button to display two **Additional Criteria** fields: **Data Source** and **Sort By**.

Data Source


Select the Live or Archive Database



The screenshot shows the 'Additional Criteria' dialog box. It has two dropdown menus: 'Data Source' and 'Sort By'. The 'Data Source' dropdown is set to 'Live Database'. The 'Sort By' dropdown is also set to 'Live Database'. To the right of the 'Sort By' dropdown is a button labeled '<< Less Options'. At the bottom of the dialog box, there are five buttons: 'View', 'Print', 'E-Mail', 'Save', and 'Cancel'. Above the 'View' button, there are navigation buttons: '<<', '1 of 1', and '>>'.

Sort By

Select option from the dropdown list



The screenshot shows the 'Additional Criteria' dialog box with the 'Sort By' dropdown menu open. The 'Data Source' dropdown is still set to 'Live Database'. The 'Sort By' dropdown menu is open, showing a list of options: 'Test Date/Time', 'Operator ID', 'Patient ID', 'Test Type', 'Reagent Lot', 'Level', and 'Comment Code'. The 'Test Date/Time' option is highlighted. To the right of the dropdown menu is a button labeled '<< Less Options'. At the bottom of the dialog box, there are five buttons: 'View', 'Print', 'E-Mail', 'Save', and 'Cancel'. Above the 'View' button, there are navigation buttons: '<<', '1 of 1', and '>>'.

Click the **<< Less Options** button to remove the Additional Criteria fields.

Audit Trail Log

The Audit Trail Log provides a tabular display of all changes made to test result record(s).

EDIT REPORT:

Audit Trail Log

Location Choice All Facilities » All Departments » All Locations

Instrument

PCx

Comment Code

Patient ID

☐ Suppress

Patient Last Name

User ID

...

Operator ID

☐ Suppress

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

View

Print


E-Mail

Save

Cancel

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Description |
|--------------|---|
| Instrument | Select an instrument type from the dropdown list. |
| Comment Code | Enter a comment code (only one comment code can be entered per report). |
| Patient ID | Enter a Patient ID. |
| Patient Name | Enter a Patient Name, filtered by last name. |
| User ID | Click on the  button to display the pop-up browser window with the User ID Lookup list. Select a user and click on the OK button in the browser window; click the Cancel button to return to Edit Report screen. You can also directly enter a User ID into the field. |
| Operator ID | Click the Suppress checkbox to exclude the Operator ID from the printed report. |
| Patient ID | Click the Suppress checkbox to exclude the Patient ID from the printed report. |

Comment Code Log

The Comment Code Log provides a tabular display of all test results for an instrument type tagged with a comment code.

EDIT REPORT:
Comment Code Log

Location Choice All Facilities » All Departments » All Locations

Instrument ...

Comment Code

Operator ID ... ☐ Suppress

Operator Name ☐ Suppress

Result Type All

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks



More Options >>

<< 2 of 18 >>

View Print E-Mail Save Cancel

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Description |
|---------------------|--|
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list. Select an instrument type and click on the OK button in the browser window; click the Cancel button to return to Edit Report screen. |
| Comment Code | Enter a comment code (only one comment code can be entered per report). |
| Operator ID | Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |

Click the **Suppress** checkbox to exclude the Operator ID from the printed report.

| Options | Description |
|----------------------|---|
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Result type | System default displays All result types. Select a single result type from the dropdown list to customize the report. |

Control Data Statistics by Instrument

The Control Data Statistics by Instrument report displays a monthly summary of control test results, including results above or below the range, mean, SD and %CV, by instrument.

EDIT REPORT:

Control Data Statistics by Instrument

Location Choice All Facilities » All Departments » All Locations

Instrument

Instrument Serial No. ...

Reagent Lot ...

Control Lot ...

Test

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks



[More Options >>](#)


<< 1 of 1 >>

[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Description |
|------------------------------|--|
| Instrument | Select an instrument type from the dropdown list. |
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |

| Options | Description |
|--------------------|---|
| Control Lot | Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Test | <p>System default displays GLU test.</p> <p>Select a test type from the dropdown list to customize the report.</p> |

Control Data Statistics by Reagent Lot

The Control Data Statistics by Reagent Lot report displays a monthly summary of control test results, including results above or below the range, mean, SD and %CV, by Reagent Lot.

EDIT REPORT:

Control Data Statistics by Reagent Lot

Location Choice All Facilities » All Departments » All Locations

Reagent Lot ...

Control Lot ...

Test ▼

Test Date Range

Range ☐ From To

Period ☒ ▼

Print Options

☐ Suppress Cover Page



☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Description |
|--------------------|---|
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Control Lot | Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Test | System default displays GLU test. Select a test type from the dropdown list to customize the report. |

Data Listing

The Data Listing report provides a tabular display of information associated to individual test results for one or all test result types (i.e. patient, control, etc).

EDIT REPORT:

Data Listing

Location Choice All Facilities » All Departments » All Locations

Instrument

Instrument Serial No. ...

Operator ID ... ☐ Suppress

Operator Name ☐ Suppress

Patient ID ☐ Suppress

Patient Last Name

Reagent Lot ...

Control Lot ...

Result Type

Test

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page


☐ Suppress Page Breaks




More Options >>

<< 1 of 1 >>

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Description |
|------------------------------|--|
| Instrument | Select an instrument type from the dropdown list. |
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |

| Options | Description |
|----------------------|---|
| Operator ID | <p>Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen.</p> <p>Click the Suppress checkbox to exclude the Operator ID from the printed report.</p> |
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Patient ID | <p>Enter a Patient ID.</p> <p>Click the Suppress checkbox to exclude the Patient ID from the printed report.</p> |
| Patient Name | Enter a Patient Name, filtered by last name. |
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Control Lot | Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Result type | <p>System default displays All result types.</p> <p>Select a single result type from the dropdown list to customize the report.</p> |
| Test | <p>System default displays GLU test.</p> <p>Select a test type from the dropdown list to customize the report.</p> |

Glucose Meter Nurse Manager Report

The Glucose Meter Nurse Manager report displays the total number of Patient and QC tests performed, the percentage of in-range and out-of-range tests performed, the total number of outliers without a comment and total number of instrument alarms generated by an Operator. There is also space document action and follow-up.

The report also indicates if an Operator is compliant. An Operator is considered compliant if the Operator has performed a successful test. A test is considered successful if the test does not have any of the following alarms associated with it: Out of Range QC, Out of Range-LIS, Invalid Operator-LIS, Invalid Patient-LIS, Repeat Test-LIS, Invalid Operator-Inst, Invalid QC-Inst and Invalid Reagent-Inst.

If any of the alarms listed are associated with a test completed by the Operator, "1" will appear under the *# Instrument Alarms* column for every instance that the alarm(s) is generated for a test performed by the Operator (i.e. two Out-of-Range QC and three Out-of-Range Result alarms would yield a value of 5 in the *# Instrument Alarms* column). If an Operator test generates an alarm, but the alarm is not among those listed above there is no impact on the "compliance" status.

The information is presented as follows:

Operator Comment Documentation of Out of Range QC

- *% of QC in Range*
- *# QC Outliers*
- *# QC Outliers without Comments*

Operator Comment Documentation of Alert Values

- *# Instrument Alarms*
- *# Instrument Alarms without Comments*

Compliant Operator

- *Yes or No*

There is also space document action and follow-up.

This report is specific to Abbott glucose monitors only and will not report data for other instrument types.

EDIT REPORT:

Glucose Meter Nurse Manager Report

Location Choice All Facilities » All Departments » All Locations

Operator ID ☐ Suppress

Operator Name ☐ Suppress

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks

[More Options >>](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|----------------------|--|
| Operator ID | Click the Suppress checkbox to exclude the Operator ID from the printed report. |
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |

Glucose Meter QI Representative Report

The Glucose Meter QI Representative Report displays location specific quality indicator data. The report includes six indicator categories: Meter Usage, QC Within Range, Documentation of Out of Range QC, Proficiency Testing, Documentation of Alert Values and Compliance.

The report displays the total number of Patient, QC tests performed, the percentage of in-range and out-of-range tests performed, the total number of outliers without a comment, total number of Proficiency tests performed and the total number of instrument alarms generated by an Operator. There is also space document action and follow-up.

The report also indicates if an Operator is compliant. An Operator is considered compliant if the Operator has performed a successful test. A test is considered successful if the test does not have any of the following alarms associated with it: Out of Range QC, Out of Range-LIS, Invalid Operator-LIS, Invalid Patient-LIS, Repeat Test-LIS, Invalid Operator-Inst, Invalid QC-Inst and Invalid Reagent-Inst.

If any of the alarms listed are associated with a test completed by the Operator a 1 will appear under the *# Instrument Alarms* column for every instance that the alarm(s) is generated for a test performed by the Operator (i.e. two Out-of-Range QC and three Out-of-Range-Result alarms would yield a value of 5 in the *# Instrument Alarms* column). If an Operator test generates an alarm, but the alarm is not among those listed above there is no impact on the "compliance" status.

The information is presented as follows:

Meter Usage (Indicator 1)

- Total # Patient Sample
- Total QC performed

QC Within Range (Indicator 2)

- % of QC in Range

Documentation of Out of Range QC (Indicator 3)

- # QC Outliers
- # QC Outliers without Comments

Proficiency Testing (Indicator 4)

- Proficiency testing performed during this period

Documentation of Alert Values (Indicator 5)

- # Instrument Alarms
- # Instrument Alarms without Comments

Compliant Operator (Indicator 6)

- Yes or No

This report is specific to Abbott glucose monitors only and will not report data for other instrument types.

EDIT REPORT:
Glucose Meter QI Representative Report

Location Choice All Facilities » All Departments » All Locations

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

View Print E-Mail Save Cancel

Select a Facility Department or Location in the Tree.

Select the filter options for the report.

Instrument Configuration Report

The Instrument Configuration report displays the configuration settings of the specified Abbott glucose monitor by location.

The instrument configuration settings must be saved (new or default settings) prior to running this report, otherwise the report will be blank.

EDIT REPORT:
Instrument Configuration

Location Choice All Facilities » All Departments » All Locations

Instrument

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

View Print E-Mail Save Cancel

Select a Facility Department or Location in the Tree.

Select the filter options for the report:


Instrument Select an instrument type from the dropdown list.

Instrument Error Log

The Instrument Error Log provides a tabular display of error codes uploaded from Abbott glucose monitors only.

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|------------------------------|--|
| Instrument | Select an instrument type from the dropdown list. |
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial Number, multiple Instrument Serial Numbers or all Instrument Serial Numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Error Code | Enter an error code (only one error code can be entered per report). |

Instrument Workload Statistics

The Instrument Workload Statistic report displays a tabular summary of workload statistics per instrument including the number of patient, control, proficiency, linearity and electronic QC (when applicable) tests.

EDIT REPORT:

Instrument Workload Statistics

Location Choice All Facilities » All Departments » All Locations

Instrument ...

Test Date Range

Range ☐ From To

Period ☒ Today

Print Options

☐ Suppress Cover Page


☐ Suppress Page Breaks

More Options >>

<< 1 of 1 >>

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|-------------------|---|
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list, select an Instrument, multiple Instruments or all Instruments, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |



Levey-Jennings Graph



The Levey-Jennings Graph provides a graphical and tabular display per instrument of control data trends and statistics by reagent lot.

EDIT REPORT:
Levey-Jennings Graph
 Location Choice All Facilities » All Departments » All Locations
 Instrument ...
 Instrument Serial No. ...
 Reagent Lot ...
 Control Lot ...
 Test ▼
Test Date Range
 Range ☐ From To
 Period ☒ ▼
Print Options
☐ Suppress Cover Page
☐ Suppress Page Breaks

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|------------------------------|--|
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list, select an Instrument, multiple Instruments or all Instruments, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |

| Options | Descriptions |
|--------------------|---|
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Control Lot | Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Test | System default displays GLU test. Select a test type from the dropdown list to customize the report. |

Levey-Jennings Graph by Reagent Lot

The Levey-Jennings Graph by Reagent Lot provides a graphical and tabular display of control data trends and statistics by reagent lot only.

EDIT REPORT:

Levey-Jennings Graph by Reagent Lot

Location Choice All Facilities » All Departments » All Locations

Reagent Lot ...

Control Lot ...

Test

Data Display

Test Date Range

Range ☐ From To

Period ☒

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks



[More Options >>](#)

<< 1 of 1 >>

[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|---------------------|---|
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Control Lot | Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Test | System default displays GLU test. Select a test type from the dropdown list to customize the report. |
| Data Display | System default displays All Data . Select data display type from the dropdown list to customize the report. |

Linearity Graph

The Linearity Graph provides a graphical and tabular display per instrument of calibration verification and linearity test results across the reagent assay range.

EDIT REPORT:
Linearity Graph

Location Choice All Facilities » All Departments » All Locations

Instrument

Instrument Serial No. ...

Linearity Lot ...

Test

Test Date Range

Range ☐ From To

Period ☒ Previous Month

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks



[More Options >>](#)

<< 1 of 1 >>

[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:


| Options | Descriptions |
|------------------------------|--|
| Instrument | Select an instrument type from the dropdown list. |
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Linearity Lot | Click on the  button to display the pop-up browser window with the Linearity Lot Lookup list, select a Linearity Lot, multiple Linearity Lots or All Linearity Lots, and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |
| Test | System default displays GLU test. Select a test type from the dropdown list to customize the report. |

Operator Certification Expiration Report

The Operator Certification Expiration report provides a tabular display of operators whose certification will expire during a specified date range.

Select a Facility Department or Location in the Tree.



Select the filter options for the report:

| Options | Descriptions |
|----------------------|---|
| Operator ID | <p>Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen.</p> <p>Click the Suppress checkbox to exclude the Operator ID from the printed report.</p> |
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Instrument | Select an instrument type from the dropdown list. |

Expiration Date Range**Range** button

- Click on the **Range** button to display the calendar icon in the **From** and **To** entry fields (this option disables the **Period** button).
- The **From** and **To** date fields automatically display today's date unless changed.

To **change** the date range:

- Click on the  calendar icon next to the From field to display the pop-up calendar and select an appropriate date.
- Click on the  calendar icon next to the **To** field to display the pop-up calendar and select an appropriate date.

Period button


- Click on the **Period** button and select a period from the dropdown list. This option disables the **Range** button.

Operator Competency Report

The Operator Competency report provides a tabular display of an operator's performance (QC and patient test statistics).

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|----------------------|---|
| Operator ID | <p>Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen.</p> <p>Click the Suppress checkbox to exclude the Operator ID from the printed report.</p> |
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Instrument | Select an instrument type from the dropdown list. |
| Test | <p>System default displays GLU test.</p> <p>Select a test type from the dropdown list to customize the report.</p> |


| Options | Descriptions |
|--------------------------------|--|
| Performance Level | System default displays All . Select a performance level from the dropdown list to customize the report. |
| Minimum No. of QC Tests | Entry field for the minimum number of QC tests required. Only operators who have performed at least the number of QC tests specified here will be included in the report. |



Patient Outlier Log

The Patient Outlier Log provides a tabular display of patient test results outside the action range specified at the time the patient test was performed.

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|---------------------|---|
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list, select an Instrument, multiple Instruments or all Instruments, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Patient ID | Enter a Patient ID. Click the Suppress checkbox to exclude the Patient ID from the printed report. |
| Patient Name | Enter a Patient Name, filtered by last name. Click the Suppress checkbox to exclude the Patient Name from the printed report. |



| Options | Descriptions |
|----------------------|---|
| Operator ID | <p>Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen.</p> <p>Click the Suppress checkbox to exclude the Operator ID from the printed report.</p> |
| Operator Name | <p>Click the Suppress checkbox to exclude the Operator Name from the printed report.</p> |
| Test | <p>System default displays GLU test.</p> <p>Select a test type from the dropdown list to customize the report.</p> |
| Reagent Lot | <p>Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen.</p> |



Proficiency Data Listing Report

The Proficiency Data Listing report provides a tabular display per instrument of information associated to proficiency test results.

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|--------------------|---|
| Operator ID | Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| | Click the Suppress checkbox to exclude the Operator ID from the printed report. |
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list, select an Instrument, multiple Instruments or all Instruments, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |

| Options | Descriptions |
|------------------------------|--|
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Proficiency Lot | Click on the  button to display the pop-up browser window with the Proficiency Lot Lookup list, select a Proficiency lot, multiple Proficiency lots or all Proficiency lots, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |



QC Outlier Log



The QC Outlier Log report provides a tabular display of control test results outside the reagent lot range specified in QCM3 (at the time of the report).

EDIT REPORT:
QC Outlier Log
 Location Choice All Facilities » All Departments » All Locations
 Instrument ...
 Reagent Lot ...
 Control Lot ...
 Control Level All ▾
 Operator ID ... ☐ Suppress
 Operator Name ☐ Suppress
 Test Glu ▾
Test Date Range
 Range ☐ From To
 Period ☒ Today ▾
Print Options
☐ Suppress Cover Page
☐ Suppress Page Breaks
[More Options >>](#)
 << 1 of 1 >>
[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|--------------------|--|
| Instrument | Click on the  button to display the pop-up browser window with the Instrument Lookup list, select an Instrument, multiple Instruments or all Instruments, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Reagent Lot | Click on the  button to display the pop-up browser window with the Reagent Lot Lookup list, select a Reagent Lot, multiple Reagent Lots or all Reagent Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen. |

| Options | Descriptions |
|----------------------|---|
| Control Lot | <p>Click on the  button to display the pop-up browser window with the Control Lot Lookup list and select a Control Lot, multiple Control Lots or all Control Lots and click on the OK button in the pop-up browser window; click on the Cancel button to return to the Edit Report screen.</p> |
| Control Level | <p>System default displays All.</p> <p>Select All, Low, Normal or high from the dropdown list to customize the Report.</p> |
| Operator ID | <p>Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen.</p> <p>Click the Suppress checkbox to exclude the Operator ID from the printed report.</p> |
| Operator Name | <p>Click the Suppress checkbox to exclude the Operator Name from the printed report.</p> |
| Test | <p>System default displays GLU test.</p> <p>Select a test type from the dropdown list to customize the report.</p> |

i-STAT Data Listing

The i-STAT Data Listing report provides a tabular display per instrument of information associated to i-STAT test results, including test attributes, according to the selected filter options.

EDIT REPORT:

i-STAT Data Listing

Location Choice All Facilities » All Departments » All Locations

Instrument Serial No. ...

Operator ID ... ☐ Suppress

Operator Name ☐ Suppress

Patient ID ☐ Suppress

Result Type ▼

Test Date Range

Range ☐ From To

Period ☒ Today ▼

Print Options

☐ Suppress Cover Page

☐ Suppress Page Breaks



More Options >>

<< 1 of 1 >>

View Print E-Mail Save Cancel

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|------------------------------|--|
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Operator ID | Click on the  button to display the pop-up browser window with the Operator ID Lookup list. Select an operator, multiple operators, or all operators and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| | Click the Suppress checkbox to exclude the Operator ID from the printed report. |


| Options | Descriptions |
|----------------------|---|
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Patient ID | Enter a Patient ID. Click the Suppress checkbox to exclude the Patient ID from the printed report. |
| Result type | System default displays All result types. Select a single result type from the dropdown list to customize the report. |

i-STAT Instrument Event

The i-STAT Instrument Event report provides a tabular display of i-STAT Event Codes per instrument.

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|------------------------------|--|
| Instrument Serial No. | Click on the  button to display the pop-up browser window with the Instrument Serial No. Lookup list, select an Instrument Serial number, multiple Instrument Serial numbers or all Instrument Serial numbers, and click on the OK button in the browser window; click on the Cancel button to return to the Edit Report screen. |
| Event Code | Enter an event code (only one event code can be entered per report). |

Quality Check Codes

The Quality Check Codes report provides a tabular display per location of quality check codes by type.

EDIT REPORT:
Quality Check Codes
 Location Choice All Facilities » All Departments » All Locations
 Quality Codes Lower Limit
 Quality Codes Upper Limit
 Quality Ratio Lower Limit
 Quality Ratio Upper Limit
Test Date Range
 Range ☐ From To
 Period ☒ Today
Print Options
☐ Suppress Cover Page
☐ Suppress Page Breaks
[More Options >>](#)
 << 1 of 1 >>
[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|----------------------------------|--|
| Quality Codes Lower Limit | Entry field for Quality Codes Lower Limit. |
| Quality Codes Upper Limit | Entry field for Quality Codes Upper Limit. |
| Quality Ratio Lower Limit | Entry field for Quality Ratio Lower Limit. |
| Quality Ratio Upper Limit | Entry field for Quality Ratio Upper Limit. |

Quality Check Codes by Operator

The Quality Check Codes by Operator report provides a tabular display of quality check codes by type per operator.

EDIT REPORT:
Quality Check Codes by Operator
 Location Choice All Facilities » All Departments » All Locations
 Operator Name ☐ Suppress
 Quality Codes Lower Limit
 Quality Codes Upper Limit
 Quality Ratio Lower Limit
 Quality Ratio Upper Limit
Test Date Range
 Range ☐ From To
 Period ☒ Today
Print Options
☐ Suppress Cover Page
☐ Suppress Page Breaks
[More Options >>](#)
 << 1 of 1 >>
[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report:

| Options | Descriptions |
|----------------------------------|--|
| Operator Name | Click the Suppress checkbox to exclude the Operator Name from the printed report. |
| Quality Codes Lower Limit | Entry field for Quality Codes Lower Limit. |
| Quality Codes Upper Limit | Entry field for Quality Codes Upper Limit. |
| Quality Ratio Lower Limit | Entry field for Quality Ratio Lower Limit. |
| Quality Ratio Upper Limit | Entry field for Quality Ratio Upper Limit. |

Reagent Usage

The Reagent Usage report provides a summary of i-STAT cartridge usage by type per location.

EDIT REPORT:
Reagent Usage
Location Choice All Facilities » All Departments » All Locations
Test Date Range
Range ☐ From To
Period ☒
Print Options
☐ Suppress Cover Page
☐ Suppress Page Breaks
[More Options >>](#)

<< 1 of 1 >>

[View](#) [Print](#) [E-Mail](#) [Save](#) [Cancel](#)

Select a Facility Department or Location in the Tree.

Select the filter options for the report.

6. Operators

Operators Overview

The **Operators** menu option provides access to three screens for use in managing operators and instrument certifications. The **New** screen allows you to add a new operator and instrument certification. The **Certification** screen allows you to view and update instrument certifications. The **Search** screen provides a simple view of operators. Each screen is available only if your account was granted this permission by the Administrator or POCC.

New Operator Screen

Operators>>New

The **New Operator** screen allows the user to add and certify a new operator. Depending on the user's assigned Function Permissions the user may only be able to add new Operators, but not add instrument certification. Separate permission must be assigned to add and certify an Operator. Applicable Operator Menu Permissions are as follows;

- **Edit Operator Certification** - Allows user permission to add or edit Operator information, including home Location. User may also add or edit (ReCertify) instrument certification.
- **Edit Operator Information** - Allows user permission only to add or edit Operator information, including home Location.

Each Operator is assigned to a home Location within a Facility. Each Facility in the Organization maintains a separate operator list.

Add and Certify New Operator

To **add** and certify a new operator:

1. Move the mouse cursor over **Operators** to display the menu options.
2. Click on the **New** option.
 - The **New Operator** screen displays the entry fields for new Operator Information (on the left of the screen) and Operator Certification information (on the right).
 - Each time you select the **New Operator** screen from the menu, the table will be empty.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

Refresh Hide Tree x

NEW OPERATOR: F2 » D2 » L2 No:175 LIS:123 Op:0

Operator Information

* First Name * Operator ID

Middle Name E-Mail Address

* Last Name * Home Department D2

* Home Location L2

Operator Certification

Instrument PCx

Original Cert. Date 02/15/2007 31

* Certification Date 02/15/2007 31

* Recert Interval Expired

Expiration Date

Automatically Recertify With Good QC ☐ Yes ☒ No

Add Operator

0 Records Found

| Edit | Delete | (Re)Certify | E-Mail | Customize | | | | | |
|------------------------------|-------------|-------------|----------|-----------------|---------------------|------------|-----------------|-----------|-------------|
| Operator Name | Operator ID | Department | Location | Instrument Type | Original Cert. Date | Cert. Date | Recert Interval | Exp. Date | Auto Recert |
| <input type="checkbox"/> All | | | | | | | | | |

No data is available for this selection. Please change selection criteria or tree selection.

Edit Delete (Re)Certify E-Mail Customize

0 Records Found

*First Name, Last Name, Operator ID, Home Department and Home Location are **required** entry fields under Operator Information.

*Certification Date and Recert Interval are **required** entry fields under Operator Certification.

- Fill in the First Name and Last Name and Operator ID fields under Operator Information.
 - The Operator ID entry field accepts up to a **maximum** of 30 characters.
- Click on a **Location** in the Tree to assign the operator to a Home Department and Location.
 - The Home Location field and the associated Home Department field automatically display the selected location.
- Select the following operator certification options:

| Options | Description |
|----------------------------|---|
| Instrument | Select instrument type for certification |
| Original Cert. Date | <p>The default is today's date.</p> <p>Enter the date the operator was first certified.</p> <p>Calendar allows you to select correct date. Selected value will appear in the Original Certification Date field.</p> |
| *Certification Date | <p>The default is today's date.</p> <p>Enter the date of the current (re)certification.</p> <p>Calendar allows you to select correct date. Selected value will appear in the Certification Date field.</p> |
| *Recert Interval | <p>Select from:</p> <p>3 months</p> |

| Options | Description |
|--|--|
| | 6 months |
| | 12 months |
| | Never Expire |
| | Specify |
| | Expired (default) |
| | Choosing "Specify" will activate the Expiration Date field below. |
| Expiration Date | Enter the expiration date. Calendar allows you to select correct date. Selected value will appear in the Expiration Date field. The default is today's date. |
| Automatically Re-certify With Good QC | Select Yes or No (default). These buttons become active only when a specific interval is selected from the Expiration Date. Selecting Automatically Re-certify With Good QC will change an Operator's Certification Date from the current date to the date of the passing QC test (low and high), thus extending the Operator's certification date. The defined Recertification Interval remains unchanged. Each time the Operator completed passing QC the Certification Date will be updated. |

6. Click the **Add Operator** button to save the data.

- After the **Add Operator** button is pressed, the certification criteria entered will remain but the Operator Information fields will be cleared for the next new operator (until the screen is closed).

7. The saved Operator will appear in the New Operator table.

To add an additional instrument certification to one or more Operators in the table:

1. Select one or more records from the table by clicking the box(es) in the **Operator Name** column.
2. Select the desired instrument type certification from the Instrument dropdown list in the Operator Certification section.
3. Click the (Re)Certify action button.

4. Select the following operator certification options:

| Options | Description |
|--|--|
| Original Cert. Date | Enter the date the operator was first certified. Calendar allows you to select correct date. Selected value will appear in the Original Certification Date field. |
| *Certification Date | Enter the date of the current (re)certification. Calendar allows you to select correct date. Selected value will appear in the Certification Date field. |
| Expiration Date | Select from: 3 months 6 months 12 months Never Expire Specify Expired (default) Choosing "Specify" will activate the box to the right of the Expiration Date field. |
| Automatically Re-certify With Good QC | Select Yes or No (default). These buttons become active only when a specific interval is selected from the Expiration Date. Selecting Automatically Re-certify With Good QC will change an Operator's Certification Date from the current date to the date of the passing QC test (low and high), thus extending the Operator's certification date. The defined Recertification Interval remains unchanged. Each time the Operator completed passing QC the Certification Date will be updated. |

5. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.
6. The instrument specific certification will appear in the New Operator table as a new row in the table.

Note: For each instrument type certification a new row will appear in the table for every Operator with the specified certification type. However, each Operator ID is unique and is noted in the database only once. As such deleting an Operator from the system will delete all instrument type certifications for the selected Operator.

Add a New Operator (no certification permission)

The certification fields in the Operator Certification section are only displayed if you have permission to certify operators. If you do not have permission to certify Operators, the New Operator screen will appear only with the Operator Information fields.

To **add** a new operator:


1. Move the mouse cursor over **Operators** to display the menu options.
2. Click on the **New** option.
 - The **New Operator** screen displays the entry fields for new Operator Information.
 - Each time you select the **New Operator** screen from the menu, the table will be empty.

*First Name, Last Name, Operator ID, Home Department and Home Location are **required** entry fields under Operator Information.

3. Fill in the First Name and Last Name and Operator ID fields under Operator Information.
 - The Operator ID entry field accepts up to a **maximum** of 30 characters.
4. Click on a **Location** in the Tree to assign the operator to a Home Department and Location.
 - The Home Location field and the associated Home Department field automatically display the selected location.
5. Click the **Add Operator** button to save the data.
 - After the **Add Operator** button is pressed, the Operator Information fields will be cleared for the next new operator (until the screen is closed).
6. The saved Operator will appear in the New Operator table.


Warning Message

If you attempt to add an operator who is already in the system, the following message is displayed (No two operators in the same Facility can have the same Operator ID):




- Click **OK** to enter alternative data for the new operator.

If you re-add a deleted operator, the following message is displayed:



- Click **OK** re-add the operator.
- Click **Cancel** to return to the New Operator screen.

If a different menu option is selected before saving the new or edited operator information, the following warning message is displayed:



- Click the **OK** button to cancel the data entered.
- Click the **Cancel** button to clear the message and return to the **Edit Operator** screen to save the data.

If a **required** field is missing when the new operator information is saved, the following warning message is displayed indicating the specific required field (First Name, Last Name, and Operator ID) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the **required** field for Home Location is not entered for a **New Operator** the following warning message is displayed:




- Click the **OK** button and select a Location from the Tree.

If data entered in the Certification Date or Expiration Date fields is not valid, the following warning message is displayed, indicating the field that requires valid data:



- Click the **OK** button to clear the message and enter the valid information.

If the Expiration Date selected is earlier than yesterday's date, the following error message is displayed:



- Click on the **OK** button to clear the message.

New Operator table

The added operators will be displayed in a table in the bottom half of the **New Operator** screen.

The Action Bar will have the following buttons: **Edit**, **Delete**, **(Re)Certify** and **E-Mail**. These actions are described in the Operator Certification screen section.

Operator Certification Screen

Operators>>Certification

The **Operator Certification** screen allows you to search operators by Last Name, ID or certification information and to re-certify the operators.

The data in the Expiration Certification Date column will be highlighted (white text on red background) for expired operator certifications.

To **view** the Operator Certification screen:

1. Move the mouse cursor over **Operators** to display the menu options.
2. Click on the **Certification** option.

The **Operator Certification** screen displays the Operator Certification table with operators based on the location highlighted in the Tree. The default Instrument type is **All**, and the default Status is **Never Certified**.

To display operators at a *different location* using the same search criteria, click on a **Facility**, **Department** or **Location** in the Tree.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

OPERATOR CERTIFICATION: Eastern Medical Center » Pediatrics » Ward A No: 134 LIS: 58 Op: 30

Search Criteria

Operator ID: Instrument: All From:

Last Name: * Status: Never Certified To:

Display Records

0 Records Found

| Operator Name | Operator ID | Department | Location | Email | Instrument Type | Original Cert. Date | Cert. Date | Recert Interval | Exp. Date | Auto Recert | Last Good QC | Last Pat Test |
|--|-------------|------------|----------|-------|-----------------|---------------------|------------|-----------------|-----------|-------------|--------------|---------------|
| No data is available for this selection. Please change selection criteria or tree selection. | | | | | | | | | | | | |

Edit Delete (Re)Certify E-Mail Customize

0 Records Found

Using the Search Criteria



Use the search criteria to find a specific operator or to find all operators with an instrument certification set to expire on a specific date.

To **find** an operator, fill in the following search criteria and click the **Display Records** button:

| Criteria | Options | Description |
|--------------------|---------|---|
| Operator ID | | Enter ID number |
| Last Name | | Enter operator name |
| Instrument | | Default is All . Selecting All will display operators for any instrument. Or, select a specific instrument type. |

Notes:


- *There may be multiple rows per operator, one per instrument.*
- *The screen query will be re-executed and the screen will be refreshed automatically whenever the Instrument dropdown is reselected*

| | | |
|----------------|-----------------|--|
| *Status | Never Certified | Use to select new operators that have never been certified. (Default) |
| | Expiration | Both date fields are disabled. Both date fields are enabled. Enter an expiration date to select operators with expired certifications, or Enter or select a date range to select operators whose certifications have expired between specific dates using calendar pop-ups in From and To fields. |
| | No Selection | Select all operators (except deleted operators) regardless of certification status. Both date fields are disabled. |
| From | | Click on  to select beginning of date range. Only enabled if Expiration Status is selected, above. |
| To | | Click on  to select end of date range. Only enabled if Expiration Status is selected, above. |

- The operators that match (or are similar to) the search criteria will be displayed in the Operator Certification table.
- To display operators at a different location using the same search criteria, click on a **Facility**, **Department** or **Location** in the Tree.

Warning Message

If you click the **Edit** button when **All** is selected from the Instruments dropdown list in the Search Criteria, the following message is displayed:



- Click OK to return to the Operator Certification screen and select a specific Instrument from the dropdown list in the Instrument field.

Edit Operator Records

To **edit** a single, multiple or all operator records:

1. Select one or more records from the table by clicking the box(es) in the **Operator Name** column. Alternatively, click the box marked **All** in the Operator Name column header to edit all records.
2. Click the **Edit** action button.

- The **Edit Operator** screen appears with the first selected record.
 - If multiple records are selected, each record is displayed on a separate page.
3. Edit the Operator information.
 - * *First Name, Last Name, Operator ID, Home Department and Home Location are **required** entry fields.*
 - Operator ID is displayed but is not an editable field.
 - To assign the operator to a different location, click on a Location in the Tree. This will automatically update the Home Location and associated Home Department fields


Note:

An operator can be moved to a different Home Location in the same Facility. An operator cannot be moved to a different Facility. To move an operator to a different Facility a New Operator record must be added under the new Facility. The operator may be added to a new Facility with the same ID number, provided the ID number is not currently assigned to another operator.

4. Click the **Save** button to save the data or Click the **Cancel** button to cancel the edit action without saving.

Warning Message

If a different Facility is selected in the Edit Operator screen, the following error message is displayed:




- Click the **OK** button to clear the message.
- Select a Home Location in the operator's current Facility.

Delete an Operator

The **delete** button will hide all records for the selected operator(s), including the operator's identity and all certification information for the operator(s).

To **delete** a single operator record or multiple operator records:

1. Click the box(es) in the Operator Name column.
2. Click the **Delete** action button.
3. The following dialogue box appears:



4. Click **OK** to delete the record(s) and return to the **Operator Certification** screen, or click **Cancel** to return to the **Operator Certification** screen without deleting the record(s).

Note: *If you delete an operator record, this operator and all associated certifications will be flagged as deleted in the database. Test results completed by the deleted Operator are still accessible after the Operator has been removed. The Operator information will still be visible as part a of a historical test record. However, the Operator information will no longer be visible on the Operator list and the Operator ID number is no longer an active ID number.*

Re(Certify) an Operator

To certify or re-certify a single operator or multiple operators in a batch:

1. Select one or more operators from the table by clicking the box(es) in the **Operator Name** column.
2. Select the desired instrument type certification from the Instrument dropdown list in the Operator Certification section.
3. Click the **Re(Certify)** action button.

4. Select the following operator certification options:

| Options | Description |
|--|--|
| Original Cert. Date | Enter the date the operator was first certified. Calendar allows you to select correct date. Selected value will appear in the Original Certification Date field. |
| *Certification Date | Enter the date of the current (re)certification. Calendar allows you to select correct date. Selected value will appear in the Certification Date field. |
| Expiration Date | Select from: 3 months 6 months 12 months Never Expire Specify Expired (default) Choosing "Specify" will activate the box to the right of the Expiration Date field. |
| Automatically Re-certify With Good QC | Select Yes or No (default). These buttons become active only when a specific interval is selected from the Expiration Date. |


| Options | Description |
|---------|---|
| | Selecting Automatically Re-certify With Good QC will change an Operator's Certification Date from the current date to the date of the passing QC test (low and high), thus extending the Operator's certification date. The defined Recertification Interval remains unchanged. Each time the Operator completed passing QC the Certification Date will be updated. |

- Click the **Save** button to certify the selected operator(s) for the selected instrument and return to the **Operator Certification** screen or Click the **Cancel** button to return to the **Operator Certification** screen without saving.

Note: In order to certify an operator for another instrument, you need to change the selection in the Instrument dropdown field.

Warning Message

If you click the **Re(Certify)** button when **All** is selected from the Instruments dropdown list in the Search Criteria, the following message is displayed:



- Click **OK** to return to the **Operator Certification** screen and select a specific Instrument from the dropdown list in the Instrument field.

E-Mail Operator Certification Records

To **e-mail** a single operator record or multiple operator records:


- Select one or more records in the table by click on the box(es) in the Operator Name column.
- Click the **E-Mail** action button.

To **e-mail** all Operator records:

- Click on the box marked **All** in the Operator Name column header.
- Click the **E-Mail** action button.

The default e-mail client opens a new mail message and includes the menu option title in the subject field.

- The body of the message includes a hyperlink to an .htm file that contains the selected records. Linked files will expire after 30 days.



Operator Search Screen

Operators>>Search

The **Operator Search** screen allows the user to find, edit, delete and e-mail operator information only. These action buttons are described in the Operator Certification section. This screen also allows the user to assign an operator to a different Home Location. Certification information (Operator ID, Instrument Certification) cannot be updated from this screen.

To **view** operator(s):

1. Move the mouse cursor over **Operators** to display the menu options.
2. Click on the **Search** option.

The **Operator Search** screen displays the Operator table based on the location highlighted in the Tree.

To display operators at a different location using the same search criteria, click on a **Facility**, **Department** or **Location** in the Tree.

Using the Search Criteria

To **find** an operator:

1. Enter a **Last Name** and/or an **Operator ID** in the entry fields.
2. Click the **Display Records** button.


The **Operator Search** screen displays the Operator table based on results that match or are similar to the Last Name and/or Operator ID entered and the location highlighted in the Tree.

To find operators at a different location using the same search criteria, click on a **Facility**, **Department** or **Location** in the Tree.

Operator Certification Information is not displayed on the operator search screen.

Warning Messages

If no match is found for the Last Name and/or Operator ID entered in the Operator Search screen, the following message is displayed:




- Click the **OK** button to add the operator.
- Click the **Cancel** button to clear the message.

If the user clicks the **OK** button to add the operator, the application opens a **New Operator** screen with the Last Name and/or Operator ID entered in the **Operator Search** screen. The Location highlighted in the Tree is displayed in the Home Location field.

- Enter the required fields.
- Click the **Add Operator** button.
- The **Operator Search** screen appears again with the new operator displayed in the Operator table.


Note: When you are presented with the **New Operator** screen as a result of searching for a non-existent operator in the **Operator Search** screen the New Operator screen will not display a Table. After you have entered the new operator's information and hit the **Add** button, you will return to the **Operator Search** screen rather than seeing the new operator presented in the New Operator screen.

If records are not selected from the table before clicking the **Edit**, **Delete** or **E-Mail** action button the following warning message is displayed:



- Click the **OK** button to clear the message and make a selection.

If an operator record is selected for deletion, the following warning message is displayed:



- Click the **OK** button to delete the operator records.
- Click the **Cancel** button to cancel the delete action.

Note: The delete action will flag the operator record as deleted and the operator will not appear on the Operator table. The test results associated with the operator will remain intact for data integrity and reports.

7. Instruments

Instruments Overview

The **Instruments** menu option provides access to three screens for use in managing instrument settings and two screens with log files. The **Configuration** screen allows you to configure instrument settings. The **Comments** screen allows you to define comment codes by instrument type. The **Ports** screen allows you to manage instrument communications. The **Tracking Log** screen allows you to review events relating to instrument location assignments. The **Upload Log** screen allows you to review instrument communication session results. Each screen is available only if your account was granted this permission by the Administrator or POCC.


Moving, Renaming and Deleting Instruments

The first time a **new** instrument uploads data to the system the instrument serial number appears in the Tree under an **instrument type icon** (🖨️) in the **Location** from which the instrument was uploaded.

Instruments are grouped by **Location** under 🖨️ **instrument type icons**.

For example in the illustration below the **Department**, Pediatrics, includes:

Two **Locations** (Ward A and Ward B) each with a PCx instrument icon (🖨️) and one instrument assigned to each location (M08990068 and M06210023, respectively).




Instruments visible on the Tree can be moved, renamed and deleted.

Only users assigned the Instrument Operations Function Permission can move, rename or delete an instrument.

Moving Instruments

- An instrument can be moved to a different Location in the same Facility via drag and drop.
- An instrument can be moved to a Location within a different Facility; however, this cannot be done via drag and drop.

To **move** an instrument within the same Facility:

1. Place the mouse cursor over the  icon to the left of the instrument **serial number** or **name**.
2. Click and hold the left mouse button to drag the instrument to a new Location.
 - The cursor will appear as a four-arrow symbol.
 - When the mouse cursor changes from a four-arrow symbol to an arrow with a tag box, the instrument **can** be moved to that Location.
3. Release the mouse button to display the instrument in its new Location.
 - When the mouse cursor changes to a circle with a line through it, the instrument **cannot** be moved to that part of the Tree. An instrument can only be assigned to the Location level of the Tree.
 - An instrument cannot be assigned to the Facility or Department level of the Tree.
 - An instrument cannot be moved from a Location within one Facility to a Location within a different Facility via drag and drop.
4. Release the mouse button to display the instrument in its original Location.

To **move** an instrument to a different Facility:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Configuration** option.
3. Select the current instrument Location in the Tree.
4. Change the Configuration Download setting to **Download Location** and click the **Save** button.
5. Select the new instrument Location in the Tree.
6. Change the Configuration Download Setting to **Download Location** and click the **Save** button.
7. Delete the instrument from its current Location in the Tree by highlighting the instrument in the Tree and selecting the Delete key on your keyboard.
8. Perform a control test and upload the instrument from the desired Location in the new Facility. This should be the same Location as identified in step 5 above.

- Upon upload the instrument will appear in Tree under the Location from which the instrument was uploaded in the new Facility.
9. Select the current instrument Location in the Tree. This should be the same Location as identified in step 5 above.
 10. Click the Default button to return the Location back to the Assigned Location setting. **Do not click the Save button.**
 11. Select the instrument's previous Location in the Tree. This should be the same Location as identified in step 3 above.
 12. Click the Default button to return the Location back to the Assigned Location setting. **Do not click the Save button.**


Notes:

- *When an instrument is moved to a Location that does not have any other instruments of that type, the system creates a new or additional **instrument type icon** for that Location.*
- *When an instrument is moved from a Location where it is the only instrument of that type, the system removes the **instrument type icon** from that Location.*
- *Saving a configuration setting at the Location level of the Tree will break the Parent-Child relationship. Therefore, it is important to select the Default button for each Location (previous and new) after the instrument appears in the new Facility. See [Re-establishing the Parent-Child Relationship](#) for more information.*

Renaming Instruments

To **rename** an instrument:

1. Place the mouse cursor over the instrument's **serial number** or **name** text and double click the left mouse button.
 - The serial number or name text changes from read-only to an editable text box.



2. Type the new instrument name and press the keyboard **Enter** key or click the left mouse button anywhere outside the text box to display the new instrument name in the Tree.
3. To cancel the rename action, press the keyboard **Esc** key on your keyboard to turn the editable text box back to read-only text. The instrument name remains unchanged in the Tree.


Note: If the serial number is replaced with text (you can leave the serial number and add text) the instrument serial number is visible in the Instrument Serial No. column with all test results. The Instrument Name column will display the name entered in the Tree (serial number and/or text).

Deleting Instruments

To **delete** an instrument:

1. Place the mouse cursor on the **serial number** or **name** text, click the left mouse button to highlight the instrument, and press the keyboard **Delete** key.

The following **warning message** is displayed:




2. Click the **OK** button to delete the instrument.
3. Click the **Cancel** button to cancel the delete action.

Notes:

The instrument record is flagged as deleted and the instrument is removed from the Tree. If the instrument has data in the live database, the test result data remains intact for data integrity and reports. If the instrument flagged as deleted uploads data at a later time, the instrument will reappear in the Tree in the Location from which it was uploaded. The Instrument name and serial number are both saved in the database.

Warning Messages

1. If the instrument name entered already exists for an instrument of the same type at the same location and it has not been flagged as deleted, the following error message is displayed:



- Click the **OK** button to clear the message.
2. If the instrument name field is cleared during the rename action the following error message is displayed:



- Click the **OK** button to clear the message.

Instrument Configuration Screens

Instruments>>Configuration

The **Instrument Configuration** screens allow the POCC to view and edit instrument configuration data for a specific location. Instruments can be configured for the entire Organization or an individual Facility, Department or Location.

This section of the user manual contains full information for configuring the Precision PCx Blood Glucose Monitoring System and limited information on Test ID mapping for the i-Stat 1 Analyzer. Please refer to the documentation provided with the i-Stat 1 Analyzer or other installed instrument type for full instrument configuration options.

Location Tree Parent-Child Relationship

The Location Tree is made up of one or more Facilities and associated Departments and Locations. Each level of the Tree is linked to the level above. This is known as the Parent-Child relationship. The Organization level is a Parent Level. The Location level is a child of the Department level. Facilities and Departments can be either a Parent or Child.

| | |
|--------------|---|
| Organization | Parent to all Facilities within the Organization |
| Facility: | Parent to all Departments within the Facility, Child to the Organization |
| Department: | Parent to all Locations within the Department, Child to the associated Facility |
| Location: | Child to the associated Department |

Instrument Configuration and the Parent-Child Relationship

Instrument configuration settings saved at the top level of the Tree affect all sub-levels. For example, configuration settings saved at the Organization (parent) level will be applied to all Facilities and associated Departments and Locations.

Instrument configuration settings saved at a lower level of the Tree affect only the levels below the selected level of the Tree. For example, configuration settings saved at the Department level will be applied to all Locations (child) associated to the Department (parent). Settings saved at the Location level affect only the selected Location.

It is recommended that instrument configuration settings be made at the Facility level when possible. Department level settings should then be used to vary the instrument configurations only as needed. The Organization level settings then remain available for reference as default values.

Note: *It is important to note the selected Tree level prior to saving instrument configuration settings to ensure the settings are applied to the appropriate locations.*

Breaking the Parent-Child Relationship

Saving instrument configuration for a lower level of the Tree (i.e. Department or Location) will sever the relationship between the child and its parent level on the Tree.

Breaking the Parent-Child relationship is necessary to enable configuration settings different from the global Organization or Facility settings (i.e. certain locations may require an Action Range different from the majority of the Organization or Facility).

Once a location is disassociated to its parent location, any changes to Instrument Configuration will need to be made independent of the parent location.

Identifying a Broken Parent-Child Relationship

There is no visual indication in the Tree that the relationship between a child location and its parent has been broken. To identify the locations that are no longer linked its parent run the Instrument Configuration report from the QCM default report suite. This reports Organization, Facility, Department or Location specific instrument settings. The first page of the report identifies instrument settings applied to the entire Organization (associated Facilities, Departments and Locations). Subsequent locations listed on this report are no longer linked to its parent location in the Tree.

Re-establishing the Parent-Child Relationship

If the Parent-Child relationship was severed in error or the configuration settings for a location should be the same as its parent, the child location can be reattached to its parent location by highlighting the child location in the Tree and selecting the Default button (located to the right of the Save button) on the Instrument Configuration screen.

To re-establish the Parent-Child relationship:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Configuration** option.
3. Select the child location in the Tree.
4. Click the **Default** button on the Instrument Configuration screen.

Note: Clicking the Save button will break the parent-child relationship. Do not click the Save button after clicking the Default button.

Configuring Instrument Settings for Precision PCx

To **view** the Instrument configuration screens for Precision PCx:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Configuration** option.
 - The **Instrument Configuration** screen displays the first of seven pages used to configure the instruments for a Location.
 - Each instrument configuration page displays the corresponding page number icon highlighted at the bottom of the screen.
 - The Instrument Configuration pages are titled as follows:

| | |
|--------|----------------------|
| Page 1 | All Test Options |
| Page 2 | Patient Test Options |
| Page 3 | Control Test Options |
| Page 4 | Security Options |
| Page 5 | Profile Options |
| Page 6 | Test Type Options |
| Page 7 | System Options |

INSTRUMENT CONFIGURATION: F1 » All Departments » All Locations No:460 LIS:15 Op:559

Instrument:

Configuration Download ☐ Download Location ☐ Assigned Location ☒

All Tests Options

Upload Interval

Expired: * After: ☒ Hours ☐ Days

Operator ID

Options: * Min Digits: * Max Digits:

Barcode Type

☒ 1 2 of 5

☒ Code 128

☒ Codabar

☒ Code 93

☒ Code 39

Truncate Digits

* First: * Last:

The **Configuration Download** option appears on all seven Configuration screens. The Download Location or Assigned Location option selected on any of the seven screens applies to all configuration screens.

| Setting | Options | Description |
|-------------------------------|---|---|
| Configuration Download | <input type="radio"/> Download Location | Click on the Download Location button to automatically reconfigure an instrument to match the instruments at the Location from |

| Setting | Options | Description |
|---------|---|--|
| | | <p>which the instrument is uploaded. The application will download the instrument configuration settings and lists (Operator and Reagent) assigned to the Location from which it is uploaded.</p> <p>All test results uploaded from the instrument will be assigned to the instrument's download Location in the Tree.</p> |
| | <input type="radio"/> Assigned Location | <p>Click on the Assigned Location button to retain the instrument's original configuration regardless of the Location from which the instrument is uploaded.</p> <p>All test results uploaded from the instrument will be assigned to the instrument's assigned Location in the Tree, regardless of the Location from which the test results were uploaded.</p> |

The **save** and **default** buttons also appear on all seven screens.

- Click the **Save** button to save the data.
- Click the **Default** button to apply the parent configuration settings for **All Tests Options** to the selected location in the Tree. This button will re-establish the Parent-Child relationship and apply the parent's configuration settings to the child location (which may not necessarily be the system default settings).

Note: You do not need to select the Save button for each configuration page. You can move from page to page within Instrument Configuration without selecting Save and the changes will be saved. Once you have entered all configuration settings select the Save button to apply the settings to the highlighted Tree Location. This behavior also applies to the Default button.

All Tests Options

Page **1** of the **Instrument Configuration** screen displays the system default settings for **All Tests Options**.

The available options on the **All tests Options** configuration page are:

| Setting/Options | Description |
|------------------------|--|
| Upload Interval | |
| Expired | Select the time interval that will activate the instrument upload expiration settings and select: Allow Test , Warn or Lockout from the dropdown list. |
| *After | <div>O Hours O Days</div> Enter a value for the time interval and click on the Hours or Days button. |
| Operator ID | |
| Options | Select Disabled , Optional or Required from the dropdown list. If Required or Optional is selected, enter the minimum and maximum digits required for Operator ID (the maximum digits for an Operator ID is 15 characters). The maximum and minimum digits fields are <i>required</i>. |
| Barcode Type | Select the Operator ID Barcode Type. |
| Truncate Digits | To ignore any beginning and/or ending values, specify the digits to truncate in the Operator ID. |

** After, Min Digits and Max Digits are required entry fields on the All Tests Options screen.*

Patient Test Options

Page 2 of the instrument configuration screen displays the system default settings for **Patient Test Options**.

The available options on the **Patient Test Options** configuration page are:

| Options | Description |
|------------------------------------|--|
| Patient ID | |
| Options | Select Disabled , Optional or Required from the dropdown list. If Required or Optional is selected, enter the minimum and maximum digits required for Patient ID (the maximum digits for Patient ID is 15 characters). The maximum and minimum digit fields are required . |
| Barcode Type | Select the Patient ID Barcode Type. |
| Truncate Digits | To ignore any beginning and/or ending values, specify the digits to truncate in the Patient ID. |
| Comment Code (In-Range) | Select an option for the In-Range Comment Code: Disabled , Optional or Required . |
| Comment Code (Out-of-Range) | Select an option for the Out-of-Range Comment Code: Disabled , Optional or Required . |
| Memory Full Option | Select an option for when the instrument's memory is full: Lockout , Prompt for Overwrite or Allow Overwrite . |

* **Min Digits**, **Max Digits**, **First** and **Last** are required entry fields on the **Patient Test Options** screen.

Control Test Options

Page **3** of the instrument configuration screen displays the system default settings for **Control Test Options**.

The screenshot shows the 'INSTRUMENT CONFIGURATION' page for 'F1' at 'All Departments' and 'All Locations'. The instrument is set to 'PCx'. Under 'Configuration Download', 'Assigned Location' is selected. The 'Control Test Options' section includes the following settings:

- Normal (Mid) Level: Disabled
- Low Level: Required
- High Level: Required
- Result Display: Numeric
- Comment Code (In-Range): Disabled
- Comment Code (Out-of-Range): Disabled
- Control Test Memory Full: Allow Overwrite

At the bottom, there are numbered tabs 1 through 7, and 'Save' and 'Default' buttons.

The available options on the **Control Test Options** configuration page are:

| Setting/Options | Description |
|--|--|
| Normal (Mid) Level | Select the Disabled or Required option for the instrument to prompt for the specified control level. |
| Low Level | Same as above. |
| High Level | Same as above. |
| Result Display | Select the Numeric or Pass/Fail option for the control test result display. |
| Note: If the Pass/Fail option is selected the instrument will display the result as Pass or Fail, however, the actual numeric value of the Pass/Fail control test result will be uploaded and displayed in the application. | |
| Comment Code (In-Range) | Select an option for the In-Range Comment Code: Disabled , Optional or Required . |
| Comment Code (Out-of-Range) | Select an option for the Out-of-Range Comment Code: Disabled , Optional or Required . |
| Control Test Memory Full | Select an option for when the instrument's control test memory is full: Lockout , Prompt for Overwrite or Allow Overwrite . |

- Click the **Save** button to save the data.
- Click the **Default** button to apply the parent configuration settings for **Control Test Options** to the selected location in the Tree. This button will re-establish the Parent-Child relationship and apply the parent's configuration settings to the child location (which may not necessarily be the system default settings).

Security Options

Page 4 of the instrument configuration screen displays the system default settings for **Security Options**.

The available options on the **Security Options** configuration page are:

| Setting/Options | Description |
|----------------------|--|
| QC Interval | |
| Time Expired | Select the Allow Test , Warn or Lockout option for the QC test interval expired security setting. |
| Time Interval | Select the Time Interval , Time of Day or No. of Patient Tests to activate the QC test interval expired security setting. <ul style="list-style-type: none"> 1. Enter the interval in hours for the Time Interval option, or 2. Enter up to three different times, such as, 06:00 AM, 02:00 PM and 10:00 PM for the Time of Day option, or 3. Enter a number between 1 and 250 for the No. of Patient Tests option |
| Operator ID | |
| Not Certified | Select the Allow Test , Warn or Lockout option for the operator ID not certified security setting. |
| Strip Lot | |
| Not On List | Select the Allow Test , Warn or Lockout option for the strip lot not on list security setting. |

Profile Options

Page 5 of the instrument configuration screen displays the system default settings for **Profile Options**.

The available options on the **Profile Options** configuration page are:

| Setting/Option | Description |
|---|---|
| Download Information | Select one download option for downloading an operator list: None, Full List or Incremental List. |
| Operator List | |
| None | Select for no download to instruments. |
| Full List | Select to download the full operator list. |
| Incremental List | Select to download changes made since the meter's last upload. |
| <input type="checkbox"/> Meter Setup | Check this option to download configuration setting to the instruments. |
| <input type="checkbox"/> Strip List | Check this option to download the approved strip lot list to instruments. |

Test Type Options

Page 6 of the instrument configuration screen displays the system default settings for patient **Test Type Options**.

The screenshot shows the 'INSTRUMENT CONFIGURATION: F1 » All Departments » All Locations' screen. At the top right, it displays 'No: 460 LIS: 15 Op: 559'. Below the header, there is a section for 'Instrument' with a dropdown menu set to 'PCx'. Underneath, there are radio buttons for 'Configuration Download', 'Download Location', and 'Assigned Location', with 'Assigned Location' being selected. The main section is titled 'Test Type Options' and contains a 'Sample Type' dropdown set to 'Prompt' and a 'Units' dropdown set to 'mg/dL'. Below these is the 'Action Range' section, which includes an 'Enabled' checkbox (unchecked), a 'Low' input field with the value '20', and a 'High' input field with the value '600'. At the bottom of the screen, there are seven numbered tabs (1-7), with tab 6 being active, and two buttons labeled 'Save' and 'Default'.


The available options on the **Test Type Options** configuration page are:

| Option | Description |
|------------------------------|--|
| Sample Type | Specify the blood sample type used for patient testing and the measurement units. |
| Capillary or Arterial | Select if patient blood sample is capillary or arterial. |
| Venous | Select if patient blood sample is venous only. |
| Prompt | Select to choose the correct blood sample type prior to each patient test. |
| Units | Select mg/dL or mmol/L . |
| Action Range | Enable Action Range to prompt operators when a patient test result is outside the user-defined range. |
| Enabled | Check this option to enable the user-defined action range and specify the limit for Low and High . |


Note: Changes made to the Low and High action ranges will only apply to test results completed since the change was downloaded to the instrument.

Warning Messages

If **Enabled** is checked and there is no entry in the **Low** field, a warning message is displayed:



If **Enabled** is checked and there is no entry in the **High** field, a warning message is displayed:



- Click the **OK** button to clear the message and enter the required information.

System Options

Page **7** of the instrument configuration screen displays the system default settings for **System Options**.

The screenshot shows the 'INSTRUMENT CONFIGURATION: F1' screen with the breadcrumb 'All Departments » All Locations'. In the top right corner, there are three red status indicators: 'No:460', 'LIS:15', and 'Op:559'. Below the breadcrumb, there is a dropdown menu for 'Instrument' set to 'PCx'. Underneath, there are three radio buttons for 'Configuration Download': 'Download Location' (unselected) and 'Assigned Location' (selected). The main section is titled 'System Options' and contains several settings: 'Language' is set to 'English'; 'Enable Beeper' is checked; 'Date' settings include 'Format' set to 'mo/dy/yr' and 'Separator' set to '/'; 'Time' settings include 'Format' set to '1-12AM/PM', 'Separator' set to ':', and 'Units' set to 'None'; 'Firmware' settings include an unchecked 'Update Firmware' checkbox, 'Firmware Version' set to 'PCx 2.2.1', and '* Instrument Password' set to '123'. At the bottom of the form, there are seven numbered tabs (1-7) and two buttons: 'Save' and 'Default'.

The available options on the **System Options** configuration page are:

- | | |
|------------------------|---|
| Language | Select a language from the dropdown list. |
| Enable Beeper | Click the checkbox to enable the beeper. |
| Date format | Select the date format from the dropdown list. |
| Separator | Select the separator to place between the date numbers. |
| Time format | Select the time format from the dropdown list. |
| Separator | Select the separator to place between the hours and minutes. |
| Units | Select the unit for the time format from the dropdown list. |
| Update Firmware | Click the checkbox to enable the firmware update option and select a Firmware Version from the dropdown list. |

Note: See the meter component instructions for specific information about upgrading instruments.

***Instrument Password** Enter the instrument password

** **Instrument Password** is a required entry field on the **System Options** screen.*

Configuring Test ID Mapping for i-STAT Portable Clinical Analyzer:

Specify the Patient ID used to identify control, proficiency and linearity tests run on the i-STAT PCA.

The Patient ID entered will be used to identify both the test type and level specified.

The Enable Automatic Conversion checkbox applies the ID setting to uploaded test results. If this box is not checked the ID's entered will not convert the test to the specified test type/level and the test result will appear as a patient test result.

To **view** the Test ID Mapping configuration screens for i-STAT PCA:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Configuration** option.
3. Select i-STAT from the Instrument dropdown list,
 - The Patient Id Mapping screen will appear. Enter the desired Patient ID for each Test Type and Level.

INSTRUMENT CONFIGURATION: City Center Hospital » All Departments » All Locations

Instrument

Configuration Download ☐ Download Location ☒ Assigned Location

Patient Id Mapping

☒ Enable Automatic Conversion

Control Test Mapping

Low(L1) Patient Id

Normal(L2) Patient Id

High(L3) Patient Id

Linearity Test Mapping

| | | | |
|---------------|---------------------------------|----------------|----------------------------------|
| L1 Patient Id | <input type="text" value="15"/> | L2 Patient Id | <input type="text" value="25"/> |
| L3 Patient Id | <input type="text" value="35"/> | L4 Patient Id | <input type="text" value="45"/> |
| L5 Patient Id | <input type="text" value="55"/> | L6 Patient Id | <input type="text" value="65"/> |
| L7 Patient Id | <input type="text" value="75"/> | L8 Patient Id | <input type="text" value="85"/> |
| L9 Patient Id | <input type="text" value="95"/> | L10 Patient Id | <input type="text" value="100"/> |

Proficiency Test Mapping

Patient Id

1

Instruments Comment Code Screens

Instruments>>Comments

The **Instrument Comment Codes** screens allow you to add, edit and delete user-defined comment codes.

- Each Facility in the Organization maintains a separate list of user-defined comment codes.
- The PCx has predefined instrument comment codes that cannot be edited.

To **view** instrument and user-defined comment codes:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Comments** option.

The **Instrument Comment Codes** screen displays the Comment Code table for the PCx with predefined instrument comment codes numbered 0 through 8 in read-only fields.

User-defined comment codes can be added, edited and deleted.

INSTRUMENT COMMENT CODES: F1

Search Criteria

Instrument

Display Records

9 Records Found

| Add | Edit | Delete | |
|------------------------------|-----------------------------|---------------|--|
| Code | Description | | |
| <input type="checkbox"/> All | | | |
| <input type="checkbox"/> 0 | No action required | | |
| <input type="checkbox"/> 1 | Repeat Test | | |
| <input type="checkbox"/> 2 | New Test Strip Box | | |
| <input type="checkbox"/> 3 | New Controls Bottle | | |
| <input type="checkbox"/> 4 | Procedure Error | | |
| <input type="checkbox"/> 5 | New Test Strip Lot | | |
| <input type="checkbox"/> 6 | New Controls Lot | | |
| <input type="checkbox"/> 7 | Changed Meters | | |
| <input type="checkbox"/> 8 | Changed Batteries | | |

Add Edit Delete

9 Records Found


To **view** comment codes at a different location:

Click on a different **Facility** in the Tree.

Add a comment code

To **add** a user-defined comment code:

1. Click the **Add** action button.
 - The **Add Instrument Comment Code** screen displays the instrument type that was selected from the dropdown list in the **Instrument Comment Codes** screen before clicking the **Add** button.



ADD INSTRUMENT COMMENT CODE: Southern Medical Center

Comment Details

Instrument Type

* Code

* Description

- Enter a Code number from 0 to 99,999 (PCx user defined comment codes start at 9) and a comment code description.
- The instrument comment code list is used to associate a description to a comment code number for reporting purposes only. The instrument comment code list does not download to the specified instrument.

* Code (number) and Description are **required** entry fields.

2. Click the **Save** button to save the data or Click the **Cancel** button to cancel the add action without saving.

Edit a comment code

To **edit** a single user-defined comment code or multiple user-defined comment codes:


Note: Default comment codes cannot be edited or deleted.

1. Click the box(es) in the **Code** column.
2. Click the **Edit** action button.

To **edit** all user-defined comment codes:

1. Click the box marked **All** in the **Code** column header.
2. Click the **Edit** action button.

The **Edit Instrument Comment Code** screen appears with the information stored in the database for the first selected user-defined comment code record.



* Code (number) and Description are **required** entry fields.

3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.


Delete a comment code

To **delete** a single user-defined comment code or multiple user-defined comment codes:

1. Click the box(es) in the **Code** column.
2. Click the **Delete** action button.

To **delete** all user-defined comment codes:


1. Click the box marked **All** in the code column header.
2. Click the **Delete** action button.
3. The following dialogue box appears:



4. Click **OK** to delete the record(s) or **Cancel** to return to the **Instrument Comment Codes** screen.

Warning Messages


1. If a different menu option is selected before saving the new or edited comment code record, the following warning message is displayed:



- Click the **OK** button to cancel the data entered.
 - Click the **Cancel** button to clear the message and return to the **Instrument Comment Codes** screen to save the data.
2. If a **required** field is missing when the comment code information is saved, the following warning message is displayed with a reference to the specific required field (Code and/or Description) that is missing:




- Click the **OK** button to clear the message and enter the required information.
3. If comment code records are not selected from the table before clicking on the **Edit** or **Delete** action button, the following warning message is displayed:




- Click the **OK** button to clear the message and make a selection.

4. If a Code is entered that is already in use, the following warning message is displayed:



- Click the **OK** button and choose a different Code that is not in use.
5. If you attempt to delete a comment with attached data, the following warning message is displayed:



Click the **OK** button to return to the Instrument Comment Code screen.

Instrument Port Assignment Screens

Instruments>>Ports

The **Instrument Port Assignment** screen allow you to view, add, edit and delete COM port assignments for instruments with a serial connection to the server.

Use the **Server Name** selection criterion to view COM port assignments for instruments with serial connections at the selected server.

To **view** instrument port assignments:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Ports** option.
3. Click on a Location in the Tree.
4. Click the **Display Records** button.

The **Instrument Port Assignment** screen displays the information stored in the database for instruments with serial connections to the server.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin LogOff

Refresh Hide Tree

INSTRUMENT PORT ASSIGNMENT:

Search Criteria
Server Name: GML TEST1
Display Records

1 Record Found

| Instrument | Serial Port | Port Speed | Facility | Department | Location | MeterCOM TCP port | MeterCOM Rev |
|------------|-------------|------------|------------------------|------------|----------|-------------------|--------------|
| PCX | COM1 | 19200 | Eastern Medical Center | Pediatrics | Ward A | 3000 | 1.0 |

1 Record Found

Add a Port Assignment

To **add** a COM port assignment for an instrument with a serial connection to the server:

1. Click the **Add** action button.
 - The **Add Instrument Port Assignment** screen displays the Server Name in a read-only field.
2. Select an instrument from the dropdown list.
3. Select an available COM Port from the dropdown list.
4. Select the baud rate from the **Port Speed** dropdown list.

* Port Speed, Facility, Department and Location fields are **required** fields.

5. Click on a **Location** in the Tree.
 - The selected Location and the associated Facility and Department fields automatically display the selected Location.
6. Click the **Save** button to save the data or Click the **Cancel** button to cancel the add action without saving.

Edit a Port Assignment

To **edit** a COM port assignment for an instrument at a specific location:

1. Click the box in the **Instrument** column.
2. Click the **Edit** action button.

The **Edit Instrument Port Assignment** screen appears with the information stored in the database for the selected COM port assignment record.

* The Port Speed, Facility, Department and Location fields are **required** fields.

To **change** the Instrument type, COM port assignment, the port speed or the location of the instrument:


1. Select a different Instrument from the dropdown list.
2. Select a different **COM port** from the dropdown list.
3. Click on a different Location in the Tree.

- The Location field and the associated Department and Facility fields automatically display the selected Location.
4. Select the baud rate from the **Port Speed** dropdown list.
 5. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

Delete a Port Assignment

To **delete** a COM port assignment for an instrument at a specific location:


1. Click the box in the **Instrument** column.
2. Click the **Delete** action button.
3. The following dialogue box is displayed.



4. Click the **OK** button to delete the COM port record(s) or Click the **Cancel** button to cancel the delete action without saving.


Warning Messages

If a different menu option is selected before saving a new or edited instrument port assignment record, the following warning message is displayed:




- Click the **OK** button to cancel the data entered.
- Click the **Cancel** button to clear the message and return to the add or edit **Port Assignment** screen to save the data.

If a **required** field is missing when the add or edit **Instrument Port Assignment** screens are saved, the following warning message is displayed with a reference to the specific field (Department and/or Location) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the COM port entered already exists in the database for the same location, the following error message is displayed:



- Click the **OK** button to clear the message.

If records are not selected from the table before clicking on the **Edit** or **Delete** action button, the following warning message is displayed:



- Click the **OK** button to clear the message and make a selection.

Tracking Log

Instruments>>Tracking Log

The **Instrument Tracking Log** screen allows the user to view events related to placing instruments in to service, moving an instruments assigned Location and taking instruments out of service.

- Enter an instrument serial number and/or name to view the log of events generated by a single instrument.
- Use the **Source** search criterion to view the log of user-generated or instrument-generated events, and use the **Date Range** search criterion to display events that occurred in a specific time period.

To **view** the instrument tracking log:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Tracking Log** option.

The **Instrument Tracking Log** screen displays the **Instrument Event** table based on the current Tree location, the system default Search Criteria and a date range of six months prior to today's date.

INSTRUMENT TRACKING LOG: F3 » D3 » L3 No:175 LIS:123 Op:9

Search Criteria

Instrument: PCx Serial No.: From: 08/15/2006 To: 02/15/2007

Source: All Name: Display Records

1 Record Found

| Date/Time | Instrument Type | Instrument Name | Serial No. | User Name | Comment |
|--------------------|-----------------|-----------------|------------|-----------|----------------|
| 01/10/2007 11:22AM | PCx | -- | M06210023 | -- | New Instrument |

1 Record Found

The available search criteria on the **Instrument Tracking Log** screen are:

| Criteria | Options | Description |
|-------------------|----------------------|---|
| Instrument | All | System displays tracking events for all instrument types. |
| | Selected Instrument | System displays results only from the selected Instrument type. |
| Source | All | System displays all tracking events. |
| | Instrument generated | System only displays instrument-generated tracking events. |
| | User Entry | System only displays user-generated tracking events. |
| Serial No. | | Entry field for an instrument serial number. |

| Criteria | Options | Description |
|-------------------|-------------|--|
| Name | | Entry field for an instrument name. System only displays tracking events for the instrument that matches, or is similar to, the serial number or name entered. |
| Date Range | | System default displays date range of six months prior to today's date. Enter or select a From and a To date to display tracking events that occurred between the selected dates. |
| | From | Click on the calendar icon to display the pop-up calendar and click on a date. The From date field automatically displays the selected date. |
| | To | Click on the calendar icon to display the pop-up calendar and click on a date. The To date field automatically displays the selected date. |

Add a User-defined Event:

To add a user-defined event:

1. Click the **Add** action button.

The **Add Instrument Tracking Log Entry** screen displays the instrument in a read-only field and the entry fields for an instrument Serial No. and a user-generated event comment.

ADD INSTRUMENT TRACKING LOG ENTRY: Eastern Medical Center

Log Details

Instrument: PCx

* Serial No. [Text Field]

* Comment [Text Field]

[Save] [Cancel]

* Serial No. and Comment are **required** entry fields.

2. Click the **Save** button to save the data or Click the **Cancel** button to cancel the add action without saving.


Edit a User-defined Event

To **edit** a single user-defined event or multiple user-defined events:

1. Click the box(es) in the Date/Time column.

2. Click the **Edit** action button.

The **Edit Instrument Tracking Log Entry** screen displays the Serial No. in a read-only field and the user-generated comment stored in the database for the first selected record.




* Serial No. and Comment are **required** entry fields.

3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action, then Click the **OK** button to clear the message.

Delete a User-defined Event

To **delete** a single user-defined event or multiple user-defined events:

1. Click the box in the **Date/Time** column.
2. Click the **Delete** action button.
3. The following dialogue box is displayed.



4. Click the **OK** button to delete the user-defined event(s) or Click the **Cancel** button to cancel the delete action without saving.

Note: If the serial number specified in the search criteria entry field of the Instrument Tracking Log screen is not in the database, the Table displays a message: "No data is available for this selection. Please change criteria or tree selection."


Upload Log

Instruments>>Upload Log

To **view** the instrument upload log:

1. Move the mouse cursor over **Instruments** to display the menu options.
2. Click on the **Upload Log** option.

The **Upload Log** screen displays information on instrument upload and data forwarding events.



8. Lots

Lots Overview

The **Lots** menu option provides access to five screens for use in managing reagent lots. The **Control** screen allows you to manage Quality Control lot information. The **Linearity** screen allows you to manage lot information for Linearity Panels. The **Proficiency** screen allows you to manage lot information for Proficiency Testing samples. The **Reagent >> Single & Multi Analyte** screens allow you to manage lot information for Test Strips and Cartridges. Each screen is available only if your account was granted this permission by the Administrator or POCC.

Control Lots Screen

Lots>>Control

The **Control Lots** screen allows you to view, add, edit and delete Quality Control lot information for a Facility.

- Each Facility in the Organization maintains a separate control lot list.
- Use the **Search Criteria** options to display lots that were put into service within a specific date range.

To **view** control lots:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the **Control** option.

Note: If the Organization is highlighted in the Tree when the **Lots** menu option is selected, no results are displayed and the **Add** action button is not displayed on the table's action bar.

The **Control Lots** screen displays the Control Lots table with lots based on the Facility highlighted in the Tree.

To **view** lots at a different location:

1. Click on a **Facility** in the Tree.
 - All other Tree options are disabled.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin: Log Off

CONTROL LOTS: Eastern Medical Center No:134 LIS:58 Op:29

Search Criteria
In-Service Date From In-Service Date To
Display Records

3 Records Found

| Lot No. | Manufacturer | Reagent Type | Level | In-Service Date | Expiration Date |
|--------------------------------|--------------|--------------------|-------|-----------------|-----------------|
| <input type="checkbox"/> All | | | | | |
| <input type="checkbox"/> 13547 | -- | Glucose Test Strip | High | -- | -- |
| <input type="checkbox"/> 13857 | -- | Glucose Test Strip | High | -- | -- |
| <input type="checkbox"/> 13552 | -- | Glucose Test Strip | Low | -- | -- |

Add Edit Delete Customize

3 Records Found

Add a Control Lot

To **add** a control lot:

1. Click on a Facility in the Tree.
2. Click the **Add** action button.

The **Add Control Lot** screen displays the entry fields for control lot information.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin: Log Off

ADD CONTROL LOT: Eastern Medical Center No:134 LIS:58 Op:29

Input Information

* Lot No.

Manufacturer

Reagent Type

In-Service Date

* Expiration Date

Level

Save Cancel

* Lot No. and Expiration Date are **required** entry fields. The Lot No. field accepts up to a **maximum** of 30 characters.

The information entered will be stored in the database.

| | |
|-------------------------|--|
| *Lot No. | Enter up to 30 characters for the lot number – required . |
| Manufacturer | Enter up to 30 characters for the manufacturer name. |
| Reagent Type | Select a Reagent Type from the dropdown list. The system default for Reagent Type is Glucose Test Strip. |
| In-Service Date | Click on the calendar icon to display the pop-up calendar and click on a date. The In-Service Date field automatically displays the selected date. |
| *Expiration Date | Click on the calendar icon to display the pop-up calendar and click on a date. |

The **Expiration Date** field automatically displays the selected date – **required**.

Level

Select the control solution level from the dropdown list:

Low
Normal
High

1. Click the **Save** button to save the data. The Control Lots screen appears with the new lot displayed in the Control Lots table.
 - Click the **Cancel** button to cancel the add action without saving.

Edit a Control Lot

To **edit** a single control lot or multiple control lots in a Facility:

1. Click the box(es) in the Lot No. column.
2. Click the **Edit** action button.

To **edit** all control lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Edit** action button.

The **Edit Control Lot** screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

The screenshot shows the 'EDIT CONTROL LOT' interface for Eastern Medical Center. The top navigation bar includes links for Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. A user status bar on the right shows 'admin: Log Off' and counts for 'No: 134', 'LIS: 58', and 'Op: 29'. On the left, a tree view shows the facility hierarchy: Medical Center Inc, Eastern Medical Center (selected), Pediatrics, Ward A, Ward B, Southern Medical Center, and Western Medical Center. The main form area is titled 'EDIT CONTROL LOT: Eastern Medical Center' and contains the following fields:

- * Lot No.: 13547 (read-only)
- Manufacturer: (empty text field)
- Reagent Type: Glucose Test Strip (dropdown menu)
- In-Service Date: 3/1 (calendar icon)
- * Expiration Date: 3/1 (calendar icon)
- Level: High (dropdown menu)

 At the bottom of the form, there are navigation buttons '<< 1 of 3 >>' and action buttons 'Save' and 'Cancel'.

* Expiration Date is a **required** entry field.


The Lot No. is displayed in a read-only field and cannot be edited. All other entry fields can be edited and will be stored in the database.

3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

Delete a Control Lot

To **delete** a single control lot or multiple control lots in a Facility:


1. Click the box(es) in the Lot No. column.
2. Click the **Delete** action button.
3. The following message will be displayed:



4. Click **OK** to delete the record.
5. Click **Cancel** to return to the Control Lots screen.

To **delete** all control lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Delete** action button.
3. The following message will be displayed:




4. Click **OK** to delete the selected records.
5. Click **Cancel** to return to the Control Lots screen.

Note: The delete action will remove the control lot only if **no** test data exists in the active database for that lot. If the selected control lot is associated with test data, the control lot cannot be deleted.

Using the Filter Options

To **find** control lots by In-Service Date:

1. Click on the calendar icon in the **Search Criteria** options to display the pop-up calendar.




2. Click on a date for the **From** and/or **To** date fields.
3. Click the **Display Records** button.

The **Control Lots** screen displays the Control Lots table based on results that match the dates entered in the From and To date fields.

Warning Messages

If a different menu option is selected before saving new or edited lot information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the add or edit screen to save the data.

If a **required** field is missing when lot information is saved, the following warning message is displayed with a reference to the specific required field (Lot No. and/or Expiration Date) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the Lot Number entered already exists in the database, the following error message is displayed:




- Click the **OK** button to clear the message.

If data entered in the Expiration Date field is not valid, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the valid information.

If records are not selected from the table before clicking the **Edit** or **Delete** action buttons, the following warning message is displayed.



- Click the **OK** button to clear the message and make a selection.

If a lot is selected for deletion with test data stored in the active database, the following error message is displayed.



- Click the **OK** button to clear the message.

Linearity Lots Screen

Lots>>Linearity

The **Linearity Lots** screen allows you to view, add, edit and delete linearity lot information for a Facility.

- Each Facility in the Organization maintains a separate linearity lot list.
- Use the **Search Criteria** options to change the Reagent Type or display lots that were put into service within a specified date range.

To **view** linearity lots:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the **Linearity** option.

The **Linearity Lots** screen displays the Linearity Lots table with lots based on the Facility highlighted in the Tree.

The default Reagent Type is Glucose Test Strip. Select the Reagent Type dropdown list to select a different reagent type.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Logout

LINEARITY LOTS: Western Medical Center

Search Criteria

Reagent Type: Glucose Test Strip
 In-Service Date From: 31
 In-Service Date To: 31
 Display Records

4 Records Found

Add Edit Delete Customize

| Lot No. | Reagent Type | Expiration Date | In-Service Date |
|------------------------------|--------------------|-----------------|-----------------|
| <input type="checkbox"/> All | | | |
| <input type="checkbox"/> 777 | Glucose Test Strip | 09/30/2004 | -- |
| <input type="checkbox"/> 789 | Glucose Test Strip | 09/30/2004 | -- |
| <input type="checkbox"/> 123 | Glucose Test Strip | 10/29/2004 | -- |

Add Edit Delete Customize

4 Records Found

Note: If the Organization is highlighted in the Tree when the **Lots** menu option is selected, then no results are displayed and the **Add** action button is not displayed on the table's action bar.

To **view** lots at a different location:

1. Click on a **Facility** in the Tree.
 - All other Tree options are disabled.

Add a Linearity Lot

To **add** a linearity lot click on a Facility in the Tree.

1. Click on a Facility in the Tree.
2. Select a reagent type from the Reagent Type dropdown list in the Search Criteria and click Display Records.
3. Click the **Add** action button and enter the linearity lot information

The **Add Linearity Lot** screen displays the entry fields for linearity lot information.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

ADD LINEARITY LOT: Western Medical Center No:134 LIS:58 Op:29

Input Information

Reagent Type: Glucose Test Strip

* Lot No.:

Manufacturer:

In-Service Date: 3/1

* Expiration Date: 3/1

Save Cancel

1 Record Found

| Analyte | L1 | L2 | L3 | L4 | L5 | L6 | L7 | L8 | L9 | L10 |
|---------|----|----|----|----|----|----|----|----|----|-----|
| GLU | -- | -- | -- | -- | -- | -- | -- | -- | -- | -- |

1 Record Found

* Lot No. and Expiration Date are **required** entry fields. The Lot No. field accepts up to a **maximum** of 30 characters.

The information entered will be stored in the database.

| Field | Description |
|-------------------------|--|
| *Lot No. | Enter up to 30 characters for the lot number – required . |
| Manufacturer | Enter up to 30 characters for the manufacturer name. |
| Reagent Type | System default for Reagent Type is Glucose Test Strip. A different reagent type may be specified from the Reagent Type dropdown list on the Linearity Lots screen. |
| In-Service Date | Click on the calendar icon to display the pop-up calendar and click on a date. The In-Service Date field automatically displays the selected date. |
| *Expiration Date | Click on the calendar icon to display the pop-up calendar and click on a date. |

| Field | Description |
|--------------|--|
| | The Expiration Date field automatically displays the selected date – required . |
| Level | Linearity levels are entered on a separate sub-screen. To access it, select the linked Analyte column header from the table in the bottom part of the screen and click the Edit button (see below). |

- Click the **Save** button to save the data. The Linearity Lots screen appears with the new lot displayed in the Linearity Lots table or Click the **Cancel** button to cancel the add action without saving.

To **add** Analyte Ranges (from the **Add Linearity Lot** screen):

- Select one or more records from the table by click the box(es) in the **Analyte** column.
- Click the **Edit** action button.

The **Edit Analyte Ranges** screen appears with information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

EDIT ANALYTE RANGES: F1
Input Information
Reagent Type
* Lot No.
Manufacturer
In-Service Date
* Expiration Date
Analyte

| | | | |
|-----|----------------------|----|----------------------|
| L1 | <input type="text"/> | To | <input type="text"/> |
| L2 | <input type="text"/> | To | <input type="text"/> |
| L3 | <input type="text"/> | To | <input type="text"/> |
| L4 | <input type="text"/> | To | <input type="text"/> |
| L5 | <input type="text"/> | To | <input type="text"/> |
| L6 | <input type="text"/> | To | <input type="text"/> |
| L7 | <input type="text"/> | To | <input type="text"/> |
| L8 | <input type="text"/> | To | <input type="text"/> |
| L9 | <input type="text"/> | To | <input type="text"/> |
| L10 | <input type="text"/> | To | <input type="text"/> |

<< 1 of 1 >>
Save Cancel

All of the fields are disabled except for the Range fields.

3. Edit the range fields (L1 - L10) as applicable.
4. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

Edit a Linearity Lot

To **edit** a single linearity lot or multiple linearity lots in a Facility:

1. Click the box(es) in the Lot No. column.
2. Click the **Edit** action button.

To **edit** all linearity lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Edit** action button.

The **Edit Linearity Lot** screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

EDIT LINEARITY LOT: Western Medical Center No: 134 LIS: 58 Op: 29

Input Information

Reagent Type: Glucose Test Strip

* Lot No.: 1

Manufacturer:

In-Service Date: 3/1

* Expiration Date: 09/30/2004

<< 1 of 4 >>

Save Cancel

1 Record Found

| Edit | Analyte | L1 | L2 | L3 | L4 | L5 | L6 | L7 | L8 | L9 | L10 | Customize |
|------------------------------|---------|----|----|----|----|----|----|----|----|----|-----|-----------|
| <input type="checkbox"/> All | | | | | | | | | | | | |
| <input type="checkbox"/> GLU | -- | -- | -- | -- | -- | -- | -- | -- | -- | -- | -- | |
| Edit | | | | | | | | | | | | Customize |

1 Record Found


* Expiration Date is a **required** entry field.

- The **Lot No.** is displayed in a read-only field and cannot be edited. All other entry fields can be edited and will be stored in the database (to edit levels, see Section To Edit Analyte Ranges, above).
3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

Delete a Linearity Lot

To **delete** a single linearity lot or multiple linearity lots in a Facility:


1. Click the box(es) in the Lot No. column.
2. Click the **Delete** action button.
3. The following message will be displayed:



4. Click **OK** to delete the record or Click **Cancel** to return to the Linearity Lots screen without saving.

To **delete** all linearity lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Delete** action button.
3. The following message will be displayed:



4. Click **OK** to delete the records or Click **Cancel** to return to the Linearity Lots screen without saving.

Note: The delete action will remove the linearity lot only if **no** test data exists in the active database for that lot. If the selected linearity lot is associated with test data, the linearity lot cannot be deleted.

Using the Filter Options

To **find** linearity lots by Reagent Type:

1. Select from the dropdown list in the Reagent Type field.
2. Click the **Display Records** button.

The **Linearity Lots** screen displays the Linearity Lots table based on results that match the Reagent Type chosen from the dropdown list.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Logout

LINEARITY LOTS: Western Medical Center No: 154 LIS: 50 Op: 20

Search Criteria

Reagent Type: Glucose Test Strip

In-Service Date From: 31 In-Service Date To: 31

Display Records

4 Records Found

| Lot No. | Reagent Type | Expiration Date | In-Service Date |
|---------|--------------------|-----------------|-----------------|
| 1 | Glucose Test Strip | 09/30/2004 | -- |
| 777 | Glucose Test Strip | 09/30/2004 | -- |
| 789 | Glucose Test Strip | 09/30/2004 | -- |
| 123 | Glucose Test Strip | 10/29/2004 | -- |

4 Records Found


To **find** linearity lots by In-Service Date:

1. Click on the calendar icon in the **Search Criteria** options to display the pop-up calendar.
2. Click on a date for the **From** and/or **To** date fields.
3. Click the **Display Records** button.

The **Linearity Lots** screen displays the Linearity Lots table based on results that match the dates entered in the **From** and **To** date fields.

Warning Messages

If a different menu option is selected before saving new or edited lot information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the add or edit screen to save the data.

If a **required** field is missing when lot information is saved, the following warning message is displayed with a reference to the specific required field (Lot Code, Lot No. and/or Expiration Date) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the Lot Number entered already exists in the database, the following error message is displayed:




- Click the **OK** button to clear the message.

If data entered in the Expiration Date field is not valid, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the valid information.

If records are not selected from the table before clicking the **Edit** or **Delete** action buttons, the following warning message is displayed:



- Click the **OK** button to clear the message and make a selection.

If a lot is selected for deletion with test data stored in the active database, the following error message is displayed:



- Click the **OK** button to clear the message.

Proficiency Lots Screen

Lots>>Proficiency

The **Proficiency Lots** screen allows you to view, add, edit and delete proficiency lot information for a Facility.

- Each Facility in the Organization maintains a separate proficiency lot list.
- Use the **Search Criteria** options to change the Reagent Type or display lots that were put into service within a specified date range.

To **view** proficiency lots:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the **Proficiency** option.

The **Proficiency Lots** screen displays the Proficiency Lots table with the lots based on the Facility highlighted in the Tree.

The default Reagent Type is Glucose Test Strip. Select the Reagent Type dropdown list to select a different reagent type.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

Refresh Hide Tree x

PROFICIENCY LOTS:

Search Criteria

Reagent Type: Glucose Test Strip In-Service Date From: 31 In-Service Date To: 31

Display Records

0 Records Found

| Lot No. | Manufacturer | Reagent Type | Expiration Date | In-Service Date | Test ID |
|--|--------------|--------------|-----------------|-----------------|---------|
| No data is available for this selection. Please change selection criteria or tree selection. | | | | | |

Edit Delete Customize

0 Records Found

Note: If the Organization is highlighted in the Tree when the **Lots** menu option is selected, then no results are displayed and the **Add** action button is not displayed on the table's action bar.

To **view** lots at a different location:

1. Click on a **Facility** in the Tree.
 - All other Tree options are disabled.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin Log Off

PROFICIENCY LOTS: Western Medical Center No:134 LIS:58 Op:29

Search Criteria

Reagent Type: Glucose Test Strip In-Service Date From: 31 In-Service Date To: 31

Display Records

2 Records Found

| Lot No. | Manufacturer | Reagent Type | Expiration Date | In-Service Date | Test ID |
|---------|--------------|--------------------|-----------------|-----------------|---------|
| 123 | -- | Glucose Test Strip | 09/27/2004 | -- | 345 |
| 789 | -- | Glucose Test Strip | 09/27/2004 | -- | 098 |

2 Records Found

Add a Proficiency Lot

1. Click on a Facility in the Tree.
2. Select a reagent type from the Reagent Type dropdown list in the Search Criteria and click Display Records.
3. Click the **Add** action button to add the proficiency lot information.

The **Add Proficiency Lot** screen displays the entry fields for proficiency lot information.


ADD PROFICIENCY LOT: Western Medical Center


Input Information

Reagent Type

* Lot No.

Manufacturer

In-Service Date 

* Expiration Date 

* Test ID

1 Record Found

| <input type="button" value="Edit"/> | <input type="button" value="Customize"/> |
|-------------------------------------|--|
| Analyte | Target Range |
| <input type="checkbox"/> All | |
| <input type="checkbox"/> GLU | -- |
| <input type="button" value="Edit"/> | <input type="button" value="Customize"/> |

1 Record Found

* Lot No., Expiration Date and Test ID are **required** entry fields. The Lot No. field accepts up to a **maximum** of 30 characters.

The information entered will be stored in the database.

| Field | Description |
|-------------------------|---|
| *Lot No. | Enter up to 30 characters for the lot number – required . |
| Manufacturer | Enter up to 30 characters for the manufacturer name. |
| Reagent Type | System default for Reagent Type is Glucose Test Strip. A different reagent type may be specified from the Reagent Type dropdown list on the Linearity Lots screen. |
| In-Service Date | Click on the calendar icon to display the pop-up calendar and click on a date. The In-Service Date field automatically displays the selected date. |
| *Expiration Date | Click on the calendar icon to display the pop-up calendar and click on a date. The Expiration Date field automatically displays the selected date – required . |

| Field | Description |
|-----------------|---|
| *Test ID | Enter test identifier. |
| Range | The target proficiency range is entered and edited on a separate subscreen, accessible from the Add/Edit Proficiency Lot screen. |

- Click the **Save** button to save the data. The Proficiency Lots screen appears with the new lot displayed in the Proficiency Lots table or Click the **Cancel** button to cancel the add action without saving.

To **add** Analyte Ranges (from the **Add Proficiency Lot** screen):

- Select one or more records from the Table by clicking the box(es) in the **Analyte** column.
- Click the **Edit** action button.

The Edit Analyte Ranges screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

All of the fields are disabled except for the Range fields.

EDIT ANALYTE RANGES: F3

Input Information

Reagent Type:

* Lot No.:

Manufacturer:

In-Service Date:

* Expiration Date:

* Test ID:

Analyte:

Range: To:

<< 1 of 1 >>

Save Cancel

- Edit the Range fields.
- Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

Edit a Proficiency Lot

To **edit** a single proficiency lot or multiple proficiency lots in a Facility:

1. Click the box(es) in the Lot No. column.
2. Click the **Edit** action button.

To **edit** all proficiency lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Edit** action button.

* Expiration Date and Test ID are **required** entry fields.

The **Edit Proficiency Lot** screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

The Lot No. is displayed in a read-only field and cannot be edited. All other entry fields can be edited and will be stored in the database (to edit Target Ranges, see Section To **Edit** Analyte Ranges, above).

3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

EDIT PROFICIENCY LOT: Western Medical Center

Input Information

Reagent Type:

* Lot No.:

Manufacturer:

In-Service Date:

* Expiration Date:

* Test ID:

<< 1 of 3 >>

Save Cancel

1 Record Found

| Analyte | Target Range |
|------------------------------|--------------|
| <input type="checkbox"/> All | |
| <input type="checkbox"/> GLU | -- |

1 Record Found

To **Edit** Analyte Ranges (from the **Edit Proficiency Lot** screen):


1. Select one or more records from the Table by clicking the box(es) in the **Analyte** column.

2. Click the **Edit** action button.

The Edit Analyte Ranges screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected, each record is displayed on a separate page.

All of the fields are disabled except for the Range fields.

3. Fill in the Range fields.
4. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.



EDIT ANALYTE RANGES: Test Facility1

Input Information

Reagent Type

* Lot No.

Manufacturer

In-Service Date

* Expiration Date

* Test ID

Analyte

Range To

<< 1 of 1 >>

Save Cancel

Delete a Proficiency Lot

To **delete** a single proficiency lot or multiple proficiency lots in a Facility:

1. Click the box(es) in the Lot No. column.
2. Click the **Delete** action button.

To **delete** all proficiency lots in a Facility:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Delete** action button.

Note: The delete action will remove the proficiency lot only if **no** test data exists in the active database for that lot. If the selected proficiency lot is associated with test data, the proficiency lot cannot be deleted.

Using the Filter Options

To **find** proficiency lots by Reagent Type:

1. Select from the dropdown list in the Reagent Type field.
2. Click the **Display Records** button.

The **Proficiency Lots** screen displays the Proficiency Lots table based on results that match the Reagent Type chosen from the dropdown list.


To **find** proficiency lots by In-Service Date:

1. Click on the calendar icon in the **Search Criteria** Options to display the pop-up calendar.
2. Click on a date for the **From** and/or **To** date fields.
3. Click the **Display Records** button.

The **Proficiency Lots** screen displays the Proficiency Lots table based on results that match the dates entered in the **From** and **To** date fields.


Warning Messages

If a different menu option is selected before saving new or edited lot information, the following warning message is displayed.




- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the add or edit screen to save the data.

If a **required** field is missing when lot information is saved, the following warning message is displayed with a reference to the specific required field (Lot No., Expiration Date and/or Test ID) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the Lot Number entered already exists in the database, the following error message is displayed:




- Click the **OK** button to clear the message.

If data entered in the Expiration Date field is not valid, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the valid information.

If records are not selected from the table before clicking the **Edit** or **Delete** action buttons, the following warning message is displayed:



- Click the **OK** button to clear the message and make a selection.

If a lot is selected for deletion with test data stored in the active database, the following error message is displayed:



- Click the **OK** button to clear the message.

Reagent Lots Screen

Lots>>Reagent

The **Reagent Lots** screen allows you to view, add, edit and delete reagent lots assigned to a Facility, Department or Location.

- Each Facility, Department or Location maintains a separate reagent lot list.
- Use the **Search Criteria** options to display lots that were put into service within a specific date range or view lots by their download to instrument status.

Reagent Lots and the Parent-Child Relationship

Reagent Lots saved to the Download to Instruments list at the top level of the Tree affect all sub-levels. For example, a reagent lot saved to the Facility (parent) level will be applied to all associated Departments and Locations.

Reagent Lots saved to the Download to Instruments list at a lower level of the Tree affect only the levels below the selected level of the Tree. For example, a reagent lot saved at the Department level will be applied to all Locations (child) associated to the Department (parent). Reagent Lots saved at the Location level affect only the selected Location.

It is important to note the selected Tree level prior to saving Reagent Lots to ensure the Download to Instruments list is applied to the appropriate locations.

Unlike with Instrument Configuration settings, there is no way to reestablish the parent-child relationship once a reagent lot has been saved to the Department or Location level of the Tree. Once a location is disassociated from its parent location any changes to reagent lots will need to be made independent of the parent location.

To **view** reagent lots:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the **Reagent** option.
3. Choose **Single Analyte** or **Multi-Analyte**

The **Reagent Lots** screen displays the Reagent Lots table with lots based on the Location, Department or Facility highlighted in the Tree.

REAGENT LOTS: F1 » All Departments » All Locations

Search Criteria: Reagent Type: Glucose Test Strip, Download: All, In-Service: , In-Service Date To:

16 Records Found

| Lot No. | Lot Code | Manufacturer | In-Service Date | Expiration Date | Download | Low Range | Normal Range | High Range |
|---------|----------|--------------|-----------------|-----------------|----------|-----------|--------------|------------|
| -- | 156245 | -- | -- | -- | No | -- | -- | -- |
| -- | 65jk2e | -- | -- | -- | No | -- | -- | -- |
| -- | 188755 | -- | -- | -- | No | -- | -- | -- |
| -- | 65sk1e | -- | -- | -- | No | -- | -- | -- |
| -- | 65dk1e | -- | -- | -- | No | -- | -- | -- |
| -- | 65hk3e | -- | -- | -- | No | -- | -- | -- |
| -- | 66fk3e | -- | -- | -- | No | -- | -- | -- |
| -- | 166915 | -- | -- | -- | No | -- | -- | -- |
| -- | 67ek7e | -- | -- | -- | No | -- | -- | -- |
| -- | 66gl3e | -- | -- | -- | No | -- | -- | -- |
| -- | 167725 | -- | -- | -- | No | -- | -- | -- |
| -- | 158595 | -- | -- | -- | No | -- | -- | -- |
| -- | 168555 | -- | -- | -- | No | -- | -- | -- |
| -- | 168045 | -- | -- | -- | No | -- | -- | -- |
| -- | 66kl5e | -- | -- | -- | No | -- | -- | -- |
| -- | 166245 | -- | -- | -- | No | -- | -- | -- |

Note: If the Organization is highlighted in the Tree when the **Lots** menu option is selected, no results are displayed, and the **Add** action button is not displayed on the table's action bar.

To **view** lots at a different location:

1. Click on a **Facility, Department** or **Location** in the Tree.

Add a Single Analyte Reagent Lot

To **add** a single analyte reagent lot:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the Reagent option.
3. Choose Single Analyte
4. Click on a Facility, Department or Location in the Tree.
5. Select Reagent Type Select a reagent type from the Reagent Type dropdown list in the Search Criteria and click Display Records.
6. Click the Add action button.

* Lot Code, Lot No, Expiration Date and the values for the Low, Normal and High Ranges are **required** entry fields. The Lot Code field accepts up to a **maximum** of 30 characters.

- The **Add Reagent Lot** screen displays the Reagent Type in a read-only field and the entry fields for reagent lot information.
- The information entered will be stored in the database.

| Field | Description |
|-------------------------|--|
| *Lot Code | Enter the first six (6) characters from the Abbott glucose test strip reagent lot barcode – required. |
| *Lot No. | The application automatically displays the six (6) characters entered in the Lot Code field. The additional eight (8) characters from the Abbott glucose test strip reagent lot barcode are required to download the reagent lot to glucose instruments – required. |
| Manufacturer | Enter up to 30 characters for the manufacturer name. |
| In-Service Date | Click on the calendar icon to display the pop-up calendar and click on a date. The In-Service Date field automatically displays the selected date. |
| *Expiration Date | Click on the calendar icon to display the pop-up calendar and click on a date. The Expiration Date field automatically displays the selected date – required. |

- Click on the checkbox to allow the Reagent Lot to be downloaded to instruments.

☐ **Download To Instruments**

Note: The application does not restrict the number of reagent lots selected for Download To Instruments. Be sure to download only reagent lots currently in use to ensure that instruments receive the appropriate reagent download list. Any reagent lot not currently in use should be removed from the download list.

- Enter the QC ranges for the Reagent Lot:

***Low Range** Enter the lower and upper values for the Low Range – **required.**

***Normal Range** Enter the lower and upper values for the Normal Range – **required.**

***High Range** Enter the lower and upper values for the High Range – **required.**

- Click the **Save** button to save the data. The Reagent Lots screen appears with the new lot displayed in the Reagent Lots table or Click the **Cancel** button to cancel the add action without saving.

ADD REAGENT LOT: F1 » All Departments » All Locations

Input Information

Reagent Type

* Lot Code

* Lot No.

Manufacturer

In-Service Date

* Expiration Date

☐ Download to Instruments

* Low Range To

* Normal Range To

* High Range To

To **add** a multi analyte reagent lot:

1. Move the mouse cursor over **Lots** to display the menu options.
2. Click on the **Reagent** option.
3. Choose **Multi Analyte**
4. Click on a Facility, Department or Location in the Tree.
5. Select a reagent type from the Reagent Type dropdown list in the Search Criteria and click Display Records.
6. Click the **Add** action button.

The **Add Reagent Lot** screen displays the Reagent Type in a read-only field and the entry fields for reagent lot information.

The information entered will be stored in the database.

* Lot No and Expiration Date entry fields are **required**. The Lot No. field accepts up to a **maximum** of 30 characters.

| Field | Description |
|---------------------|--|
| *Lot No. | Enter up to 30 characters for the reagent lot number – required . |
| Manufacturer | Enter up to 30 characters for the manufacturer name. |

| Field | Description |
|-------------------------|---|
| In-Service Date | Click on the calendar icon to display the pop-up calendar and click on a date. The In-Service Date field automatically displays the selected date. |
| *Expiration Date | Click on the calendar icon to display the pop-up calendar and click on a date. The Expiration Date field automatically displays the selected date – required . |

7. Click on the checkbox to allow the Reagent Lot to be downloaded to instruments (if applicable).

☐ **Download To Instruments**

Note: The application does not restrict the number of reagent lots selected for Download To Instruments. Be sure to download only reagent lots currently in use to ensure that instruments receive the appropriate reagent download list. Any reagent lot not currently in use should be removed from the download list.

8. Click the **Save** button to save the data. The Reagent Lots screen appears with the new lot displayed in the Reagent Lots table or Click the **Cancel** button to cancel the add action without saving.

ADD REAGENT LOT: Test Facility1 » All Departments » All Locations

Input Information

Reagent Type:

* Lot No.:

Manufacturer:

In-Service Date:

* Expiration Date:

☐ Download to Instruments

14 Records Found

| Analyte | Low Range | Normal Range | High Range |
|--------------------------------|-----------|--------------|------------|
| <input type="checkbox"/> All | | | |
| <input type="checkbox"/> AnGap | -- | -- | -- |
| <input type="checkbox"/> BE | -- | -- | -- |
| <input type="checkbox"/> BUN | -- | -- | -- |
| <input type="checkbox"/> Cl | -- | -- | -- |
| <input type="checkbox"/> GLU | -- | -- | -- |
| <input type="checkbox"/> Hb | -- | -- | -- |
| <input type="checkbox"/> HCO3 | -- | -- | -- |
| <input type="checkbox"/> Hct | -- | -- | -- |
| <input type="checkbox"/> K | -- | -- | -- |
| <input type="checkbox"/> Na | -- | -- | -- |
| <input type="checkbox"/> PCO2 | -- | -- | -- |
| <input type="checkbox"/> pH | -- | -- | -- |
| <input type="checkbox"/> TCO2 | -- | -- | -- |
| <input type="checkbox"/> Urea | -- | -- | -- |

14 Records Found

To add Analyte Ranges (from the **Add Reagent Lot** screen):

1. Select one or more records from the Table by clicking the box(es) in the **Analyte** column.
2. Click the **Edit** action button.

The Edit Analyte Ranges screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

All of the fields are disabled except for the Range fields.

3. Edit the Range fields.
4. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

EDIT ANALYTE RANGES: F1 » All Departments » All Locations

Input Information

Reagent Type

Lot No.

Analyte

* Low Range To

* Normal Range To

* High Range To

<< 1 of 1 >>

Save Cancel

Edit a Reagent Lot

To **edit** a single reagent lot or multiple reagent lots in a Facility, Department or Location:

1. Click the box(es) in the Lot No. column.
2. Click the **Edit** action button.

To **edit** all reagent lots in a Facility, Department or Location:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Edit** action button.

The **Edit Reagent Lot** screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected each record is displayed on a separate page.

All other entry fields can be edited and will be stored in the database.

EDIT REAGENT LOT: Western Medical Center » All Departments » All Locations

Input Information

Reagent Type

* Lot Code

* Lot No.

Manufacturer

In-Service Date

* Expiration Date

☐ Download to Instruments

* Low Range To

* Normal Range To

* High Range To

<< 1 of 7 >>

Save Cancel

* Reagent Type and Lot Code are displayed in read-only fields and cannot be edited. The Expiration Date and the values for the Low, Normal and High Ranges are **required** entry fields.

3. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

To **Edit** Multi Analyte Ranges (from the **Edit Reagent Lot** screen):

1. Select one or more records from the Table by clicking the box(es) in the **Analyte** column.
2. Click the **Edit** action button.

The Edit Analyte Ranges screen appears with the information stored in the database for the first selected lot record. If multiple lot records are selected, each record is displayed on a separate page.

All of the fields are disabled except for the Range fields.

3. Fill in the Range fields.
4. Click the **Save** button to save the edited data or Click the **Cancel** button to cancel the edit action without saving.

EDIT ANALYTE RANGES: Test Facility1 » All Departments » All Locations

Input Information

Reagent Type

Lot No.

Analyte

* Low Range To

* Normal Range To

* High Range To

<< 1 of 1 >>

Save Cancel

Delete a Reagent Lot

To **delete** a single reagent lot or multiple reagent lots in a Facility, Department or Location:

1. Click the box(es) in the Lot No. column.
2. Click the **Delete** action button.

To **delete** all reagent lots in a Facility, Department or Location:

1. Click the box marked **All** in the Lot No. column header.
2. Click the **Delete** action button.

Note: The delete action will remove the reagent lot only if **no** test data exists in the active database for that lot. If the selected reagent lot is associated with test data, the reagent lot cannot be deleted.

Using the Filter Options

To **find** reagent lots by Reagent Type:

1. Select from the dropdown list in the Reagent Type field.
2. Click the **Display Records** button.

The **Reagent Lots** screen displays the Reagent Lots table based on results that match the Reagent Type chosen from the dropdown list.

PrecisionWeb Dashboard | Review | Reports | Operators | Instruments | Lots | Administrative | Help | admin | Log Off

REAGENT LOTS: Western Medical Center » All Departments » All Locations | No: 134 | LIS: 58 | Up: 29

Search Criteria

Reagent Type: **Glucose Test Strip** (dropdown menu open showing: Glucose Test Strip, ACT-4c, ACT-C, Crea, cTri)

In-Service Date From: [31] | In-Service Date To: [31]

Display Records

7 Records Found

| Lot No. | Lot Code | Manufacturer | In-Service Date | Expiration Date | Download | Low Range | Normal Range | High Range |
|------------------------------|----------|--------------|-----------------|-----------------|----------|-----------|--------------|------------|
| <input type="checkbox"/> All | | | | | | | | |
| <input type="checkbox"/> -- | 625123 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 625124 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 625125 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 625126 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 625127 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 149735 | -- | -- | -- | No | -- | -- | -- |
| <input type="checkbox"/> -- | 155255 | -- | -- | -- | No | -- | -- | -- |

7 Records Found

To **find** reagent lots by In-Service Date:

1. Click on the calendar icon to display the pop-up calendar.
2. Click on a date for the **From** and/or **To** date fields.
3. Click the **Display Records** button.

The **Reagent Lots** screen displays the Reagent Lots table based on results that match the dates entered in the From and To date fields.


To **find** reagent lots by Download status:

1. Select from the Download dropdown list.
2. Click the **Display Records** button.

| Criteria | Options | Description |
|----------|---------|---|
| Download | All | Displays reagent lots regardless of their download status. |
| | Yes | Only displays reagent lots selected for download to instruments. |
| | No | Only displays reagent lots not currently selected to download to instruments. |

Warning Messages

If a different menu option is selected before saving new or edited lot information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the add or edit screen to save the data.

If a required field is missing when lot information is saved, the following warning message is displayed with a reference to the specific required field (Lot Code, Lot No., Expiration Date and/or value for Low, Normal and High Range) that is missing:




- Click the **OK** button to clear the message and enter the required information.

If the Lot Number entered already exists in the database, the following error message is displayed:




- Click the **OK** button to clear the message.

If data entered in the Expiration Date field is not valid, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the valid information.

If records are not selected from the table before clicking the **Edit** or **Delete** action buttons, the following warning message is displayed:



- Click the **OK** button to clear the message and make a selection.

If a lot is selected for deletion with test data stored in the active database, the following error message is displayed:



- Click the **OK** button to clear the message.

9. Administrative

Administrative Overview

The **Administrative** menu option provides access to eleven screens for use in managing the Location Tree, User Accounts, the Database and Instrument Test Panels. The four **Organization Setup** screens allow you add or edit Organization, Facility, Department and Location information to the Tree. The **Facility Settings** screen allows you to manage facility-specific settings such as LIS result forwarding, notifications, warnings and review comments. The two **User Setup** screens allow you to create a new user or find and edit existing QCM user accounts and permissions. The **My Profile** screen allows you to set your preferred startup screen, adjust which alarm types to display or change your password. The two **Database** screens allow you to import operator information or check the database status. The **Panels** screen allows you to manage multi analyte instrument test panels. Each screen is available only if your account was granted this permission by the Administrator or POCC.

Organization Setup Screen

Administrative: Organization Setup >> Organization

The **Organization Setup** screens allow the Administrator or POCC to setup and maintain the information for the Facilities, Departments and Locations within the Organization.

- Additions, edits and deletions made to a Facility, Department or Location name are reflected in the **Location Tree**.
- When the system is first installed, the default for the Organization Name appears in the Tree and in the entry field as **Organization**.
- Only one organization is supported by the system.
- The Organization name and associated details can be edited.

Name an Organization

To **name** an Organization:

1. Move the mouse cursor over **Administrative** to display the menu options.

The screenshot shows the PrecisionWeb application interface. At the top, there is a navigation bar with tabs: Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. The 'Administrative' tab is active. Below the navigation bar, there is a user status bar showing 'admin', 'Log Off', and counts for 'No: 134', 'LIS: 55', and 'Op: 29'. On the left, a tree view shows the organizational structure: Medical Center Inc. (selected), Eastern Medical Center, Pediatrics, Ward A, Ward B, Southern Medical Center, and Western Medical Center. The main content area is titled 'ORGANIZATION SETUP: Organization Details'. It contains the following fields: * Organization Name (Medical Center Inc.), Address 1, Address 2, City, State, Zip, Country, Phone Number, Fax Number, Administrator, Phone Number, Ext., and E-Mail Address. A 'Save' button is located at the bottom of the form.

2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Organization** option.

* Organization Name is a **required** entry field.

3. Enter the name of the organization in the **Organization Name** field in the **Organization Details** table. Enter additional organization details in the associated fields.
4. Click the **Save** button to save the data.
 - The information entered will be stored in the database and the **Organization** name will appear in the Tree.

Edit Organization Name and Details

To **edit** information for an Organization:


1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Organization** option.
3. Click on the **Organization** name in the Tree so the details for the organization appear on the screen.
4. Add to or edit the fields in the **Organization Details** table.
5. Click the **Save** button to save the new or edited data.

The **Organization Setup** screen displays the information currently in the database.

* Organization Name is a **required** entry field.

Warning Messages

If a different menu option is selected before saving the **Organization Setup** information, the following warning message is displayed:



- Click the **OK** button to cancel the data entered and move to the selected menu option, or
- Click the **Cancel** button to clear the message and return to the **Organization Setup** screen.

If the **required Organization Name** field is empty in the **Organization Setup** screen, the following warning message is displayed:



- Click the **OK** button and enter the required information.

Facility Setup Screen

Administrative: Organization Setup >> Facility

Add a Facility Name and Details

To **add** a Facility to the Organization:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Facility** option.

The screenshot displays the PrecisionWeb application interface. On the left is a navigation tree with 'Medical Center Inc' highlighted. The main content area is titled 'FACILITY SETUP: Facility Details'. It contains the following fields:

- Organization Name: Medical Center Inc (read-only)
- * Facility Name: (required text input)
- Address 1: (text input)
- Address 2: (text input)
- City: (text input)
- State: (text input)
- Zip: (text input)
- Country: (text input)
- Time Zone: (No Selection) (dropdown menu)
- Phone Number: (text input)
- Fax Number: (text input)
- Administrator: (text input)
- Phone Number: (text input)
- Ext.: (text input)
- E-Mail Address: (text input)

At the bottom of the form are 'Save' and 'Delete' buttons.

3. Click on the **Organization** name in the Tree.
 - * Facility Name is a **required** entry field.
4. Enter the Facility name and associated information in the fields in the **Facility Details** table.
5. Click the **Save** button to save the new data.
 - The **Facility Setup** screen displays the Organization name in a read-only field.
 - The information entered will be stored in the database and the **Facility** name will appear in the Tree.

Note: The **Time Zone** dropdown list appears in Facility, Department and Location Setup screens and allows you to select the time zone for a specific facility, department or location, which will be displayed with test times associated with those locations. If no setting is specified, then "(No Selection)" appears in the field and no time zone is appended to the test times.

Edit Facility Name and Details

To **edit** a Facility in the Organization:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Facility** option.
3. Click on a **Facility** name in the Tree.
4. Add to or edit the fields in the **Facility Details** table.
5. Click the **Save** button to save the new or edited data.

The **Facility Setup** screen displays the Organization name in a read-only field and the Facility information currently in the database.

The screenshot displays the PrecisionWeb application interface. The top navigation bar includes links for Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. The left sidebar shows a tree view of the organization structure, with 'Eastern Medical Center' highlighted. The main content area is titled 'FACILITY SETUP: Facility Details' and contains the following fields:

- Organization Name: Medical Center Inc (read-only)
- * Facility Name: Eastern Medical Center (required)
- Address 1
- Address 2
- City
- State
- Zip
- Country
- Time Zone: CST (dropdown menu)
- Phone Number
- Fax Number
- Administrator: Mary Smith
- Phone Number
- Ext.
- E-Mail Address: msmith@emc.org

At the bottom of the form are two buttons: 'Save' and 'Delete'.

* The Facility Name is a **required** entry field.

Delete a Facility

To **delete** a Facility from the Organization:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Facility** option.


3. Click on a **Facility** name in the Tree.
4. Click the **Delete** button to remove the Facility name from the database and the Tree.

The **Facility Setup** screen displays the Organization name in a read-only field and the Facility information currently in the database.

Note: The delete action will remove the Facility Name only if there is **no** data (instruments, test results, terminal servers, operators or lots) associated with the specified Facility in the active database. If a Department or Location is associated with Facility, even if there is no data, the Facility cannot be deleted.

Warning Messages

If a different menu option is selected before saving the Facility Setup information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered.
- Click the **Cancel** button to clear the message and return to the **Facility Setup** screen to save the data.

If the **required Facility Name** field is empty in the Facility Setup screen, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the required information.

If the Facility name entered already exists in the database, the following error message is displayed:




- Click the **OK** button to clear the message and edit the Facility name.

If a Facility name is selected for deletion, the following warning message is displayed:



- Click the **OK** button to delete the facility.
- Click the **Cancel** button to cancel the delete action.

If a Facility name is selected for deletion with data stored in the active database, the following error message is displayed:



- Click the **OK** button to clear the message.

Department Setup Screen

Administrative: Organization Setup >> Department

Add a Department Name and Details

To **add** a Department to a Facility:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Department** option.
3. Click on a **Facility** name in the Tree.
4. Enter the Department name and associated information in the fields in the **Department Details** table.
5. Click the **Save** button to save the new data.

The screenshot displays the PrecisionWeb application interface. At the top, there is a navigation bar with tabs: Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. Below this, a left-hand navigation tree shows a hierarchy starting with 'Medical Center Inc', followed by 'Eastern Medical Center', 'Pediatrics', 'Ward A', 'Ward B', 'Southern Medical Center', and 'Western Medical Center' (which is highlighted). The main content area is titled 'DEPARTMENT SETUP: Department Details'. It contains the following fields: 'Organization Name' (pre-filled with 'Medical Center Inc'), 'Facility Name' (pre-filled with 'Western Medical Center'), '* Department Name' (marked as required), 'Manager', 'Phone Number', 'Fax Number', 'E-Mail Address', and 'Time Zone' (set to '(No Selection)'). There are 'Save' and 'Delete' buttons at the bottom of the form.

* Department Name is a **required** entry field.

- The **Department Setup** screen appears with the Organization and Facility names displayed in read-only fields.
- The information entered will be stored in the database and the **Department** name will appear in the Tree.

Edit Department Name and Details

To **edit** a Department in a Facility:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Department** option.
3. Click on a **Department** name in the Tree.
4. Add to or edit the fields in the **Department Details** table.
5. Click the **Save** button to save the new or edited data.

The **Department Setup** screen displays the Organization and Facility names in read-only fields and the Department information currently in the database.

The screenshot displays the PrecisionWeb interface for editing department details. On the left, a tree view shows the organizational structure, with 'Intensive Care' selected under 'Western Medical Center'. The main panel, titled 'DEPARTMENT SETUP: Department Details', contains the following fields:

- Organization Name: Medical Center Inc
- Facility Name: Western Medical Center
- * Department Name: Intensive Care
- Manager: Raymond Walker
- Phone Number: [empty]
- Ext.: [empty]
- Fax Number: [empty]
- E-Mail Address: rwalker@wmc.org
- Time Zone: PST

At the bottom of the form are 'Save' and 'Delete' buttons.

* Department Name is a **required** entry field.

Delete a Department Name

To **delete** a Department in a Facility:


1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Department** option.
3. Click on a **Department** name in the Tree.
4. Click the **Delete** button to remove the Department name from the database and the Tree.

The **Department Setup** screen displays the Organization and Facility names in read-only fields and the Department information currently in the database.

Note: The delete action will remove the Department Name only if there is **no** data (instruments, test results, terminal servers, operators or lots) associated with the specified Department in the active database. If a Location is associated with Department, even if there is no data, the Department cannot be deleted.

Warning Messages

If a different menu option is selected before saving the Department Setup information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered.
- Click the **Cancel** button to clear the message and return to the **Department Setup** screen to save the data.

If the **required Department Name** field is empty in the Department Setup screen, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the required information.

If the Department name entered already exists in the database in the same Facility, the following error message is displayed:




- Click the **OK** button to clear the message and edit the Department name.

If a Department name is selected for deletion, the following warning message is displayed:




- Click the **OK** button to delete the department.
- Click the **Cancel** button to cancel the delete action.

If a Department name is selected for deletion with test data stored in the active database, the following warning message is displayed:



- Click the **OK** button to clear the message.

If the Department Setup submenu option is selected and the **Organization** is highlighted in the Tree, the following warning message is displayed:



- Click the **OK** button to clear the message and select a Facility or a Department from the Tree.

Location Setup Screen

Administrative: Organization Setup >> Location

Add a Location Name and Details

To **add** a Location to a Department:

- Move the mouse cursor over **Administrative** to display the menu options.
- Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Location** option.
- Click on a **Department** name in the Tree.
- Enter the Location name and associated information in the in the fields in the **Location Details** table.
- Click the **Save** button to save the new data.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin: LogOff

Refresh Hide Tree

LOCATION SETUP: Location Details

Organization Name

Facility Name

Department Name

* Location Name

Manager

Phone Number Ext.

Fax Number

E-Mail Address

Time Zone (No Selection)

Save Delete

0 Records Found

| Edit | Delete | Hostname | IP Address | MAC Address | COM Port | Comment |
|--------------------------|--------|--|------------|-------------|----------|---------|
| <input type="checkbox"/> | All | No data is available for this selection. Please change selection criteria or tree selection. | | | | |

0 Records Found

* The Location Name is a **required** entry field.

- The **Location Setup** screen appears with the Organization, Facility and Department names displayed in read-only fields.
- The information entered will be stored in the database and the **Location** name will appear in the Tree.

Edit Location Name and Details

To **edit** a Location:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Location** option.
3. Click on a **Location** name in the Tree.
4. Add to or edit the fields in the **Location Details** table.
5. Click the **Save** button to save the new or edited data.

PrecisionWeb Dashboard Review Reports Operators Instruments Lots Administrative Help admin | Log Off

Refresh Hide Tree x

LOCATION SETUP: Location Details No:134 LIS:58 Op:29

Location Details

Organization Name: Medical Center Inc

Facility Name: Eastern Medical Center

Department Name: Pediatrics

* Location Name: Ward B

Manager: Viola Allen

Phone Number: Ext.:

Fax Number:

E-Mail Address: vallen@emc.org

Time Zone: EST

Save Delete

0 Records Found

| Add | Edit | Delete | Hostname | IP Address | MAC Address | COM Port | Comment |
|--------------------------|------|--------|----------|------------|-------------|----------|---------|
| <input type="checkbox"/> | All | | | | | | |

No data is available for this selection. Please change selection criteria or tree selection.

| Add | Edit | Delete |
|-----|------|--------|
|-----|------|--------|

* Location Name is a **required** entry field.

The **Location Setup** screen displays the Organization, Facility and Department names in read-only fields and the Location information currently in the database.

Delete a Location Name and Information

To **delete** a Location in a Department:

1. Move the mouse cursor over **Administrative** to display the menu options.

2. Move the mouse cursor over the **Organization Setup** option until the submenu appears and click on the **Location** option.
3. Click on a **Location** name in the Tree.
4. Click the **Delete** button to remove the Location name from the database and the Tree.
 - The **Location Setup** screen displays the Organization, Facility and Department names in read-only fields and the Location information currently in the database.

Note: *The delete action will remove the Location Name only if there is **no** data (instruments, test results, terminal servers, operators or lots) associated with the specified Location in the active database.*

Add a Terminal Server to a Location

To add a Terminal Server to a Location, identify the Terminal Server by adding at least one of the Terminal Server's identifiers (Hostname, IP Address, or MAC Address) to the Location:

From the **Location Setup** screen,

1. Click on a **Location** name in the Tree.
2. Click the **Add** button on the **Action Bar** of the **Hostname/IP Address** table to display the **Terminal Server Details** screen.
3. Enter a Hostname, IP Address or MAC Address for the specified Location.
4. The **Comment** field allows the Administrator to record a comment and is an optional field.
5. Click the **Save** button to save the new data and return to the **Location Setup** screen, or Click the **Cancel** button to return to the **Location Setup** screen without saving the new data.

Notes:

- *In a DHCP environment, you only should add the Hostname. The IP Address should update after approximately 5 – 10 minutes after Hostname is initially populated or after any IP Address changes.*
- *If DHCP is being used, network topology and configuration will determine if the Hostname must be added as a Fully Qualified Name (FQN).*
- *If IP Addresses are statically assigned the IP Address and Hostname should be added to the Location.*
- *The next time the instrument uploads data from this location to the database, the remaining fields are detected and automatically refreshed.*

The **Location Setup** screen displays the information stored in the database for the Location and the action button to add a Hostname, IP Address and/or a MAC address.

* You must fill in at least one of the specified fields: Hostname, IP Address or MAC-Address.

Edit Terminal Server Details

To edit terminal server details for a Location:

From the **Location Setup** screen,

1. Click on a **Location** name in the Tree.
2. Click on the box in the Hostname column to select the record. Click the **Edit** action button to display the **Edit Terminal Server** screen.
3. Edit the Hostname, IP Address and/or MAC Address.
4. The **Comment** field allows the Administrator to record a comment and is an optional field.
5. Click the **Save** button to save the edited data and return to the **Location Setup** screen, or Click the **Cancel** button to return to the **Location Setup** screen without saving the new data without saving.

The **Location Setup** screen displays the information stored in the database for the Location and the current Hostname/IP Address and/or MAC Address information.

* Either a Hostname, an IP Address or MAC Address is **required**.

Delete a Terminal Server

To delete a terminal server for a location:


From the **Location Setup** screen,

1. Click on a **Location** name in the Tree.
2. Click on the box in the Hostname column to select the record and click the **Delete** action button.
3. Click the **OK** button to delete the Hostname/IP Address or MAC Address of the selected Location.

Note: *The delete action will remove the terminal server entry for the specified Location.*


Warning Messages

If a different menu option is selected before saving the Location Setup information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the **Location Setup** screen to save the data.

If the **required** field for the Location Name is missing on the **Location Setup** screen, the following warning message is displayed:




- Click the **OK** button to clear the message and enter the required information.

If the Location name entered already exists in the database in the same Department, the following error message is displayed:




- Click the **OK** button to clear the message and edit the Location name.

If a Location Name is selected for deletion, the following warning message is displayed:




- Click the **OK** button to delete the Location, or
- Click the **Cancel** button to cancel the delete action.

If a Location name is selected for deletion with test data stored in the active database, the following error message is displayed:




- Click the **OK** button to clear the message.

If the Location Setup submenu option is selected and the **Organization** or a **Facility** is highlighted in the Tree, the following warning message is displayed:



- Click the **OK** button to clear the message and select a Department or Location from the Tree.

If no Hostname/IP Address or MAC Address record is selected before clicking the **Edit** or **Delete** action button, the following warning message is displayed:



- Click the **OK** button to clear the message and select a record from the table.

Facility Settings

Administrative>>Facility Settings

The **Facility Settings** screen allows you to select the settings for a Facility, including:

- Selecting a Result Forwarding option,
- Selecting when advance warning alarms are displayed on the Dashboard, and
- Setting Operator Performance Target Levels for passing QC.
- If more than one Facility exists within the Organization, each Facility must be configured separately, as settings are Facility specific.

To select the settings for a **Facility**:

1. Move the mouse cursor over **Administrative** to display the menu options and click on the **Facility Settings** option.
2. Click on a **Facility** Name in the Tree.
3. Enter settings for the specified Facility.

The screenshot shows the PrecisionWeb interface. The top navigation bar includes links for Dashboard, Review, Reports, Operators, Instruments, Lots, Administrative, and Help. The user is logged in as 'admin' with a 'Log Off' option. The left sidebar shows a tree view of the organization structure, with 'Western Medical Center' selected under 'Southern Medical Center'. The main content area is titled 'FACILITY SETTINGS: Western Medical Center'. It contains several sections: 'Result Forwarding' with 'On' and 'Off' radio buttons; 'Expiration Date Notifications' with input fields for Operator Certification, Reagent Lot, Control Lot, and Linearity Lot, each followed by 'Days in Advance'; 'Performance Target Levels' with input fields for '* Excellent' (90%) and '* Acceptable' (80%); 'Meter Download Warning' with an input field for '* Interval' (1) followed by 'Hours'; and 'Facility-Defined Review Comment' with ten text input fields. A 'Save' button is located at the bottom right of the form.

4. Click the **Save** button to save the data.

The **Facility Settings** screen displays the Facility Name selected in the screen title.

Select settings from the following options:

| Setting | Description |
|---------------------------------------|--|
| Result Forwarding | On - forwards test results to the LIS/HIS for processing. Off - holds test results in a queue (detail not visible in the UI). |
| Expiration Date Notifications: | |
| Operator Certification | Enter the number of days for advance warning notification. |
| Reagent Lot | Enter the number of days for advance warning notification. |
| Control Lot | Enter the number of days for advance warning notification. |
| Linearity Lot | Enter the number of days for advance warning notification. |
| Performance Target Levels: | |
| *Excellent | Enter the target level for Excellent QC performance. |
| *Acceptable | Enter the target level for Acceptable QC performance. |
| Meter Download Warning | |
| *Interval | Enter the number of hours allowed between instrument uploads. |

Note: The summary screen will display red highlighted text in the Last Upload Date/Time field if the specified upload time limit is exceeded. The Meter Download Warning Interval setting is independent of instrument Upload Interval setting.

| | |
|--|---|
| Facility-Defined Review Comment | <p>These additional fields will allow you to specify up to 10 Review Comments that will then be displayed in the Predefined Review Comment list of the Edit Result screen. Only entered comments will be displayed in the list.</p> <p>You can delete a comment by emptying a field and selecting Save. You may leave fields empty. The empty fields will be ignored when displaying the list of Predefined Review Comments. If a Review Comment is deleted, however, the results to which this comment was attached will still display the comment.</p> |
|--|---|

The settings selected will be stored in the database.

* The QC Performance Target Levels of Excellent and Acceptable and the Meter Download Warning Interval are **required** entry fields.

New User Setup Screen

Administrative: User Setup >> New

- The **New User** screen allows the Administrator or POCC to add users to the system.
- A user is granted permission to access locations and menu options in the system to meet the user's responsibilities within the organization.
- This screen also allows the Administrator or POCC to permit the user to perform specific functions, such as edit records and download reagent lots.

Add a new user

To **add** a new user:

1. Move the mouse cursor over **Administrative** to display menu options.
2. Move the mouse cursor over the **User Setup** option until the submenu appears.
3. Click the **New** option.

NEW USER:

Identification

* First Name

Middle Name

* Last Name

* Employee ID

E-Mail

* User ID

Permit Access To Locations

Eastern Medical Center ☐ All ☒ None ☐ Selected ...

Southern Medical Center ☐ All ☒ None ☐ Selected ...

Western Medical Center ☐ All ☒ None ☐ Selected ...

QCM Menu Permissions

☐ Dashboard

☐ Review

☐ Reports

☐ Operators

☐ Instruments

☐ Lots

☐ Administrative

☐ Help

QCM Function Permissions

☐ Edit Patient ID

☐ Edit All Result Details

☐ Forward To LIS

☐ Edit Operator Certification

☐ Edit Operator Information

☐ Download Reagents

☐ Edit Network Settings

☐ Instrument Operations

The **New User** screen has three sections:

- Identification
- Permit Access to Locations
- QCM3 Menu Permissions and QCM3 Function Permissions

Identification

Enter the identification information for the new user.

* First Name, Last Name, Employee ID and User ID are **required** entry fields.


Note: The application automatically enters the first five characters of the user's Last Name and the initial character of the First and Middle Names for the new user's User ID. The default User ID can be edited prior to saving the new user entry.


Permit Access to Locations

The selections made in this section determine the Locations displayed in the Tree. If the user is not granted permission to access a specific Facility, Department or Location, it will not appear in the Tree.

To permit a user **Access** to Facilities in the organization:

- Click on the **All** button to allow the user access to all locations in that Facility.
- Click on the **None** button to deny the user access to a specific Facility and its Departments and Locations.
- Click on the **Selected** button to allow the user access to selected Locations within a Facility, then:


1. Click the  button. A pop-up browser window appears with a list of associated Departments and Locations for the selected Facility.



2. Click on the box in the Department Name column to allow the user access to associated Location listed in the Location Name column and click the **OK** button, or
3. Click the **Cancel** button to cancel access to the selected locations and close the pop-up browser window.

QCM3 Menu Permissions









The selections made in this section determine the options displayed in the **Menu Bar**. If the user is not granted permission to access a menu or sub-menu option it will not appear in the Menu Bar.

- To permit a user access to **all** associated menu options click on the checkbox to the left of the menu option.
- To permit a user access to **some** associated menu options click on the checkbox to the left of the menu option, then click on the  arrow to display the associated menu options. The arrow points down when the associated menu are displayed. The arrow points to the right if the associated menu options are hidden from view.

For example, the **Review** menu option includes the following sub-menu options:

- Summary
- Alarms
- All Results
- Patient Results
- LIS Patient Lookup
- Click on the checkbox to select/deselect the associated options.
- A checked box means the user has access to that menu option.
- An unchecked box means the user is denied access to that menu option.
- The Administrative menu option will remain checked to assign the My Profile screen even if all sub-menu options remain unchecked.

QCM Menu Permissions

| | | |
|---|-------------------------------------|--------------------|
|  | <input type="checkbox"/> | Dashboard |
|  | <input type="checkbox"/> | Review |
| | <input type="checkbox"/> | Summary |
| | <input type="checkbox"/> | Alarms |
| | <input type="checkbox"/> | All Results |
| | <input type="checkbox"/> | Patient Results |
| | <input type="checkbox"/> | Patient Lookup |
|  | <input type="checkbox"/> | Reports |
| | <input type="checkbox"/> | View |
| | <input type="checkbox"/> | Suites |
| | <input type="checkbox"/> | Import |
|  | <input type="checkbox"/> | Operators |
| | <input type="checkbox"/> | New |
| | <input type="checkbox"/> | Certification |
| | <input type="checkbox"/> | Search |
|  | <input type="checkbox"/> | Instruments |
| | <input type="checkbox"/> | Configuration |
| | <input type="checkbox"/> | Comments |
| | <input type="checkbox"/> | Ports |
| | <input type="checkbox"/> | Tracking Log |
| | <input type="checkbox"/> | Upload Log |
|  | <input type="checkbox"/> | Lots |
| | <input type="checkbox"/> | Control |
| | <input type="checkbox"/> | Linearity |
| | <input type="checkbox"/> | Proficiency |
| | <input type="checkbox"/> | Reagent |
|  | <input checked="" type="checkbox"/> | Administrative |
| | <input type="checkbox"/> | Organization Setup |
| | <input type="checkbox"/> | Facility Settings |
| | <input type="checkbox"/> | User Setup |
| | <input type="checkbox"/> | Database |
| | <input type="checkbox"/> | Panels |
|  | <input type="checkbox"/> | Help |

QCM3 Function Permissions

The selections made in this section allow the user to perform specific functions in the system. If the user does not have permission to perform a function, the corresponding Action Button is not displayed in the Action Bar.

QCM Function Permissions

Edit Patient ID ☐

Edit All Result Details ☐

Forward To LIS ☐

Edit Operator Certification ☐

Edit Operator Information ☐

Download Reagents ☐

Edit Network Settings ☐

Instrument Operations ☐

- Click on the checkbox to allow the user to perform that function.
- An unchecked box means the user is not allowed to perform that function.


| | | |
|------------------------------------|--------------------------|--|
| Edit Patient ID | <input type="checkbox"/> | Allows user permission to edit the Patient ID associated to a record. This option is disabled if Edit All Result Details (below) is selected. |
| Edit All Result Details | <input type="checkbox"/> | Allows user permission to edit all result details (including Patient ID) with the exception test date and time. |
| Forward to LIS | <input type="checkbox"/> | Allows user permission to forward results to the LIS/HIS and Bypass (LIS/HIS) Error Check. |
| Edit Operator Certification | <input type="checkbox"/> | Allows user permission to add or edit Operator information, including home Location. User may also add or edit (ReCertify) instrument certification. |
| Edit Operator Information | <input type="checkbox"/> | Allows user permission only to add or edit Operator information, including home Location. |
| Download Reagents | <input type="checkbox"/> | Allows user permission to add reagent lots to instrument reagent lot download list. |
| Edit Network Settings | <input type="checkbox"/> | Allows user permission to add, edit and delete terminal server information. |
| Instrument Operations | <input type="checkbox"/> | Allows user permission to rename or relocate (drag and drop) an instrument in the Tree. |

When finished selecting options for the new user in all three sections of the **New User** screen, click the **Save** button to save the new data.

The information entered and the access privileges selected will be stored in the database.


Warning Messages

If a different menu option is selected before saving the User information, the following warning message is displayed:



- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the **New User** screen to save the data.

If the Employee ID entered already exists in the database in the same Facility, the following warning message is displayed:



- Click the **OK** button to clear the message and edit the Employee ID.

If the **required** field(s) for First Name, Last Name, User ID and Employee ID are missing when the New User information is saved, the following warning message is displayed:



- Click the **OK** button to clear the message and enter the required information.

User Search Screen

Administrative: User Setup >> Search

The **User Search** screen allows the user to find, edit and delete users.

To **view** users:

1. Move the mouse cursor over **Administrative** to display menu options.
2. Move the mouse cursor over the **User Setup** option until the submenu appears and click on the **Search** option.

The **User Search** screen displays the Search Criteria entry fields and the User table.

USER SEARCH:

Search Criteria

Last Name

Employee ID

Display Records

13 Records Found

| <input type="checkbox"/> All | Name | Employee ID | E-Mail |
|-------------------------------------|----------------------|-----------------------------|------------------------|
| <input type="checkbox"/> | admin, admin | 1 | -- |
| <input type="checkbox"/> | Allen, Viola | 2696 | vallen@emc.org |
| <input type="checkbox"/> | Bailey, Lori | 3565 | lbailey@wmc.org |
| <input type="checkbox"/> | Bell, William | 3564 | wbell@wmc.org |
| <input type="checkbox"/> | Booth, Debra | 3563 | dbooth@wmc.org |
| <input type="checkbox"/> | Dunby, Charles | 3561 | cdunby@wmc.org |
| <input type="checkbox"/> | Grant, David | 3567 | dgrant@wmc.org |
| <input type="checkbox"/> | Hayes, Robert | 96544 | rhayes@mci.org |
| <input type="checkbox"/> | Keene, Brent | 2695 | bkeene@emc.org |
| <input type="checkbox"/> | Smith, Mary | 2693 | msmith@emc.org |
| <input type="checkbox"/> | Stone, Beverly | 3566 | bstone@wmc.org |
| <input type="checkbox"/> | Walker, Raymond | 3562 | rwalker@wmc.org |
| <input type="checkbox"/> | Williams, Ross | 2694 | rwilliams@emc.org |

Edit Delete

13 Records Found

To **search** for a user:

1. Enter a **Last Name** and/or an **Employee ID** and click the **Display Records** button.
2. **Results** that match or are similar to the name and/or ID entered are displayed in the User table.

Edit Users

To **edit** a single user or multiple users:

1. Select one or more records in the table by clicking on the box(es) in the Name column.
2. Click the **Edit** action button.

To **edit** all users:

1. Click on the box marked **All** in the Name column header.
2. Click the **Edit** action button.

The **Edit User** screen has the same three sections as the **New User** screen (see Section [New User Setup Screen](#), above).

Note: All User Identification fields may be modified with the exception of the User ID field. Permit Access To Locations, Menu Permissions and Function Permissions may also be modified.

My Profile Screen

Administrative: My Profile

The **My Profile** screen allows you to customize the behavior of the system for your User ID. The My Profile settings are unique to every user. All users regardless of their assigned Menu Permissions have access to the My Profile screen.

To select **My Profile** options:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Click on the **My Profile** option.

MY PROFILE:

Profile Options

Startup Screen (No Selection) ▼

Default Instrument PCx ▼

Alarm Preferences

☒ Show Total Alarm Counts

Show Alarms of Type

☒ New Alarms

☒ Expiration Warning Alarms

☒ LIS Transfer/Generated Alarms

☒ Operator Alarms

☒ Test Alarms

☒ Lot Alarms

Password

☐ Change Password

Old

New

Confirm

Save

- Your location in the **Tree** is irrelevant as users are not associated with a certain facility, department or location.
- The **My Profile** screen displays the Profile Options available to customize the system.

Select a Startup Screen

Select your choice of a Startup Screen from the dropdown list under the **Profile Options** header. Only menu options assigned to you will appear as an option in the Startup Screen dropdown list.

MY PROFILE:

Profile Options

Startup Screen (No Selection)

Alarm Preferences

Dashboard

Review > Summary

Review > Alarms

Review > All Results

Review > Patient Results

Reports > View

Reports > Suites

Reports > Import

Operators > New

Operators > Certification

☒ Test Alarms

☒ Lot Alarms

Password

☐ Change Password

Old

New

Confirm

Save

The system will display the selected Startup Screen each time you log on.

Default Instrument

Select a default instrument from the dropdown list. The default instrument will display in all instrument fields on all screens.

Alarm Preferences

- Click on the checkbox to **enable** the **Show Total Alarm Counts** option. The Total Alarm Count will appear in the top right corner of every screen (except for the log on screen).
- An unchecked box means the **Show Total Alarm Counts** option is disabled.

Show Alarms of Type

Click on the checkbox(es) to select the warning alarm types the system will display on the Dashboard. These selections will recalculate the total alarm counts based on which alarms are selected. Disabling an Alarm Type means that the associated alarms will not be displayed on the Dashboard, however, alarms will be generated as applicable. Turning on Alarm Type category at a later date will display associated alarms on the Dashboard along with all unacknowledged alarms.

- ☐ New Alarms
- ☐ Expiration Warning Alarms
- ☐ LIS Transfer/Generated Alarms
- ☐ Operator Alarms
- ☐ Test Alarms
- ☐ Lot Alarms

The following lists each My Profile Alarm Type category and the corresponding Alarms and Dashboard category:

| <u>Alarm Type</u> | <u>Alarm Description</u> | <u>Dashboard Category</u> |
|--------------------------|---------------------------------|----------------------------------|
| New | Instrument Moved | Notifications |
| New | New Reagent Lot | Notifications |
| New | New Control Lot | Notifications |
| New | New Linearity Lot | Notifications |
| New | New Proficiency Lot | Notifications |
| New | New Instrument | Notifications |
| New | New Operator | Operator |
| Expiration Warning | Reagent Exp. Warning | Notifications |
| Expiration Warning | Control Exp. Warning | Notifications |
| Expiration Warning | Linearity Exp. Warning | Notifications |
| Expiration Warning | Operator Exp. Warning | Operator |
| LIS Transfer/Generated | Out-of-Range Result - LIS | LIS/Test Alarms |
| LIS Transfer/Generated | Invalid Operator - LIS | LIS/Test Alarms |
| LIS Transfer/Generated | Invalid Patient - LIS | LIS/Test Alarms |
| LIS Transfer/Generated | Test Rejected - LIS | LIS/Test Alarms |
| LIS Transfer/Generated | Scripting Error | LIS/Test Alarms |
| LIS Transfer/Generated | Result Edit Not Forwarded | LIS/Test Alarms |
| LIS Transfer/Generated | Custom LIS Alarm 1-10 | LIS/Test Alarms |
| LIS Transfer/Generated | Repeat Test | LIS/Test Alarms |
| Operator | Invalid Operator - Inst | Notifications |
| Operator | Expired Op. Certification | Operator |
| Test | Invalid QC - Inst | LIS/Test Alarms |

| <u>Alarm Type</u> | <u>Alarm Description</u> | <u>Dashboard Category</u> |
|-------------------|---------------------------------|---------------------------|
| Test | Test Error | Notifications |
| Test | Out-of-Range QC | LIS/Test Alarms |
| Test | Out-of-Range Result | LIS/Test Alarms |
| Test | Patient ID Not Confirmed - Inst | LIS/Test Alarms |
| Lot | Expired Reagent Lot | Notifications |
| Lot | Expired Control Lot | Notifications |
| Lot | Expired Linearity Lot | Notifications |
| Lot | Invalid Reagent - Inst | Notifications |
| n/a | Rejected Result Alert | Notifications |
| n/a | Not Connected | Instrument Summary |
| n/a | Upload Required | Instrument Summary |
| n/a | Upload Location Warning | Instrument Summary |

Note: These alarms will still show in the Alarms and Dashboard screens even if the **Show Total Alarm Counts** option is disabled.

Password


- Click on the checkbox to change the current password.

| | |
|----------------|--------------------------------------|
| [add box] | Change Password. |
| Old | Enter current password. |
| New | Enter a new password. |
| Confirm | Enter new password again to confirm. |

- Click the **Save** button to save the data.


Warning Messages

If a different menu option is selected before saving the selected options, the following warning message is displayed:



- Click the **OK** button to cancel the data entered, or
- Click the **Cancel** button to clear the message and return to the My Profile screen to save the data.

If the Confirm password entered does not match the new password entered or the password does not have between five and 30 characters, the following error message is displayed:



- Click the **OK** button to clear the message and re-enter the password.

Database Screens

Database Import


Administrative: Database >> Import

Operator data files must be imported in Comma delimited (.csv) format.


- If Operators will be assigned to different Locations the Location specific Operator data must be saved in separate files, one .csv file per Operator Location.
- Valid instrument types are Abbott glucose monitors and the i-STAT analyzer.

To **save** the data from an Excel spreadsheet in Comma delimited (.csv) format:

1. Open the Excel file with the operator data.
2. Click the arrow in the **Save as type:** entry field to display the dropdown list.
3. Click on the **CSV (Comma delimited) (*.csv)** file type from the dropdown list.
4. The **File name:** entry field will display your file name with the extension .csv.
5. The following message will appear:




6. Click **OK**.
7. The following message will appear:



8. Click **Yes**.

9. Select **File, Save** from the menu to save the changes. You will be prompted with the same messages boxes, click **OK** and then **Yes**.
10. Select File, Exit from the menu (or click X in the top right corner) to close the file. You will be prompted with the following message;



11. Select No. The file will close, but the changes have been saved.
12. Right click on the saved .csv file, select Open With from the menu and then Notepad from the program list. The file will open in Notepad.
13. Hit the Enter key once on your keyboard to create an empty line just above the operator data.
14. Enter the instrument type certification (i.e. PCx) in the empty space, just above the operator data.
15. Select File, Save from the menu and close the file.

To **import** data to the system:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Move the mouse cursor over the **Database** option until the submenu appears.
3. Click on the **Import** option to display the **Database Import** screen.

A screenshot of the "DATABASE IMPORT:" screen. The title "DATABASE IMPORT:" is in red. Below it is a section titled "Input Details" with a light blue background. Inside this section, there is a "Data For" dropdown menu set to "Operators", a "Mode" section with radio buttons for "Append" (selected) and "Overwrite", and a "* File Location" text box with a "Browse..." button next to it. Below the "Input Details" section is an "Import" button. At the bottom, there is a section titled "Import Status" with a large empty text area and a scrollbar on the right.


To **import** information into the database:

1. Select the **Append** Mode to **add** the imported information to the existing data, or
2. Select the **Overwrite** Mode to **overwrite** the existing data for operators with the imported information.
3. Select the **File Location** of the operator .csv file by clicking on the **Browse** button.

A new browser window appears with a file system browser for the PC running software. Depending on your location's network, it may display only files on the local PC or files available on network drives.

4. Locate and select the file to import and click the **Open** button on the browser window.
5. The entry field on the **Database Import** screen displays the file name selected in the **File Location** entry field.
6. Click the **Import** button to begin the import of data.

If import is successful a message will display in the Import Status window indicating the number of records imported successfully and/or the number of records rejected and the reason for the error. If any records are rejected correct the error in the .csv file and attempt the import again.



The screenshot shows a window titled "Import Status" with a text area containing the following text:

```
Thu Nov 2 16:20:35 EST 2006: Importing 'H:\My Documents\POC
Training\Product Training\QC Manager 3.x\Operator Export_Import\Operator
Lists\operator list.csv'.
Operator List Import in progress...

3 record(s) imported successfully.

2 record(s) rejected due to errors (error list truncated at 100):
-- Invalid Certification date specified. [OperatorId =]
-- Invalid Expiration date specified. [OperatorId =]
-- The required Operator ID field is missing. [OperatorId =]
-- Invalid Certification date specified. [OperatorId =]
-- Invalid Expiration date specified. [OperatorId =]
-- The required Operator ID field is missing. [OperatorId =]

Thu Nov 2 16:20:38 EST 2006: Import Complete.
```

Importing Data from an Excel Spreadsheet

To **import** data from an Excel spreadsheet:

Create a table of operators using an Excel spreadsheet:

When creating the spreadsheet for operators the first row must remain blank.

Columns A through I must be populated as follows;

Operator Import List

| | |
|----------|---|
| Row 1 | Leave blank |
| Column A | Operator ID - required |
| Column B | Last Name - optional |
| Column C | First Name - optional |
| Column D | Certification Date - required |
| Column E | Expiration Date - required |
| Column F | Certification Interval - required |
| Column G | 0 (the number zero) - required |
| Column H | 1 (the number one) - required |
| Column I | Original Certification Date - required |

Valid entries for Certification Interval are:

3 (3 months), 6 (6 months), 12 (12 months), 0 (Never Expire), -1 (Specify), -2 (Expired).

Date entered in columns D, E and I must be formatted with the four digit year (see example spreadsheet below)

Note: Field I (Original Certification Date) is a required field, however, Original Certification Date information will not be populated as part of the operator certification record. Original Certification Date must be updated manually via the Operator Certification screen.


| | A | B | C | D | E | F | G | H | I | |
|----|------|----------|-----------|----------|----------|----|---|---|-----------|--|
| 1 | | | | | | | | | | |
| 2 | 1000 | Smith | Elizabeth | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 3/16/2001 | |
| 3 | 1001 | Taylor | Ann | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 3/16/2001 | |
| 4 | 1002 | Jones | Mary | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 7/1/2006 | |
| 5 | 1003 | Miller | Kathleen | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/15/2004 | |
| 6 | 1004 | Flaherty | Margaret | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/22/2002 | |
| 7 | 1005 | Martin | Jennifer | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/24/2001 | |
| 8 | 1006 | Hurt | Sally | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/11/2005 | |
| 9 | 1007 | Price | Stacey | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 3/16/2001 | |
| 10 | 1008 | Harris | Esther | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/8/1999 | |
| 11 | 1009 | Jordon | Irene | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 8/18/2002 | |
| 12 | 1010 | Horton | Barbra | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 5/12/2001 | |
| 13 | 1011 | Mills | Kathryn | 7/1/2006 | 7/1/2007 | 12 | 0 | 1 | 6/1/2004 | |

Database System Status

Administrative: Database > System Status

The **System Status Log** screen displays information on QCM managed system backup, restore and archive functions as well as any problem status recorded in the past 30 days.. *Backup activity managed via SQL Server is not noted in the System Status Log.*

- Displays system backup, restore and archive functions
i.e. 07/25/2006 08:45 Backup started. Device = Tape
- Records communication interruption with instrument
i.e. 08/29/2006 22:01 PCx Component Meter stopped responding - Please check the meter's battery: M07610053
- Reports other system issues as they occur
i.e. 09/08/2006 10:13 QCM Service Control Service Control found QCM Communication Service not running – 9/8/2002 10:13:25.968



```
SYSTEM STATUS LOG: No:151 LIS:101 Op:8

02/15/2007 02:00 PM EXPECTED ERROR
      Invalid results xml - IPAddress not in tbTerminal and MAC address not specifi
      Please resubmit results file MCSSTRERR2204BLK0000TMO21507140019TID0324ORND023

02/15/2007 02:00 PM EXPECTED ERROR
      Invalid results xml - IPAddress not in tbTerminal and MAC address not specifi
      Please resubmit results file MCSSTRERR2204BLK0000TMO21507140019TID01096RND053

02/15/2007 02:04 PM Archive started.

02/15/2007 02:04 PM Archive completed sucessfully.

02/15/2007 02:06 PM Archive started.

02/15/2007 02:06 PM Archive completed sucessfully.

02/15/2007 02:06 PM Backup started. Device = File

02/15/2007 02:07 PM Backup completed sucessfully.
```

Panels Screen

Administrative>>Panels

The **Panels** screen allows you to add, edit and delete instrument test panel information for multi analyte instruments. Creating a test panel allows the user to view multiple analytes in a summary format via the Review Patient Results screen.


The test panels added will be listed in the Panel Name dropdown list in the Patient Results screen. The test panels are accessible to the entire organization (Facilities, Departments, and Locations).

To **view** test panels:

1. Move the mouse cursor over **Administrative** to display the menu options.
2. Click on the **Panels** option.
3. In the **Search Criteria** section of the screen use the **Panel Name** filter option to display all panel names or a specific panel from the dropdown list.
4. Click the **Display Records** button to display records that fit the chosen filter option.

PANELS:

Search Criteria

Panel Name 

Display Records

1 Record Found

| Add Edit Delete | |
|-------------------------------|-------------------------------------|
| <u>Panel Name</u> ^ | <u>Tests</u> |
| <input type="checkbox"/> All | |
| <input type="checkbox"/> test | Cl, BE, aPTT WBT, aPTT PET, ACT WBT |
| Add Edit Delete | |

1 Record Found

Add a Test Panel

To **add** a test panel:

1. Click the **Add** action button.
 - The **Add Panel** screen displays the entry field for the test panel name and the test selection options.

* Panel Name is a **required** entry field.

2. Select tests by checking the appropriate boxes.
3. Click the **Save** button to save the data. The information entered will be stored in the database. The saved panel will be available in the **Panel** dropdown list in the **Patient Results** screen.

ADD PANEL:
Panel Information
 * Panel Name
Test Selection

| | |
|---|---|
| <input type="checkbox"/> ACT WBT | <input type="checkbox"/> AnGap |
| <input type="checkbox"/> aPTT PET | <input type="checkbox"/> aPTT Ratio |
| <input type="checkbox"/> aPTT WBT | <input type="checkbox"/> B/Cre |
| <input type="checkbox"/> BE | <input type="checkbox"/> BNP |
| <input type="checkbox"/> BUN | <input type="checkbox"/> CK-MB |
| <input type="checkbox"/> Cl | <input type="checkbox"/> Crea |
| <input type="checkbox"/> cTnI | <input type="checkbox"/> EQC |
| <input type="checkbox"/> Glu | <input type="checkbox"/> Hb |
| <input type="checkbox"/> HCO ₃ | <input type="checkbox"/> Hct |
| <input type="checkbox"/> iCa | <input type="checkbox"/> INR |
| <input type="checkbox"/> K | <input type="checkbox"/> Ketone |
| <input type="checkbox"/> Lac | <input type="checkbox"/> Na |
| <input type="checkbox"/> PCO ₂ | <input type="checkbox"/> PCO ₂ (patient temp.) |
| <input type="checkbox"/> pH | <input type="checkbox"/> pH (patient temp.) |
| <input type="checkbox"/> PO ₂ | <input type="checkbox"/> PO ₂ (patient temp.) |
| <input type="checkbox"/> PT | <input type="checkbox"/> PT PET |
| <input type="checkbox"/> PT Ratio | <input type="checkbox"/> sO ₂ |
| <input type="checkbox"/> TCO ₂ | <input type="checkbox"/> U/Cre |
| <input type="checkbox"/> Urea | |

The **Panels** screen appears with the new test panel displayed in the **Test Panel** table.

4. Click the **Cancel** button to cancel the Add action.

Edit a Test Panel



To **edit** a single test panel or multiple test panels:

1. Click the box(es) in the Panel Name column.
2. Click the **Edit** action button.

To **edit** all test panels:



1. Click the box marked **All** in the Panel Name column header.
2. Click the **Edit** action button.

The **Edit Panel** screen appears with the information stored in the database for the first selected test panel record.

If multiple test panel records are selected each record is displayed on a separate page. Use the  and  buttons to change pages.

EDIT PANEL:
Panel Information
 * Panel Name
Test Selection

| | |
|--|---|
| <input type="checkbox"/> ACT WBT | <input type="checkbox"/> AnGap |
| <input type="checkbox"/> aPTT PET | <input type="checkbox"/> aPTT Ratio |
| <input type="checkbox"/> aPTT WBT | <input checked="" type="checkbox"/> B/Cre |
| <input type="checkbox"/> BE | <input type="checkbox"/> BNP |
| <input checked="" type="checkbox"/> BUN | <input type="checkbox"/> CK-MB |
| <input type="checkbox"/> Cl | <input type="checkbox"/> Crea |
| <input checked="" type="checkbox"/> cTnI | <input type="checkbox"/> EQC |
| <input type="checkbox"/> Glu | <input type="checkbox"/> Hb |
| <input type="checkbox"/> HCO3 | <input type="checkbox"/> Hct |
| <input type="checkbox"/> iCa | <input type="checkbox"/> INR |
| <input type="checkbox"/> K | <input checked="" type="checkbox"/> Ketone |
| <input type="checkbox"/> Lac | <input type="checkbox"/> Na |
| <input checked="" type="checkbox"/> PCO2 | <input type="checkbox"/> PCO2 (patient temp.) |
| <input type="checkbox"/> pH | <input type="checkbox"/> pH (patient temp.) |
| <input checked="" type="checkbox"/> PO2 | <input type="checkbox"/> PO2 (patient temp.) |
| <input type="checkbox"/> PT | <input type="checkbox"/> PT PET |
| <input type="checkbox"/> PT Ratio | <input checked="" type="checkbox"/> sO2 |
| <input type="checkbox"/> TCO2 | <input type="checkbox"/> U/Cre |
| <input type="checkbox"/> Urea | |

 1 of 1 

* Panel Name is a **required** entry field.

3. Click the **Save** button to save the edited data, or Click the **Cancel** button to cancel the edit action without saving.

Delete a Test Panel

To **delete** a single test panel or multiple test panels:


1. Click the box(es) in the Panel Name column.
2. Click the **Delete** action button.

To **delete** all test panels:

1. Click the box marked **All** in the Panel Name column header.
2. Click the **Delete** action button.

Warning Messages

If a different menu option is selected before saving new or edited panel information, the following warning message is displayed:




- Click the **OK** button to cancel the data entered.
- Click the **Cancel** button to clear the message and return to the **Panels** screen to save the data.

If a **required** field (Panel Name) is missing when test panel information is saved, the following warning message is displayed:




- Click the **OK** button to clear the message.

If records are not selected from the table before clicking on the **Edit** or **Delete** action button, the following warning message is displayed:




- Click the **OK** button to clear the message and make a selection.
- Click the **OK** button to clear the message and enter the required information.

If the Panel Name entered already exists in the database, the following warning message is displayed:



- Click the **OK** button and choose a different Panel name.

If a test panel is selected for deletion, the following warning message is displayed:



- Click the **OK** button to delete the panels.
- Click the **Cancel** button to cancel the delete action.

10. Technical Support Information

Contacting Abbott Service

Abbott is committed to helping you resolve any problems with the PrecisionWeb Point of Care Data Management System. For assistance, please call Abbott Technical Support at 1-877-529-7185, twenty-four (24) hours per day. Outside the United States and Canada, contact your local office or distributor listed below.

ASIA

Australia and New Zealand

Abbott Australasia
Abbott Diabetes Care
Victoria, Australia
Australia Tel: +1-800-801-478
New Zealand Tel: +0800-106-100

China

Abbott Laboratories S.A.
Shanghai Representative Office
Abbott Diabetes Care
Shanghai, China
Tel: +800-820-3959

Hong Kong

Abbott Laboratories Ltd.
Abbott Diabetes Care
North Point, Hong Kong
Tel: +852-2806-4488

Indonesia

Abbott Indonesia
Jakarta Selatan, Indonesia
Tel: +62-21-52961529

Japan

Abbott Japan Co., Ltd.
Abbott Diabetes Care
Tokyo, Japan
Tel: +047 386 3173

Korea

Abbott Korea Ltd.
Abbott Diabetes Care
Seoul, South Korea
Tel: +080-014-5757

Malaysia

Abbott Laboratories (M) Sdn. Bhd.
Shah Alam, Selangor, Malaysia
Tel: +00-603-5569-1919

Pakistan

Abbott Laboratories (Pak) Ltd.
Abbott Diabetes Care
Karachi, Pakistan
Tel: 0092-21-5100-321

Philippines

Abbott Laboratories Philippines
Metro Manila, Philippines
Tel: +632-6890-495

Singapore

Abbott Laboratories (S) Pte Ltd
Abbott Diabetes Care
Singapore
Tel: +65-6272-2881

Taiwan

Abbott Laboratories LTD
Abbott Diabetes Care
Taipei, Taiwan
Tel: +0800-521-125

Thailand

Abbott Laboratories Ltd
Bangkok, Thailand
Tel: +66-2-252-2004

EUROPE

Austria

Abbott Ges.m.b.H
Abbott Diabetes Care
Vienna, Austria
Tel: +0800-930-093

Belgium and Luxembourg

Abbott S.A./N.V.
Abbott Diabetes Care
Louvain-la-Neuve, Belgium
Belgium Tel: +0800 167-72
Luxembourg Tel: +32 10 475 480

*EUROPE (Continued)**Czech Republic*

Abbott Laboratories s.r.o
Abbott Diabetes Care
Prague, Czech Republic
Tel: +800-189-564

Denmark

Abbott Laboratories A/S
Abbott Diabetes Care
Gentofte, Denmark
Tel: +80 81 53 54

Finland

Abbott Oy
Abbott Diabetes Care
Espoo, Finland
Tel: +0800-555-500

France

Abbott France S.A.S.
Abbott Diabetes Care
Rungis Cedex, France
Tel: +0800-10-11-56

Germany

Abbott GmbH & Co. KG
Abbott Diabetes Care
Wiesbaden, Germany
Tel: +0800-519-9519

Greece

Abbott Laboratories (Hellas) S.A.
Abbott Diabetes Care
Athens, Greece
Tel: Athens +30-210-998-5220
Thessaloniki +30-231-047-2947

Ireland

Abbott Laboratories Ireland, Ltd.
Abbott Diabetes Care
Dublin, Ireland
Tel: +1-800-776-633

Italy

Abbott S.r.l.
Abbott Diabetes Care
Roma, Italia
Tel: +800-334-216

Netherlands

Abbott B.V.
Abbott Diabetes Care
Amersfoort, Netherlands
Tel: +0800-022-8828

Norway

Abbott Norge AS
Abbott Diabetes Care
Fornebu, Norway
Tel: +800-87-100

Portugal

Abbott Laboratórios, Lda
Abbott Diabetes Care
Amadora, Portugal
Tel: +800-200-891

Spain

Abbott Científica S.A.
Abbott Diabetes Care
Madrid, Spain
Tel: +900-300-119

Sweden

Abbott Scandinavia AB
Abbott Diabetes Care
Solna, Sweden
Tel: +020-190-11-11

Switzerland

ABBOTT AG
Abbott Diabetes Care
Baar, Switzerland
Tel: +08 00 / 804 404

United Kingdom

Abbott Laboratories Ltd.
Abbott Diabetes Care
Maidenhead, Berkshire, UK
Tel: +0500-467-466

LATIN AMERICA

Argentina

Abbott Laboratorios de Argentina SA
Abbott Diabetes Care
Ciudad Autonoma De Buenos Aires,
Argentina
Tel: +0800-800-6334

Brazil

Abbott Laboratorios do Brasil Ltda
Abbott Diabetes Care
São Paulo, Brasil
Tel: +0800-703-0128

Chile

Abbott Laboratorios de Chile Ltda.
Abbott Diabetes Care
Santiago, Chile
Tel: +800-802-226

Columbia

Abbott Laboratorios de Colombia SA
Abbott Diabetes Care
Bogotá, Colombia
Tel: +01-8000-522268

Dominican Republic

Abbott Laboratorios Republica Dominicana
Abbott Diabetes Care
Santo Domingo, Dominican Republic
Tel: +809-542-7181

Guatemala

Abbott Laboratorios de Guatemala
Abbott Diabetes Care
Europlaza Torre 1 Nivel 10, Guatemala
Tel: 597-4109

Mexico

Abbott Laboratorios de Mexico S.A. de CV
Abbott Diabetes Care
Colonia del Valle, Mexico
Tel: +01-800-711-2208

Peru

Abbott Laboratorios de Peru
Abbott Diabetes Care
Lima, Peru
Tel: +0-800-1-1215

Uruguay

Abbott Laboratorios Uruguay SA
Abbott Diabetes Care
Montivideo, Uruguay
Tel: +000-405-4342

Venezuela

Abbott Laboratories C.A.
Abbott Diabetes Care
Caracas, Venezuela
Tel: +0-800-Glucosa

MIDDLE EAST AND AFRICA

Israel

ILEX Medical Ltd.
Rosh-Ha'ayin, Israel
Tel: +972- 03-938-5555

Saudi Arabia

Abbott Diabetes Care
Mediserv
Riyadh, Saudi Arabia
Tel: +800-124-1238

South Africa

ADC South Africa
Abbott Laboratories
Gauteng, South Africa
Tel: +27-11-858-2000

Turkey

Abbott Laboratuarlari
Abbott Diabetes Care
Istanbul, Turkey
Tel: +800-261-7504

NORTH AMERICA

Canada

Abbott Laboratories, Ltd.
Abbott Diabetes Care
Kanata, Ontario, Canada
Tel: +800-363-3367

Puerto Rico

Abbott Laboratories Puerto Rico Inc.
Abbott Diabetes Care
San Juan, Puerto Rico, USA
Tel: +1-787-750-5454

USA

Abbott Diabetes Care Inc.
1360 South Loop Road
Alameda, CA USA
Tel: +1-877-529-7185